

1.1x Pre-Aliquoted ReddyMix™ PCR Master Mix, 25µl Reaction

Description: ReddyMix™ PCR Master Mix is a pre-aliquoted ready-to-use master mix. It is a convenient way of amplifying DNA fragments without the need to thaw individual components, reducing the risk of contamination and pipetting errors. The Thermoprime Plus DNA Polymerase, dNTPs, reaction buffer and magnesium chloride are all present in the mix. ReddyMix™ Master Mix also contains a dye and precipitant to facilitate gel loading.

Ordering Information: Select the consumable type and size from Table 1 and MgCl₂ concentration from Table 2.

Cat. No.	Description	Quantity
AM-266-2MM/LD/a	0.2ml Thermo-Strips & Domed Caps	25µl PCR x 12 strips of 8 tubes
AM-266-2MM/LD/b	0.2ml Thermo-Strips & Domed Caps	25µl PCR x 10 x 12 strips of 8 tubes
AM-337-2MM/LD/a	Individual 0.2ml Domed Cap Thermo-Tubes	25µl x 96 tubes
AM-337-2MM/LD/b	Individual 0.2ml Domed Cap Thermo-Tubes	25µl x 10 x 96 tubes
AM-350-2MM/LD/a	Individual 0.5ml Flat Cap Thermo-Tubes	25µl x 96 tubes
AM-350-2MM/LD/b	Individual 0.5ml Flat Cap Thermo-Tubes	25µl x 10 x 96 tubes
AM-384-2MM/LD	Thermo-Fast® 384 Plates	25µl x 10 plates of 384 wells
AM-489-2MM/LD/a	Individual 0.5ml Domed Cap Thermo-Tubes	25µl x 96 tubes
AM-489-2MM/LD/b	Individual 0.5ml Domed Cap Thermo-Tubes	25µl x 10 x 96 tubes
AM-600-2MM/LD	0.2ml Thermo-Fast® 96 Plates	25µl x 10 plates of 96 wells
AM-620-2MM/LD/a	Individual 0.2ml Flat Cap Thermo-Tubes	25µl x 96 tubes
AM-620-2MM/LD/b	Individual 0.2ml Flat Cap Thermo-Tubes	25µl x 10 x 96 tubes
AM-700-2MM/LD	0.2ml Low Profile Thermo-Fast® 96 Plates	25µl x 10 plates of 96 wells
AM-772-2MM/LD/a	0.2ml Low Profile Thermo-Strips & Domed Caps	25µl PCR x 12 strips of 8 tubes
AM-772-2MM/LD/b	0.2ml Low Profile Thermo-Strips & Domed Caps	25µl PCR x 10 x 12 strips of 8 tubes
AM-773-2MM/LD/a	0.2ml Low Profile Thermo-Strips & Flat Caps	25µl PCR x 12 strips of 8 tubes
AM-773-2MM/LD/b	0.2ml Low Profile Thermo-Strips & Flat Caps	25µl PCR x 10 x 12 strips of 8 tubes
AM-800-2MM/LD	0.2ml Skirted Thermo-Fast® 96 Plates	25µl x 10 plates of 96 wells
AM-848-2MM/LD/a	0.2ml Low Profile Thermo-Strips & Domed Caps	25µl PCR x 8 strips of 12 tubes
AM-848-2MM/LD/b	0.2ml Low Profile Thermo-Strips & Domed Caps	25µl PCR x 10 x 8 strips of 12 tubes
AM-849-2MM/LD/a	0.2ml Low Profile Thermo-Strips & Flat Caps	25µl PCR x 8 strips of 12 tubes
AM-849-2MM/LD/b	0.2ml Low Profile Thermo-Strips & Flat Caps	25µl PCR x 10 x 8 strips of 12 tubes
AM-900-2MM/LD	0.2ml Semi Skirted Thermo-Fast® 96 Plates	25µl x 10 plates of 96 wells
AM-937-2MM/LD	Thermo-Fast® 384 Mark II Plates	25µl x 10 plates of 384 wells
AM-990-2MM/LD	0.2ml Ultra Rigid Semi Skirted Thermo-Fast® 96 Plates	25µl x 10 plates of 96 wells
AM-1000-2MM/LD	0.2ml Ultra Rigid Skirted Thermo-Fast® 96 Plates	25µl x 10 plates of 96 wells
AM-1100-2MM/LD	0.2ml Thermo-Fast® 96 PCR Detection Plate	25µl x 10 plates of 96 wells
AM-1300-2MM/LD	Thermo-Fast® 96 Robotic PCR Plate, Capped	25µl x 10 plates of 96 wells
AM-1400-2MM/LD	Thermo-Fast® 96 PCR Detection Plate Mark II, Capped	25µl x 10 plates of 96 wells

For research purposes only

ABgene House, Blenheim Road, Epsom, Surrey, KT19 9AP, UK

Tel: + 44 (0) 1372 723456, Fax: + 44 (0) 1372 741414, Email: info@abgene.com, Website: www.abgene.com

ABgene® is a registered trademark of Advanced Biotechnologies Ltd.

Table 1 (continued):

Cat. No.	Description	Quantity
AM-1111-2MM/LD	Thermo-Fast [®] 384 Diamond PCR Plate	25µl x 10 plates of 384 wells
AM-1113-2MM/LD/a	0.2ml Thermo-Strips with Domed Caps	25µl PCR x 8 strips of 12 tubes
AM-1113-2MM/LD/b	0.2ml Thermo-Strips with Domed Caps	25µl PCR x 10 x 8 strips of 12 tubes
AM-1114-2MM/LD/a	0.2ml Thermo-Strips with Flat Caps	25µl PCR x 8 strips of 12 tubes
AM-1114-2MM/LD/b	0.2ml Thermo-Strips with Flat Caps	25µl PCR x 10 x 8 strips of 12 tubes
AM-1182-2MM/LD/a	0.2ml Thermo-Strips with Flat Caps	25µl PCR x 12 strips of 8 tubes
AM-1182-2MM/LD/b	0.2ml Thermo-Strips with Flat Caps	25µl PCR x 10 x 12 strips of 8 tubes

Note: MM denotes Magnesium Chloride concentration

Table 2:

Cat. No.	MgCl ₂ Concentration
AM-xxx-215/LD	1.5mM
AM-xxx-220/LD	2.0mM
AM-xxx-225/LD	2.5mM
AM-xxx-230/LD	3.0mM
AM-xxx-235/LD	3.5mM
AM-xxx-240/LD	4.0mM

Note: xxx denotes plastic consumable type

Example: **AM-337-225/LD/a**
1.1x ReddyMix™ PCR Master Mix with 2.5mM MgCl₂ (Reaction Volume 25µl) Pre-aliquoted into 96 x 0.2ml Domed Cap Thermo-Tubes

Kit Components: Each well contains 22.5µl of a 1.1x working concentration ReddyMix™ PCR Master Mix. The addition of the template and primers (in a volume of 2.5µl) results in a final reaction volume of 25µl, containing:

0.625 units Thermoprime Plus DNA Polymerase
 75mM Tris-HCl (pH 8.8 at 25°C)
 20mM (NH₄)₂SO₄
 MgCl₂ (see table)
 0.01% (v/v) Tween[®] 20
 0.2mM each of dATP, dCTP, dGTP, dTTP

Thermo-Fast[®] 96 plates are provided capped with Domed Cap Strips. An extra set of caps for application after the addition of template and primers is included with both plates and Thermo-Strips. Thermo-Fast[®] 384 plates are supplied heat sealed with Easy Peel (or Diamond Foil Seal for Diamond plates).

For research purposes only



Enzyme Source: *Thermus aquaticus*

Storage Conditions: Store at -20°C until ready for use for up to 1 year. Shipped on ice.

Associated Activities: The enzyme has a 5' to 3' polymerisation-dependent exonuclease replacement activity but lacks a 3' to 5' exonuclease activity.

Tip: The gel precipitant in ReddyMix™ Master Mix causes a slight increase in the thermal mass of the reaction mix. In a small number of cases this may necessitate some minor re-optimisation of the thermal cycler programme. If this is the case we suggest increasing the temperature of the denaturation step by 1–2°C and decreasing the temperature of the annealing step by 1–2°C. Alternatively, increase the duration of each step by 5–10 seconds.

A license under the foreign counterparts of U.S. Patents Nos. 4,683,202, 4,683,195 and 4,965,188, owned by F. Hoffman-La Roche Ltd, for use in research and development, has an up-front fee component and a running royalty component. The purchase price of this product includes limited, non-transferable rights under the running-royalty component to use only this amount of the product to practice the Polymerase Chain Reaction ("PCR") and related processes described in said patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license. These rights under the up-front fee component may be purchased from Perkin-Elmer or obtained by purchasing an Authorized Thermal Cycler. No right to perform or offer commercial services of any kind using PCR, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at the Perkin-Elmer Corporation, 850 Lincoln Centre Drive, Foster City, California 94404

For research purposes only