Part 3 LC (LC/MS) カラム

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

For more information on the HPLC analysis of biomolecules, please request our Technical Guide TG20026.

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| Phase | Particle Sizes | Pore Sizes | Stationary Phase Chemistries | General Description |
|--|--------------------|----------------------------------|---|---|
| Hypersil GOLD™ | 1.9 μm, 3 μm, 5 μm | 175 Å | C18 selectivity | Outstanding peak shape using generic gradients with C18 selectivity, providing increased peak capacity, improved resolution, sensitivity and signal to noise. |
| Hypersil GOLD C8 | 1.9 μm, 3 μm, 5 μm | 175 Å | C8 | Offers similar selectivity to Hypersil GOLD but with less retention. |
| Hypersil GOLD aQ™ | 1.9 μm, 3 μm, 5 μm | 175 Å | C18 polar endcapped | Excellent for polar compounds. Good results with low buffer concentrations. Can be used for challenging reverse phase separations employing highly aqueous mobile phases. |
| Hypersil GOLD PFP | 1.9 μm, 3 μm, 5 μm | 175 Å | Perfluorinated phenyl | Offers alternative selectivity in reverse phase applications, particularly for halogenated analytes. |
| Hypersil GOLD CN | 1.9 μm, 3 μm, 5 μm | 175 Å | Cyano | Offers alternative selectivity. Can be used for both reversed and normal phase separations. |
| Hypersil GOLD Phenyl | 1.9 μm, 3 μm, 5 μm | 175 Å | Phenyl | Offers unique selectivity for the analysis of aromatic and moderately polar compounds. |
| BioBasic™ Reversed Phase | 5 μm | 300 Å | C18, C8, C4, CN, Phenyl | Based on 300 Å pore size silica specifically designed for the separation of protein and peptides with increasing hydrophobicity: cyano and phenyl phases to provide alternative selectivity where required. |
| BioBasic Ion Exchange Phase | 5 μm | 300 Å | SCX, AX | Large pore size for biomolecules. BioBasic SCX and AX stationary phases comprise silica particles coated with polymeric ion exchange ligands, which shield proteins from adsorbing to the silica surface. |
| BioBasic SEC | 5 μm | 60 Å, 120 Å, 300 Å, 1000 Å | SEC | BioBasic SEC columns are available with 60 Å, 120 Å, 300 Å and 1000 Å pore sizes, allowing separation of a wide range of molecular weights. |
| Hyperscarb™ | 3 μm, 5 μm | 250 Å | Porous Graphitic Carbon | Unique 100% porous graphitic carbon. Exceptional retention of very polar analytes. Separation of structurally similar substances. pH stable from 0 - 14. Ideal for high temperature applications. |
| BetaBasic™ | 3 μm, 5 μm | 150 Å | C18, C8, C4, CN, Phenyl | General purpose packing. Superb high pH stability. |
| BETASIL™ | 3 μm, 5 μm 5 μm | 100 Å | C18, C8, C6, Phenyl/Hexyl, Cyano, Silica C1, Diol | General purpose packing with high surface coverage for strong retention and use with high organic mobile phases. |

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| USP Code | Description | Recommended Phase |
|----------|--|---|
| L1 | Octadecyl silane (C18) chemically bonded to porous or ceramic microparticles, 1.5 - 10 µm in diameter, or a monolithic rod | Hypersil GOLD™ |
| | | Hypersil GOLD aQ™ |
| | | Hypersil GOLD 1.9 µm |
| | | Hypersil GOLD aQ 1.9 µm |
| | | BioBasic™ 18 |
| | | Hypersil BDS C18 |
| L2 | Octadecyl silane chemically bonded to a silica gel of controlled surface porosity bonded to a solid spherical core, 30 - 50 µm in diameter | Pellicular ODS |
| L3 | Porous silica microparticles, 5 - 10 µm in diameter | BETASIL™ Silica |
| | | Hypersil Silica |
| L4 | Silica gel of controlled surface porosity bonded to a solid spherical core, 30 - 50 µm in diameter | Pellicular Silica |
| L7 | Octyl silane (C8) chemically bonded to totally porous silica particles, 1.5 - 10 µm in diameter | Hypersil GOLD C8 |
| | | Hypersil GOLD C8 1.9 µm |
| | | BioBasic 8 |
| | | Hypersil BDS C8 |
| | | Hypersil MOS |
| | | Hypersil MOS-2 |
| L8 | An essentially monomolecular layer of aminopropylsilane chemically bonded to totally porous silica gel support, 3 - 10 µm in diameter | Hypersil APS-2 |
| L9 | Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation exchange coating (SCX), 3 - 10 µm in diameter | Partisil™ SCX |
| L10 | Nitrile groups (CN) chemically bonded to porous silica particles, 3 - 10 µm in diameter | Hypersil GOLD CN |
| | | BioBasic CN |
| | | Hypersil BDS Cyano |
| | | Hypersil CPS |
| | | Hypersil CPS-2 |
| L11 | Phenyl groups chemically bonded to porous silica particles, 1.5 - 10 µm in diameter | Hypersil GOLD Phenyl |
| | | Hypersil GOLD Phenyl 1.9 µm |
| | | BioBasic Phenyl |
| | | Hypersil BDS Phenyl |
| | | Hypersil Phenyl |
| | | Hypersil Phenyl-2 |
| L13 | Trimethylsilane chemically bonded to porous silica particles, 3 - 10 µm in diameter | Hypersil SAS |
| | | BETASIL C1 |
| L14 | Silica gel having a chemically bonded, strongly basic quaternary ammonium anion exchange (SAX) coating, 5 - 10 µm in diameter | Partisil SAX |
| | | Hypersil SAX (5 µm) |
| L15 | Hexylsilane (C6) chemically bonded to totally porous silica particles, 3 - 10 µm in diameter | BETASIL C6 |
| L17 | Strong cation exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 7 - 11 µm in diameter | HyperREZ™ XP Carbohydrate H ⁺ HyperREZ XP Organic Acids |
| L19 | Strong cation exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the calcium form, about 9 µm diameter | HyperREZ XP Carbohydrate Ca ²⁺ |
| | | HyperREZ XP Sugar Alcohols |
| L20 | Dihydroxypropane groups chemically bonded to porous silica particles, 5 - 10 µm in diameter | BETASIL Diol |
| L21 | A rigid spherical styrene-divinylbenzene copolymer, 5 - 10 µm in diameter | HyperREZ XP RP100 |
| | | HyperREZ XP RP300 |
| L22 | A cation exchange resin made of porous polystyrene gel with sulfonic acid groups, about 10 µm in size | HyperREZ XP SCX |
| L25 | Packing having the capacity to separate compounds with a molecular weight range 100 - 5000 (as determined by polyethylene oxide) applied to neutral, anionic and cationic water-soluble polymers. | HyperGEL™ AP |
| L26 | Butyl silane (C4) chemically bonded to totally porous silica particles, 3 - 10 µm in diameter | BioBasic 4 |
| | | BetaBasic™ 4 |
| L27 | Porous silica particles, 30 - 50 µm in diameter | HyperPrep™ Silica |
| L33 | Packing having the capacity to separate dextrans by molecular size over a range of 4,000 to 500,000 daltons. It is spherical, silica-based, and processed to provide pH stability | BioBasic SEC 120 |
| | | BioBasic SEC 300 |
| | | BioBasic SEC 1000 |
| L34 | Strong cation exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the lead form, about 9 µm in diameter | HyperREZ XP Carbohydrate Pb ²⁺ |
| L43 | Pentafluorophenyl groups chemically bonded to silica particles by a propyl spacer, 5 - 10 µm in diameter | Hypersil GOLD PFP |
| | | Fluophase™ PFP |
| L52 | A strong cation exchange resin made of porous silica with sulfopropyl groups by a propyl spacer, 5 - 10 µm in diameter | BioBasic SCX |
| L58 | Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the sodium form, about 7 to 11 µm in diameter | HyperREZ XP Carbohydrate Na ⁺ |
| L59 | Packing having the capacity to separate proteins by molecular weight over the range of 10 to 500 kDa. It is spherical (10 µm), silica-based, and processed to provide hydrophilic characteristics and pH stability | BioBasic SEC 300 (5 µm) |
| L60 | Spherical, porous silica gel, 10 µm or less in diameter, the surface of which has been covalently modified with alkyl amide groups and endcapped | HyPURITY ADVANCE™ |

* These are the recommended Thermo Scientific HPLC columns for various USP categories although other columns for each category are also available.

Thermo Scientific 製 LCカラム

The tables below list Thermo Scientific HPLC sorbents offered. Please also refer to the Advanced User Graphic (AUG) for each HPLC phase on the pages indicated.

| Phase | Particle Type | Particle Size (µm) | Pore Size (Å) | Nominal Surface Area (m ² /g) | % Carbon | Endcapping | USP Code | Phase Code |
|-------------------------------------|------------------------------------|--------------------|---------------|--|----------|------------|----------|------------|
| AQUASIL | | | | | | | | |
| C18 | spherical, silica | 3, 5 | 100 | 310 | 12 | polar | L1 | 775 |
| BetaBasic™ | | | | | | | | |
| 18 | spherical, silica | 3, 5 | 150 | 200 | 13 | Yes | L1 | 715 |
| 8 | spherical, silica | 3, 5 | 150 | 200 | 7 | Yes | L7 | 714 |
| 4 | spherical, silica | 3, 5 | 150 | 200 | 6 | Yes | L26 | 716 |
| Phenyl | spherical, silica | 3, 5 | 150 | 200 | 7 | Yes | L11 | 718 |
| CN | spherical, silica | 3, 5 | 150 | 200 | 5 | Yes | L10 | 717 |
| BetaMax™ | | | | | | | | |
| Neutral | spherical, silica | 5 | 60 | 540 | 29 | Yes | L1 | 950 |
| Acid | spherical, silica | 5 | 60 | 540 | 15 | Yes | L28 | 952 |
| Base | spherical, silica | 5 | 60 | 540 | 9 | Yes | L10 | 951 |
| BETASIL™ | | | | | | | | |
| C18 | spherical, silica | 3, 5, 10 | 100 | 310 | 20 | Yes | L1 | 701 |
| C8 | spherical, silica | 3, 5, 10 | 100 | 310 | 12 | Yes | L7 | 702 |
| C6 | spherical, silica | 3, 5 | 100 | 310 | 11 | Yes | L15 | 703 |
| C1 | spherical, silica | 5 | 100 | 310 | 4 | Yes | L13 | 705 |
| Phenyl | spherical, silica | 3, 5 | 100 | 310 | 11 | Yes | L11 | 706 |
| Phenyl-Hexyl | spherical, silica | 3, 5 | 100 | 310 | 11 | Yes | L11 | 730 |
| CN | spherical, silica | 3, 5 | 100 | 310 | 6 | Yes | L10 | 708 |
| Silica 100 | spherical, silica | 3, 5 | 100 | 310 | – | – | L3 | 700 |
| Diol 100 | spherical, silica | 5 | 100 | 310 | 6 | – | L20 | 726 |
| BioBasic™ | | | | | | | | |
| 18 | spherical, silica | 5 | 300 | 100 | 9 | Yes | L1 | 721 |
| 8 | spherical, silica | 5 | 300 | 100 | 5 | Yes | L7 | 722 |
| 4 | spherical, silica | 5 | 300 | 100 | 4 | Yes | L26 | 723 |
| Phenyl | spherical, silica | 5 | 300 | 100 | 3 | Yes | L11 | 724 |
| CN | spherical, silica | 5 | 300 | 100 | 3.5 | Yes | L10 | 729 |
| AX | spherical, silica | 5 | 300 | 100 | 3 | No | – | 731 |
| SCX | spherical, silica | 5 | 300 | 100 | 3 | – | L52 | 732 |
| DELTABOND™ | | | | | | | | |
| Resolution AK | spherical, polymer coated silica | 5 | 300 | 100 | 12 | No | – | 322 |
| Fast AK | spherical, polymer coated silica | 5 | 300 | 100 | 12 | No | – | 323 |
| Fluophase™ | | | | | | | | |
| RP | spherical, silica | 5 | 100 | 310 | 10 | Yes | – | 825 |
| PFP | spherical, silica | 5 | 100 | 310 | 12 | Yes | L43 | 827 |
| WP | spherical, silica | 5 | 300 | 100 | 5 | Yes | – | 826 |
| Hypercarb™ | | | | | | | | |
| Hypercarb | spherical, porous graphitic carbon | 3, 5, 7 | 250 | 120 | 100 | – | – | 350 |
| HyperREZ™ XP | | | | | | | | |
| Carbohydrate H⁺ | spherical, polymer | 8 | – | – | – | – | L17 | 690 |
| Carbohydrate Pb²⁺ | spherical, polymer | 8 | – | – | – | – | L34 | 691 |
| Carbohydrate Ca²⁺ | spherical, polymer | 8 | – | – | – | – | L19 | 692 |
| Carbohydrate Na⁺ | spherical, polymer | 10 | – | – | – | – | – | 693 |
| Organic Acid | spherical, polymer | 8 | – | – | – | – | L17 | 696 |
| Sugar Alcohol | spherical, polymer | 8 | – | – | – | – | L19 | 697 |

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| Phase | Manufacturer | Pore Size (Å) | Area (m ² /g) | % C | Recommended Thermo Scientific Alternative |
|-------------------------|--------------|---------------|--------------------------|------|--|
| ACE C18 | ACT | 100 | 300 | 15.5 | Hypersil GOLD™ |
| ACE C8 | ACT | 100 | 300 | 9 | Hypersil GOLD C8 |
| ACE C4 | ACT | 100 | 300 | 5.5 | HyPURITY™ C4 |
| ACE CN | ACT | 100 | 300 | 5.5 | Hypersil GOLD CN |
| ACE Phenyl | ACT | 100 | 300 | 9.5 | Hypersil GOLD Phenyl |
| ACE AQ | ACT | 100 | 300 | 14 | Hypersil GOLD aQ™ |
| ACE C18-300 | ACT | 300 | 100 | 9 | BioBasic™ 18 |
| ACE C8-300 | ACT | 300 | 100 | 5 | BioBasic 8 |
| ACE C4-300 | ACT | 300 | 100 | 2.6 | BioBasic 4 |
| ACE CN-300 | ACT | 300 | 100 | 2.6 | BioBasic CN |
| ACE Phenyl-300 | ACT | 300 | 100 | 5.3 | BioBasic Phenyl |
| ACQUITY UPLC BEH C18 | Waters | 130 | 185 | - | Hypersil GOLD (1.9 µm) |
| ACQUITY UPLC BEH C8 | Waters | 130 | 185 | - | Hypersil GOLD C8 (1.9 µm) |
| ACQUITY UPLC BEH Phenyl | Waters | 130 | 185 | - | Hypersil GOLD Phenyl (1.9 µm) |
| ACQUITY UPLC HSS T3 | Waters | 100 | 230 | - | Hypersil GOLD aQ (1.9 µm) |
| Alltima™ HP C18 | Grace | 190 | 200 | 12 | Hypersil GOLD |
| Alltima HP C18 Amide | Grace | 190 | 200 | 12 | HyPURITY ADVANCE™ |
| Alltima HP C18 AQ | Grace | 100 | 450 | 20 | Hypersil GOLD aQ |
| Alltima HP C18 HiLoad | Grace | 100 | 450 | 24 | BETASIL™ C18 |
| Alltima HP C18-Amide | Grace | 190 | 200 | 12 | HyPURITY ADVANCE |
| Alltima HP C8 | Grace | 190 | 200 | 8 | Hypersil GOLD C8 |
| Alltima HP CN | Grace | 190 | 200 | 4 | Hypersil GOLD CN |
| Alltima HP Silica | Grace | 190 | 200 | - | BETASIL Silica |
| Aminex™ HPX42C | Bio-Rad | - | - | - | HyperREZ™ XP Carbohydrate Ca ²⁺ |
| Aminex HPX72S | Bio-Rad | - | - | - | HyperREZ XP Carbohydrate H ⁺ |
| Aminex HPX87C | Bio-Rad | - | - | - | HyperREZ XP Carbohydrate Ca ²⁺ |
| Aminex HPX87H | Bio-Rad | - | - | - | HyperREZ XP Carbohydrate H ⁺ |
| Aminex HPX87N | Bio-Rad | - | - | - | HyperREZ XP Carbohydrate Na ⁺ |
| Aminex HPX87P | Bio-Rad | - | - | - | HyperREZ XP Carbohydrate Pb ²⁺ |
| AQUA™ C18 | Phenomenex | 125 | 320 | 15 | Hypersil GOLD aQ |
| Ascentis C18 | Supelco | 100 | 450 | 25 | Hypersil GOLD BETASIL C18 |
| Ascentis C8 | Supelco | 100 | 450 | 15 | Hypersil GOLD C8 BETASIL C8 |
| Ascentis RP-Amide | Supelco | 100 | 450 | 19.5 | HyPURITY ADVANCE BetaMax™ Acid |
| Ascentis Phenyl | Supelco | 100 | 450 | 19 | Hypersil GOLD Phenyl BETASIL Phenyl |
| Atlantis™ dC18 | Waters | 100 | 330 | 12 | Hypersil GOLD aQ |
| Columbus™ C18 | Phenomenex | 110 | 375 | 19 | Hypersil GOLD BETASIL C18 |
| Columbus C8 | Phenomenex | 110 | 375 | 13 | Hypersil GOLD C8 |

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