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Nutrition Facts	
Serving Size 1 cup (228g)	
Servings Per Container 2	
Amount Per Serving	
Calories 260	Calories from Fat 120
%Daily Value*	
Total Fat 13g	20%
Saturated Fat 5g	25%
Trans Fat 2g	
Cholesterol 30mg	10%
Sodium 660mg	28%
Total Carbohydrate 31	10%
Dietary Fiber 0g	0%
Sugars 5g	
Protein 5g	
Vitamin A 4%	Vitamin C 2%
Calcium 15%	Iron 4%
*Percent Daily Values are based on a 2,000 calorie diet. Your Daily Values may be higher or lower depending on your calorie needs:	
	Calories: 2,000 2,500
Total Fat	Less than 65g 80g
Sat Fat	Less than 20g 25g
Cholesterol	Less than 300g 300g
Sodium	Less than 2,400mg 2,400 mg
Total Carbohydrate	300g 375g
Dietary Fiber	25g 30g
Calories per gram:	
Fat 9	Carbohydrate 4 Protein 4

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1. FDA Web site: <http://www.fda.gov/oc/initiatives/transfat/>
2. Copenhagen Post On-Line, Denmark adopts limits for trans-fatty acid, 11-03-2003, <http://www.cphpost.dk/get/65922.html>
3. Rapid Determination of Isolated *trans* fat Geometric Isomers in Fats and Oils by Attenuated Total Reflectance, AOCS Official Method Cd 14d-99, 1999
4. Simultaneous Determination of Iodine Value and *trans* Content of Fats and Oils by Single-Bounce Horizontal Attenuated Total Reflectance Fourier Transform Infrared Spectroscopy, Jacqueline Sedman, F.R. van de Voort, and Ashraf A. Ismail, JAOCS, Vol. 77 no. 4, 2000
5. Official Methods for the Determination of *Trans* Fat, Magdi Mossoba, John K.G. Kramer, Pierluigi Delmonte, Martin P. Yurawecz, Jeanne I. Rader, AOCS Press, 2003

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E-mail
info-jp@thermo.com

www.thermofisher.co.jp
(日本)
www.thermo.com
(グローバル)

©2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specification, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.