

# Sigma Steroids FT-IR Spectral Library

## *Introduction*

The Sigma Steroids Library contains 3011 FT-IR spectra of a wide variety of natural and synthetic steroids. The steroids and spectra were obtained and provided by Sigma Chemical Company. The spectra of steroids not available from Sigma were obtained from samples provided to Sigma by the United Kingdom Medical Research Council.

The spectra are listed in this document in alphabetical order by compound name. The index number of each spectrum appears to the left of the name. If you are using the library with Nicolet's OMNIC spectroscopy software, the molecular formula for each compound is available on-line.

The spectra were measured from KBr pellets at  $2\text{ cm}^{-1}$  digital resolution. The spectra were mathematically deresolved when converted to Nicolet library format.

The following table gives the skeletal classifications for the index number ranges.

<i>Index # range</i>	<i>Skeletal Classification</i>
1-104	Estrane
105-1074	Androstane
1075-1182	Etianic acid
1183-2099	Pregnane
2100-2125	Cardane, Cardenolide
2126-2321	Cholane, Cholanic acid
2322-2682	Cholestane, Spirostane
2683-2949	Ergostane, Lumistane, Campestone
2950-2976	Stigmastane, Poriferastane
2977-3006	Lupane, Ursane, Oleanane, etc.
3007-3011	Steroid analogs

Copyright © 1999 by Nicolet Instrument Corporation, Madison WI 53711. Printed in the United States of America. All world rights reserved. No part of this publication may be stored in a retrieval system, transmitted, or reproduced in any way, including but not limited to photocopy, photograph, magnetic or other record, without the prior written permission of Nicolet Instrument Corporation.

# Sigma Steroids FT-IR Spectral Library

## Index of Compounds

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2547	(24 R,S)-24,25-Dihydroxycholesterol-[26,27-d6]	2199	11-Amino-12-oxo-5b-cholan-24-oic acid methyl ester hydrochloride
2546	(24 R,S)-24,25-Dihydroxycholesterol; (24 R,S)-Cholest-5-ene-3b,24,25-triol	1870	11-Dehydrocorticosterone; 21-Hydroxypregn-4-ene-3,11,20-trione; Kendall's A
2544	(24R)-24,25-Dihydroxycholesterol-[26,27-d6]	626	11-Ketoetiocholanolone sulfate, sodium salt; 3a-Hydroxy-5b-androstane-11,17-dione 3-sulfate, sodium salt
2545	(24S)-24,25-Dihydroxycholesterol; (24S)-Cholest-5-ene-3b,24,25-triol	869	11-Ketotestosterone; 17b-Hydroxy-4-androstene-3,11-dione
2549	(25 R,S)-25,26-Dihydroxycholesterol; (25 R,S)-Cholest-5-ene-3b,25,26-triol	2195	11-Oxo-5b-cholan-24-oic acid methyl ester
2548	(25S)-25,26-Dihydroxycholesterol; (25S)-Cholest-5-ene-3b,25,26-triol	1730	11-Pregnene-3,20-dione
95	1,3,5(10),16-Estratetraene-3,17-diol diacetate	1038	11a,15a-Dihydroxy-3a,5-cyclo-5a-androstane-6,17-dione
97	1,3,5(10),6-Estratetraene-3,17b-diol 3-glucuronide	5	11a,15a-Dihydroxy-5a-estran-3-one
96	1,3,5(10),6-Estratetraene-3,17b-diol diacetate; 6-Dehydroestradiol diacetate	1046	11a,15a-Dihydroxy-A-nor(5->6b)abeo-androstan-5-one
93	1,3,5(10),7-Estratetraene-3,17a-diol; 17 a-Dihydroequilin	415	11a,16a-Dihydroxy-5a-androstan-3-one
94	1,3,5(10),7-Estratetraene-3,17b-diol; 17 b-Dihydroequilin	469	11a,16b-Diacetoxy-5a-androstan-3-one; 11a,16b-Dihydroxy-5a-androstan-3-one 11,16-diacetate
24	1,3,5(10)-Estratrien-3-ol	227	11a,16b-Diacetoxy-5a-androstan-3b-ol; 5a-Androstane-3b,11a,16b-triol 11,16-diacetate
69	1,3,5(10)-Estratriene-2,3,16a,17b-tetrol; 2-Hydroxyestriol	416	11a,16b-Dihydroxy-5a-androstan-3-one
43	1,3,5(10)-Estratriene-2,3,17b-triol; 2-Hydroxyestradiol	419	11a,16b-Dihydroxy-5a-androstan-4-one
71	1,3,5(10)-Estratriene-3,15a,16a,17b-tetrol; 15a-Hydroxyestriol	1823	11a,17a,21-Trihydroxy-4-pregnene-3,20-dione
53	1,3,5(10)-Estratriene-3,16a,17a-triol 3-methyl ether; Epimestrol	2072	11a,17a-Dihydroxy-19-nor-4-pregnene-3,20-dione
64	1,3,5(10)-Estratriene-3,16a,17b-triol 16-glucuronide, sodium salt; Estril 16a-(b-D-glucuronide), sodium salt	1782	11a,17a-Dihydroxy-4-pregnene-3,20-dione
62	1,3,5(10)-Estratriene-3,16a,17b-triol 3-glucuronide, sodium salt; Estril 3-(b-D-glucuronide), sodium salt	22	11a,17b-Diacetoxy-4-estren-3-one
47	1,3,5(10)-Estratriene-3,16a,17b-triol; Estril	951	11a,17b-Dihydroxy-17-methyl-5a-androstan-3-one
48	1,3,5(10)-Estratriene-3,16b,17a-triol; 16,17-Epiestriol	952	11a,17b-Dihydroxy-17-methyl-5b-androstan-3-one
55	1,3,5(10)-Estratriene-3,16b,17b-triol 3-methyl ether	1014	11a,17b-Dihydroxy-17a-methyl-4-androsten-3-one; 11a-Hydroxy-17a-methyltestosterone (approx.95-98%)
25	1,3,5(10)-Estratriene-3,17a-diol; 17a-Estradiol	953	11a,17b-Dihydroxy-4,4-dimethyl-5a-androstan-3-one
40	1,3,5(10)-Estratriene-3,17b-diol 3-benzoate; Estradiol benzoate	1061	11a,17b-Dihydroxy-4-oxa-A-homo-5a-androstan-3-one
70	1,3,5(10)-Estratriene-3,6a,16a,17b-tetrol; 6a-Hydroxyestriol	430	11a,17b-Dihydroxy-5a-androstan-7-one
45	1,3,5(10)-Estratriene-3,6a,17b-triol; 6a-Hydroxyestradiol	1171	11a,19-Diacetoxy-14-hydroxy-3-oxo-14b-eti-4-enic acid ethyl ester
923	1,4-Androstadiene-3,11,17-trione	368	11a-(p-Toluenesulfonyloxy)-5a-androstan-3-one; 11a-Hydroxy-5a-androstan-3-one tosylate
915	1,4-Androstadiene-3,17-dione	391	11a-Acetoxy-16,16-difluoro-5a-androstan-3-one; 16,16-Difluoro-11a-hydroxy-5a-androstan-3-one acetate
1964	1-Dehydroaldosterone; 11b,18-Epoxy-18,21-dihydroxy-1,4-pregnadiene-3,20-dione	1548	11a-Acetoxy-4a,5-epoxy-5a-pregnane-3,20-dione
911	1-Dehydrotestosterone benzoate; 17 b-Benzoyloxyandrosta-1,4-dien-3-one	2717	11a-Acetoxy-5a-ergostan-3-one
2815	10a-Ergosterol; 5a,10a-Ergosta-5,7,22-trien-3b-ol; Pyrocalciferol	1755	11a-Acetoxy-5a-pregn-1-ene-3,20-dione
2141	11,12-Dibromo-5b-cholan-24-oic acid methyl ester	2914	11a-Benzoyloxy-3b-methoxy-9a-methyl-5a-ergost-7-ene; 3b-Methoxy-9a-methyl-5a-ergost-7-en-11a-ol benzoate
2239	11,12-Dioxo-3a-hydroxy-5b-cholan-24-oic acid	2200	11a-Bromo-12-oxo-5b-cholan-24-oic acid
2320	11,12-Seco-5b-cholane-11,12,24-trioic acid	920	11a-Hydroxy-1,4-androstadiene-3,17-dione
1645	11-Acetoxy-3a,20b-dibenzoyloxy-5b-pregn-9(11)-ene; 5b-Pregn-9(11)-ene-3a,11,20b-triol 11-acetate 3,20-dibenzoate	1036	11a-Hydroxy-3a,5-cyclo-5a-androstane-6,17-dione
1920	11-Acetoxy-3a,20b-dibenzoyloxy-5b-pregna-7,9(11)-diene; 5b-Pregna-7,9(11)-diene-3a,11,20b-triol 11-acetate 3,20-dibenzoate	836	11a-Hydroxy-3b-methoxy-5-androsten-17-one; 3b,11a-Dihydroxy-5-androsten-17-one 3-methyl ether
1921	11-Acetoxy-3a,20b-dibenzoyloxy-5b-pregna-8(14),9(11)-diene; 5b-Pregna-8(14),9(11)-diene-3a,11,20b-triol 11-acetate dibenzoate	770	11a-Hydroxy-4-androsten-3-one; 4-Androsten-11a-ol-3-one
2198	11-Amino-12-oxo-5b-cholan-24-oic acid 24->11 lactam	2093	11a-Hydroxy-4-oxa-A-homo-5a-pregnane-3,20-dione
		872	11a-Hydroxy-5a-androst-1-ene-3,17-dione
		600	11a-Hydroxy-5a-androstane-2,17-dione
		613	11a-Hydroxy-5a-androstane-3,16-dione
		601	11a-Hydroxy-5a-androstane-3,6-dione
		604	11a-Hydroxy-5a-androstane-3,7-dione
		623	11a-Hydroxy-5a-androstane-6,17-dione

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2763	11a-Hydroxy-5a-ergost-1-en-3-one	1815	11b,21-Dihydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester 3,20-disemicarbazone
1740	11a-Hydroxy-5a-preg-1-ene-3,20-dione	1794	11b,21-Dihydroxy-4-pregnene-3,20-dione 21-sulfate, potassium salt; Corticosterone 21-sulfate
1306	11a-Hydroxy-5a-pregnan-20-one	1784	11b,21-Dihydroxy-4-pregnene-3,20-dione; Corticosterone; Kendall's B; Reichstein's H
1496	11a-Hydroxy-5a-pregnane-2,20-dione	1607	11b,21-Dihydroxy-5a-pregnan-18-al-3,20-dione; 5a-Dihydroaldosterone
1758	11a-Hydroxyprogesterone 11-hemisuccinate	1608	11b,21-Dihydroxy-5b-pregnan-18-al-3,20-dione
1771	11a-Hydroxyprogesterone 11-hemisuccinate: BSA conjugate	126	11b-Acetoxy-5a-androstane
2084	11a-Methoxy-17a-ethynylestradiol; 11a-Methoxy-19-nor-17a-pregna-1,3,5(10)-trien-20-yne-3,17-diol	2718	11b-Acetoxy-5a-ergostan-3-one
1185	11b,12a-Dibromo-5b-pregnane	2696	11b-Acetoxy-5a-ergostan-3b-ol; 5a-Ergostane-3b, 11b-diol 11-acetate
2180	11b,12b-Epoxy-3a-hydroxy-5b-cholan-24-oic acid methyl ester	2201	11b-Bromo-12-oxo-5b-cholan-24-oic acid
417	11b,16b-Dihydroxy-5a-androstan-3-one	848	11b-Hydroxy-17b-propionoxy-4-androsten-3-one; 11b-17b-Dihydroxy-4-androsten-3-one 17-propionate
1720	11b,17a,20a,21-Tetrahydroxy-4-pregnen-3-one	1739	11b-Hydroxy-4-pregnene-3,20-dione; 11b-Hydroxyprogesterone
1931	11b,17a,20b,21-Tetrahydroxy-1,4-pregnadien-3-one	775	11b-Hydroxy-5a-androst-2-en-17-one; 5a-Androst-2-en-11b-ol-17-one
1721	11b,17a,20b,21-Tetrahydroxy-4-pregnen-3-one; Reichstein's compound E	308	11b-Hydroxy-5a-androstan-3-one; 5a-Androstan-11b-ol-3-one
1946	11b,17a,21-Trihydroxy-1,4-pregnadiene-3,20-dione 21-(disodium phosphate)	614	11b-Hydroxy-5a-androstane-3,16-dione
1952	11b,17a,21-Trihydroxy-1,4-pregnadiene-3,20-dione 21-glucuronide methyl ester	1307	11b-Hydroxy-5a-pregnan-20-one
1951	11b,17a,21-Trihydroxy-1,4-pregnadiene-3,20-dione 21-glucuronide; Prednisolone 21-glucuronide	2147	11b-Hydroxy-5b-cholan-24-oic acid methyl ester
1843	11b,17a,21-Trihydroxy-4-pregnene-3,20-dione 21-glucuronide methyl ester	453	11b-Hydroxyandrosterone; 3a,11b-Dihydroxy-5a-androstan-17-one
1844	11b,17a,21-Trihydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester 3,20-disemicarbazone	84	11b-Hydroxyestrone; 3,11b-Dihydroxy-1,3,5(10)-estratrien-17-one
1842	11b,17a,21-Trihydroxy-4-pregnene-3,20-dione 21-glucuronide, ammonium salt	2085	11b-Methoxy-17a-ethynylestradiol; 11b-Methoxy-19-nor-17a-pregna-1,3,5(10)-trien-20-yne-3,17-diol; Moxestrol
1841	11b,17a,21-Trihydroxy-4-pregnene-3,20-dione 21-glucuronide; Cortisol 21-glucuronide	2176	11x,12x-Diacetoxy-5b-cholan-24-oic acid methyl ester
1829	11b,17a,21-Trihydroxy-4-pregnene-3,20-dione 21-sulfate, sodium salt; Cortisol 21-sulfate, potassium salt	2165	11x,12x-Dihydroxy-5b-cholan-24-oic acid methyl ester
1824	11b,17a,21-Trihydroxy-4-pregnene-3,20-dione; Hydrocortisone; Kendall's F; Reichstein's M	2156	11x,12x-Epoxy-5b-cholan-24-oic acid methyl ester
2034	11b,17a-Dihydroxy-6a-methyl-1,4-pregnadiene-3,20-dione	2215	11x-Bromo-3a-hydroxy-12-oxo-5x-cholan-24-oic acid
1825	11b,18,21-Trihydroxy-4-pregnene-3,20-dione	495	12,12-Difluoro-1b,6a-dihydroxy-5a-androstan-17-one
1858	11b,18-Epoxy-18,21-dihydroxy-17a-preg-4-ene-3,20-dione 18-a-D-glucuronide	285	12,12-Difluoro-5a-androstan-17-one
1859	11b,18-Epoxy-18,21-dihydroxy-4-pregnene-3,20-dione 18-a-D-glucuronide triacetyl methyl ester	282	12,12-Difluoro-5a-androstan-3-one (crude)
1123	11b,18-Epoxy-18-hydroxy-3-oxo-5b-androstane-17b-carboxylic acid 17,18-lactone; 5b-Dihydroaldosterone g-etiolactone	689	12,12-Difluoro-5a-androstane-1,6,17-trione
1860	11b,18-Epoxy-21-acetoxy-18-hydroxy-4-pregnene-3,20-dione 18-a-D-glucuronide triacetyl methyl ester	236	12,12-Difluoro-5a-androstane-1b,6a,17b-triol
1558	11b,18-Epoxy-3a,18,21-trihydroxy-5b-pregnan-20-one; 3a,5b-Tetrahydroaldosterone	593	12,12-Difluoro-5a-androstane-7,17-dione
1700	11b,20b,21-Trihydroxy-4-pregnen-3-one	2518	12-(Trinitrophenylamino)lauryl cholesterol
1671	11b,20b-Dihydroxy-4-pregnen-3-one	2321	12-Oxo-3,4-seco-5b-cholane-3,4,24-trioic acid trimethyl ester
1886	11b,21-Diacetoxy-4-pregnen-18-al-3,20-dione; Aldosterone diacetate	2193	12-Oxo-5b-cholan-24-oic acid
1881	11b,21-Dihydroxy-17a-preg-4-ene-3,18,20-trione; 17a-Isoaldosterone	2196	12-Oxo-5b-cholan-24-oic acid methyl ester
1883	11b,21-Dihydroxy-4-pregnen-18-al-3,20-dione 21-monosulfate ammonium salt	418	12a,16a-Dihydroxy-5a-androstan-3-one
1813	11b,21-Dihydroxy-4-pregnene-3,20-dione 11,21-di(glucuronide triacetyl methyl ester)	1413	12a,21-Diacetoxy-3a-hydroxy-5b-pregnan-20-one
1807	11b,21-Dihydroxy-4-pregnene-3,20-dione 11-glucuronide	1530	12a,21-Diacetoxy-5b-pregnane-3,20-dione
1809	11b,21-Dihydroxy-4-pregnene-3,20-dione 21-glucuronide (1:1 sodium salt + free acid)	2214	12a-Acetoxy-3,3-dimethoxy-5b-cholan-24-oic acid methyl ester
1810	11b,21-Dihydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester	2301	12a-Acetoxy-3-oxo-5b,20R-bisnorcholan-22-oic acid methyl ester
		2300	12a-Acetoxy-3-oxo-5b-bisnorcholan-22-oic acid methyl ester
		2209	12a-Acetoxy-3-oxo-5b-cholan-24-oic acid methyl ester
		1752	12a-Acetoxy-4-pregnene-3,20-dione; 12a-Acetoxyprogesterone
		2216	12a-Acetoxy-4x-bromo-3-oxo-5b-cholan-24-oic acid methyl ester
		127	12a-Acetoxy-5a-androstane
		2259	12a-Acetoxy-5b-chol-3-en-24-oic acid methyl ester
		1508	12a-Acetoxy-5b-pregnane-3,20-dione
		1603	12a-Acetoxy-5b-pregnane-3,7,20-trione
		1099	12a-Acetoxy-6b-hydroxy-5a-etianic acid methyl ester

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1201	12a-Benzoyloxy-5b-pregnane; 5b-Pregnan-12a-ol benzoate	174	15b-Methoxy-5a-androstan-17b-ol
1545	12a-Bromo-3a,21-diacetoxy-5b-pregnane-11,20-dione	192	15b-Methoxy-5a-androstan-17b-ol tosylate
392	12a-Bromo-3a-hydroxy-5b-androstan-11-one	142	16,16-Difluoro-5a-androstan-3b-ol
594	12a-Bromo-5b-androstane-3,11-dione	237	16,16-Difluoro-5a-androstane-3b,7b,12b-triol
1741	12a-Hydroxy-4-pregnene-3,20-dione; 12a-Hydroxyprogesterone	494	16,16-Difluoro-7b,12b-dihydroxy-5a-androstan-3-one
2254	12a-Hydroxy-5b-cholesterol-3-en-24-oic acid methyl ester	270	16,16-Ethylenedioxy-5a-androstane
2144	12a-Hydroxy-5b-cholan-24-oic acid 24->12 lactone	1019	16-Anisylidene-3b-(p-toluenesulfonyloxy)-5a-androstane-11,17-dione
2148	12a-Hydroxy-5b-cholan-24-oic acid methyl ester	1024	16-Anisylidene-3b-hydroxy-5-androsten-17-one
1498	12a-Hydroxy-5b-pregnane-3,20-dione	992	16-Anisylidene-5a-androstan-11-one
2022	12a-Methyl-4-pregnene-3,11,20-trione	993	16-Anisylidene-5a-androstan-17-one
441	12b,14b-Dihydroxy-5a-androstan-15-one	1016	16-Anisylidene-5a-androstane-11,17-dione
852	12b,15a-Dihydroxy-17-fluoro-5a-androst-16-en-3-one	1023	16-Benzylidene-3b-hydroxy-5-androsten-17-one
526	12b,15a-Dihydroxy-3b-ethoxy-5a-androstan-17-one	1005	16-Benzylidene-3b-hydroxy-5a-androstan-17-one
645	12b,15a-Dihydroxy-5a-androstane-3,17-dione	1007	16-Benzylidene-3b-hydroxy-5a-androstan-17-one tetrahydropyran ether
6	12b,15a-Dihydroxy-5a-estran-11-one	994	16-Benzylidene-5a-androstan-17-one semicarbazone
128	12b-Acetoxy-5a-androstane	974	16-Benzylidene-5a-androstan-17b-ol
1665	12b-Acetoxy-5a-pregn-16-en-20-one	973	16-Benzylidene-5a-androstan-3b-ol
327	12b-Hydroxy-5a-androstan-17-one; 5a-Androstan-12b-ol-17-one	971	16-Benzylidene-5a-androstane
309	12b-Hydroxy-5a-androstan-3-one; 5a-Androstan-12b-ol-3-one	1927	16-Dehydropregnenolone oxime; 20-Hydroximinopregna-5,16-dien-3b-ol
690	12b-Hydroxy-5a-androstane-2,6,16-trione	1935	16-Dehydroprogesterone; Pregna-4,16-diene-3,20-dione
1503	12b-Hydroxy-5a-pregnane-7,20-dione	61	16-Epiestriol triacetate; 3,16b,17b-Triacetoxy-1,3,5(10)-estratriene
2220	12x-Bromo-3a,9a-epoxy-11-oxo-5b-cholan-24-oic acid methyl ester	49	16-Epiestriol; 1,3,5(10)-Estratriene-3,16b,17b-triol
1827	14,19,21-Trihydroxy-14b,17a-pregn-4-ene-3,20-dione	885	16-Hydroximino-3b-hydroxy-5-androsten-17-one
1826	14,19,21-Trihydroxy-14b-pregn-4-ene-3,20-dione	586	16-Hydroximino-5a-androstan-17-one; 5a-Androstane-16,17-dione 16-oxime
1785	14,19-Dihydroxy-14b-pregn-4-ene-3,20-dione	999	16-Hydroxymethylene-5a-androstan-17-one
2123	14,19-Dihydroxy-3-oxo-14b,17a-carda-4,20(22)-dienolide	870	16-Ketotestosterone; 17b-Hydroxy-4-androstene-3,16-dione
2100	14,21-Epoxy-3b,5,19-trihydroxy-5b,14b-cardanolide	2026	16-Methyl-16-dehydropregnenolone acetate; 3b-Acetoxy-16-methyl-5,16-pregnadien-20-one
1087	14-Hydroxy-5b,14b-etianic acid methyl ester	1989	16-Methyl-5a-pregn-16-en-20-one
148	14a,15a-Epoxy-5a-androstane	2030	16-Methylene-17a-hydroxypregnenolone 3-acetate; 3b-Acetoxy-17a-hydroxy-16-methylene-5-pregnen-20-one
643	14a,16a-Dihydroxy-5a-androstane-3,11-dione	2029	16-Methylene-17a-hydroxypregnenolone; 3b,17a-Dihydroxy-16-methylene-5-pregnen-20-one
819	14a,17b-Dihydroxy-4-androsten-3-one	1786	16a,17a-Dihydroxy-4-pregnene-3,20-dione
876	14a-Hydroxy-4-androstene-3,17-dione	1987	16a,17a-Epoxy-3b-hydroxy-16b-methyl-5a-pregnan-20-one
1742	14a-Hydroxy-4-pregnene-3,20-dione; 14a-Hydroxyprogesterone	1698	16a,17a-Epoxy-3b-hydroxy-5-pregnen-20-one
605	14a-Hydroxy-5a-androstane-3,7-dione	1379	16a,17a-Epoxy-3b-hydroxy-5a-pregnan-20-one
986	14b-Methyl-5a-androst-16-en-15-one	810	16a,17a-Epoxy-4-androsten-3-one
927	14b-Methyl-5a-androstane	1329	16a,17a-Epoxy-5a-pregnan-20-one
987	14b-Methyl[d3]-5a-androst-16-en-15-one	1266	16a,17a-Epoxy-5a-pregnane-3b,20b-diol
1906	14b-Pregn-4-ene-3,12,15,20-tetrone	820	16a,17b-Dihydroxy-4-androsten-3-one; 16a-Hydroxytestosterone
273	15,15-Ethylenedithio-5a-androstane	849	16a,17b-Dipropionyloxy-4-androsten-3-one; 16a,17b-Dihydroxy-4-androsten-3-one 16,17-dipropionate
400	15a,16a-Epoxy-5a,13a-androstan-17-one	1232	16a,20a-Diacetoxy-5b-pregnane; 5b-Pregnane-16a,20a-diol diacetate
1037	15a-Hydroxy-3a,5-cyclo-5a-androstane-6,17-dione	1233	16a,20b-Diacetoxy-5a-pregnane; 5a-Pregnane-16a,20b-diol diacetate
877	15a-Hydroxy-4-androstene-3,17-dione	185	16a-Acetoxy-17b-mesyloxy-5a-androstane; 5a-Androstane-16a,17b-diol 16-acetate 17-mesylate
608	15a-Hydroxy-5a-androstane-3,11-dione	354	16a-Acetoxy-5a-androstan-17-one; 16a-Hydroxy-5a-androstan-17-one acetate
611	15a-Hydroxy-5a-androstane-3,12-dione	180	16a-Acetoxy-5a-androstan-17b-ol; 5a-Androstane-16a,17b-diol 16-acetate
619	15a-Hydroxy-5a-androstane-6,11-dione	1317	16a-Acetoxy-5a-pregnan-20-one
401	15b,16b-Epoxy-5a-androstan-17-one		
356	15b-Allyloxy-5a-androstan-17-one; 15b-Hydroxy-5a-androstan-17-one allyl ether		
175	15b-Allyloxy-5a-androstan-17b-ol		
741	15b-Benzoyloxy-5-androstene-3b,17b-diol; 5-Androstene-3b,15b,17b-triol 15-benzoate		
938	15b-Cyano-5a-androstan-17-one		
383	15b-Hydroxy-5a-androstan-17-one ethylenethioacetal		
328	15b-Hydroxy-5a-androstan-17-one; 5a-Androstan-15b-ol-17-one		
342	15b-Methoxy-5a-androstan-17-one		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1223	16a-Acetoxy-5a-pregnan-20b-ol; 5a-Pregnane-16a, 20b-diol 16-acetate	2006	16b-Methyl-16a,17a-epoxypregnenolone; 16a,17a-Epoxy-3b-hydroxy-16b-methyl-5-pregnen-20-one (approx.95%)
1197	16a-Acetoxy-5a-pregnane; 5a-Pregnan-16a-ol acetate	2010	16b-Methyl-4-pregnene-3,20-dione
1222	16a-Acetoxy-5b-pregnan-20a-ol; 5b-Pregnane-16a, 20a-diol 16-acetate	1980	16b-Methyl-5a-pregnan-20-one
1321	16a-Benzoyloxy-5a-pregnan-20-one	1997	16b-Methylpregnenolone acetate; 3b-Acetoxy-16b-methyl-5-pregnen-20-one
294	16a-Bromo-5a,13a-androstan-17-one	1992	16b-Methylpregnenolone; 3b-Hydroxy-16b-methyl-5-pregnen-20-one
295	16a-Bromo-5a-androstan-17-one	1058	17(a)-Oxa-D-homo-5a-androstan-17-one
289	16a-Chloro-5a-androstan-17-one	637	17,17-Difluoro-16b-hydroxy-5a-androstane-3, 11-dione
587	16a-Fluoro-5a-androstane-2,7-dione	239	17,17-Difluoro-1b,11a,16b-triacetoxy-5a-androstane
1743	16a-Hydroxy-4-pregnene-3,20-dione; 16a-Hydroxyprogesterone	106	17,17-Dinitro-5a-androstane
329	16a-Hydroxy-5a-androstan-17-one; 5a-Androstan-16a-ol-17-one	926	17,17-Ethylenedioxy-1,4,6-androstatrien-3-one; 1,4,6-Androstatriene-3,17-dione 17-ethylene ketal
602	16a-Hydroxy-5a-androstane-3,6-dione	652	17,17-Ethylenedioxy-12b,15a-dihydroxy-5a-androstan-3-one
1308	16a-Hydroxy-5b-pregnan-20-one	296	17,17-Ethylenedioxy-16b-bromo-5a-androstane
458	16a-Hydroxyandrostereone; 3a,16a-Dihydroxy-5a-androstan-17-one	1017	17,17-Ethylenedioxy-4,4-dimethyl-5-androsten-3-one; 4,4-Dimethyl-5-androstene-3,17-dione 17-ethylene ketal
300	16a-Iodo-5a-androstan-17-one	895	17,17-Ethylenedioxy-4-androstene-3,11-dione; 4-Androstene-3,11,17-trione 17-ethylene ketal
364	16a-Lauryloxy-5a-androstan-17-one; 16a-Hydroxy-5a-androstan-17-one laurate	760	17,17-Ethylenedioxy-5-androstene
187	16a-Lauryloxy-5a-androstan-17b-ol; 5a-Androstane-16a,17b-diol 16-laurate	761	17,17-Ethylenedioxy-5a-androst-14-ene
1983	16a-Methoxy-21-methyl-5a-pregnan-20-one	762	17,17-Ethylenedioxy-5a-androst-15-ene
929	16a-Methyl-5a-androstan-17b-ol	580	17,17-Ethylenedioxy-5a-androstan-11-one; 5a-Androstane-11,17-dione 17-ethylene ketal
1979	16a-Methyl-5a-pregnan-20-one	578	17,17-Ethylenedioxy-5a-androstan-3-one; 5a-Androstane-3,17-dione 17-ethylene ketal
1976	16a-Methyl-5a-pregnan-20b-ol	546	17,17-Ethylenedioxy-5a-androstan-3a,12b,15a-triol
1996	16a-Methylpregnenolone acetate; 3b-Acetoxy-16a-methyl-5-pregnen-20-one	381	17,17-Ethylenedioxy-5a-androstan-3b-ol; 5a-Androstan-3b-ol-17-one ethylene ketal (crude)
2009	16a-Methylprogesterone; 16a-Methyl-4-pregnene-3,20-dione	271	17,17-Ethylenedioxy-5a-androstane
648	16b,18-Diacetoxy-5a-androstane-3,7-dione	633	17,17-Ethylenedioxy-7a-hydroxy-5a-androstan-3-one; 7a-Hydroxy-5a-androstane-3,17-dione 17-ethylene ketal
1234	16b,20a-Diacetoxy-5a-pregnane; 5a-Pregnane-16b, 20a-diol diacetate	634	17,17-Ethylenedioxy-7b-hydroxy-5a-androstan-3-one; 7b-Hydroxy-5a-androstane-3,17-dione 17-ethylene ketal
363	16b-(p-Nitrobenzoyloxy)-5a-androstan-17-one; 16b-Hydroxy-5a-androstan-17-one p-nitrobenzoate	939	17,17-Ethylenedithio-15b-cyano-5a-androstane
1328	16b-Acetoxy-17a-bromo-5a-pregnan-20-one	1808	17,21-Dihydroxy-4-pregnene-3,20-dione 21-glucuronide (free acid)
2629	16b-Acetoxy-3a,11-dihydroxy-18-nor-4a,8a,14-trimethyl-5a-cholesta-17,24-dien-21-oic acid, sodium salt; Fusidic acid,sodium salt	58	17-Acetoxy-2,3-dimethoxy-1,3,5(10)-estratriene; 2,3-Dimethoxy-1,3,5(10)-estratrien-17-ol acetate
355	16b-Acetoxy-5a-androstan-17-one; 16b-Hydroxy-5a-androstan-17-one acetate	712	17-Acetoxy-5a,13a-androst-16-ene
181	16b-Acetoxy-5a-androstan-17a-ol; 5a-Androstane-16b,17a-diol 16-acetate	46	17-Epietriol; 1,3,5(10)-Estratriene-3,16a,17a-triol
129	16b-Acetoxy-5a-androstane	767	17-Fluoro-5a-androst-16-en-3-one
359	16b-Acetyloxyacetoxy-5a-androstan-17-one	742	17-Fluoro-5a-androst-16-ene-3b,12b,15a-triol
937	16b-Benzyl-5a-androstan-17-one	265	17-Hydroximino-5a-androstane; 5a-Androstan-17-one oxime
393	16b-Bromoandrostereone; 16b-Bromo-3a-hydroxy-5a-androstan-17-one	1181	17-Methyl-3a,12a-diacetoxy-5b-etianic acid methyl ester
290	16b-Chloro-5a-androstan-17-one	970	17-Methylene-5a-androstane
145	16b-Chloro-5a-androstan-17b-ol	1900	17a,18,21-Trihydroxy-4-pregnene-3,11,20-trione; 18-Hydroxycortisone
1004	16b-Cyclohexyl-3b-hydroxy-5-androsten-17-one	1701	17a,20a,21-Trihydroxy-4-pregnen-3-one
130	16b-Ethoxy-5a-androstane	1672	17a,20a-Dihydroxy-4-pregnen-3-one
918	16b-Hydroxy-1,4-androstadiene-3,11-dione	1944	17a,20b,21-Trihydroxy-1,4-pregnadiene-3,11-dione
310	16b-Hydroxy-5a-androstan-3-one; 5a-Androstan-16b-ol-3-one	1673	17a,20b-Dihydroxy-4-pregnen-3-one
609	16b-Hydroxy-5a-androstane-3,11-dione	1330	17a,20b-Dihydroxy-5a-pregnan-3-one
606	16b-Hydroxy-5a-androstane-3,7-dione	1714	17a,21-Diacetoxy-3b-formyloxy-5-pregnen-20-one
1309	16b-Hydroxy-5a-pregnan-20-one	1968	17a,21-Dihydroxy-1,4-pregnadiene-3,11,20-trione 21-glucuronide methyl ester; Prednisone 21-glucuronide methyl ester
1310	16b-Hydroxy-5b-pregnan-20-one	1969	17a,21-Dihydroxy-1,4-pregnadiene-3,11,20-trione 21-glucuronide triacetyl methyl ester; Prednisone 21-G.A.M.E.
86	16b-Hydroxyestrone; 3,16b-Dihydroxy-1,3,5(10)-estratrien-17-one	1967	17a,21-Dihydroxy-1,4-pregnadiene-3,11,20-trione 21-glucuronide; Prednisone 21-glucuronide
2008	16b-Methyl-16a,17a-epoxypregnenolone acetate; 3b-Acetoxy-16a,17a-epoxy-16b-methyl-5-pregnen-20-one (approx.95%)		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1962	17a,21-Dihydroxy-1,4-pregnadiene-3,11,20-trione; Prednisone	1463	17a-Hydroxy-3a,20a,21-triacetoxy-5b-pregnan-11-one
1940	17a,21-Dihydroxy-1,4-pregnadiene-3,20-dione	1744	17a-Hydroxy-4-pregnene-3,20-dione; 17a-Hydroxyprogesterone
1963	17a,21-Dihydroxy-4,6-pregnadiene-3,11,20-trione	302	17a-Hydroxy-5a-androstan-2-one; 5a-Androstan-17a-ol-2-one
1891	17a,21-Dihydroxy-4-pregnene-3,11,20-trione 21-glucuronamide; Cortisone 21-glucuronamide	330	17a-Hydroxy-5a-androstan-3-one 17-sulfate, sodium salt
1890	17a,21-Dihydroxy-4-pregnene-3,11,20-trione 21-glucuronide, ammonium salt; Cortisone 21-glucuronide, ammonium salt	607	17a-Hydroxy-5a-androstane-3,7-dione
1889	17a,21-Dihydroxy-4-pregnene-3,11,20-trione 21-glucuronide; Cortisone 21-glucuronide	1311	17a-Hydroxy-5a-pregnan-20-one
1811	17a,21-Dihydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester	1600	17a-Hydroxy-5a-pregnene-3,11,20-trione
1814	17a,21-Dihydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester 3,20-disemicarbazone	1499	17a-Hydroxy-5a-pregnene-3,20-dione
2052	17a,21-Dihydroxy-6a-methyl-1,4-pregnadiene-3,11,20-trione; 6a-Methylprednisone	331	17a-Hydroxy-5b-androstan-3-one 17-sulfate, pyridinium salt
2059	17a-Acetoxy-3a,5-cyclo-5a-pregnan-20-one	311	17a-Hydroxy-5b-androstan-3-one; 5b-Androstan-17a-ol-3-one
1688	17a-Acetoxy-3b-hydroxy-5-pregnen-20-one; 17a-Hydroxypregnenolone 17-acetate	1297	17a-Hydroxy-5b-pregnan-3-one
782	17a-Acetoxy-4-androsten-3-one; 4-Androsten-17a-ol-3-one acetate	2015	17a-Hydroxy-6a-methyl-4-pregnene-3,20-dione
386	17a-Acetoxy-5a-androstan-2-one ethylene ketal; 17a-Acetoxy-2,2-ethylenedioxy-5a-androstane	2017	17a-Hydroxy-6a-methylprogesterone 17-butyrate
346	17a-Acetoxy-5a-androstan-3-one; 5a-Androstan-17a-ol-3-one acetate	2018	17a-Hydroxy-6a-methylprogesterone 17-pivalate
2054	17a-Acetoxy-6-chloro-16-methylene-4,6-pregnadiene-3,20-dione	1937	17a-Hydroxypregna-4,14-diene-3,20-dione
2055	17a-Acetoxy-6-methoxy-16-methylene-4,6-pregnadiene-3,20-dione	1764	17a-Hydroxyprogesterone 3-(O-carboxymethyl)oxime (approx.97%)
2033	17a-Acetoxy-6-methyl-4,6-pregnadiene-3,20-dione; Megestrol acetate	1773	17a-Hydroxyprogesterone 3-(O-carboxymethyl) oxime:BSA conjugate
2032	17a-Acetoxy-6a-methyl-1,4-pregnadiene-3,20-dione; 1-Dehydroprovera	339	17a-Methoxy-5a-androstan-3-one
1003	17a-Allyl-17b-hydroxy-4-androsten-3-one; 17 a-Allyltestosterone	1020	17a-Methyl-1,5-androstadiene-3b,17b-diol
103	17a-Allyl-4-estren-17b-ol	2073	17a-Methyl-19-norpregna-4,9(10)-diene-3,20-dione; Demegestone
794	17a-Benzoyloxy-4-androsten-3-one; 4-Androsten-17a-ol-3-one benzoate	2011	17a-Methyl-4-pregnene-3,20-dione (crude)
1327	17a-Bromo-16b-hydroxy-5a-pregnan-20-one	932	17a-Methyl-5b-androstane-3a,17b-diol
2051	17a-Butyryloxy-21-chloro-9a-fluoro-16b-methyl-1,4-pregnadiene-3,11,20-trione	978	17a-Methylandrost-5-ene-3b,17b-diol
771	17a-Epitestosterone; 17a-Hydroxy-4-androsten-3-one	933	17a-Methylandrost-5-ene-3b,17b-diol dipropionate
2060	17a-Ethylestr-4-en-17b-ol; Ethylestrenol; 19-Nor-4-pregnen-17b-ol	949	17a-Methylandrost-17b-ol-3-one 3-(O-carboxymethyl)oxime
1656	17a-Ethyltestosterone; 17b-Hydroxy-4-pregnen-3-one; 17a-Ethyl-17b-hydroxy-4-androsten-3-one	1009	17a-Methyltestosterone 3-(O-carboxymethyl)oxime
2079	17a-Ethynyl-17b-hydroxy-18-methyl-4-estren-3-one; 13-Ethyl-17b-hydroxy-18,19-bisnorpregn-4-en-20-yn-3-one; Norgestrel	1908	17a-Pregna-5,20-dien-3b,17b-diol
2086	17a-Ethynyl-17b-hydroxy-4,9,11-estratrien-3-one; 17b-Hydroxy-19-norpregna-4,9,11-trien-20-yn-3-one	1757	17a-Propionyloxy-4-pregnene-3,20-dione; 17a-Hydroxyprogesterone 17-propionate
2082	17a-Ethynyl-3-methoxy-1,3,5(10)-estratrien-17b-ol; 3-Methoxy-19-nor-17a-pregna-1,3,5(10)-trien-20-yn-17b-ol; Mestranol	979	17a-Propyl-5-androstene-3b,17-diol
2074	17a-Ethynyl-4-estren-17b-ol; 19-Norpregn-4-en-20-yn-17b-ol; Lynestrenol	1702	17b,20a,21-Trihydroxy-17a-pregna-4-en-3-one
1917	17a-Ethynyl-5a-androstane-3b,17b-diol diacetate; 3b,17b-Diacetoxy-5a-pregna-20-yne	796	17b-(p-Bromobenzoyloxy)-4-androsten-3-one; Testosterone p-bromobenzoate
2083	17a-Ethynylestradiol 17b-(a-glucuronide), potassium salt	369	17b-(p-Toluenesulfonyloxy)-5a-androstan-15-one; 5a-Androstan-17b-ol-15-one tosylate
2081	17a-Ethynylestradiol; 17a-Ethynyl-1,3,5(10)-estratriene-3,17b-diol; 19-Nor-1,3,5(10)-pregnatrien-20-yne-3,17b-diol	908	17b-Acetoxy-1,4-androstadien-3-one
275	17a-Fluoro-5a-androstan-3-one (crude)	144	17b-Acetoxy-12,12-difluoro-5a-androstane
589	17a-Fluoro-5a-androstane-3,11-dione	975	17b-Acetoxy-16-benzylidene-5a-androstane
2067	17a-Hydroxy-19-norprogesterone; 17a-Hydroxy-19-nor-4-pregnene-3,20-dione	2080	17b-Acetoxy-17a-ethynyl-19-nor-4-androsten-3-one; 17b-Acetoxy-19-norpregn-4-en-20-yn-3-one; Norethisterone acetate
		1913	17b-Acetoxy-17a-pregna-5,20-dien-3b-ol; 17a-Pregna-5,20-diene-3b,17b-diol 17-acetate
		387	17b-Acetoxy-2,2-ethylenedioxy-5a-androstane; 17b-Hydroxy-5a-androstan-2-one 17-acetate ethylene ketal
		57	17b-Acetoxy-2-methoxy-1,3,5(10)-estratrien-3-ol; 2-Methoxyestradiol 17-acetate
		1008	17b-Acetoxy-2-methyl-5a-androst-1-en-3-one; Stenobolone acetate
		19	17b-Acetoxy-2a-hydroxy-4-estren-3-one
		909	17b-Acetoxy-3,5-androstadien-7-one
		14	17b-Acetoxy-4-estren-3-one
		945	17b-Acetoxy-4a-methyl-5a-androstan-3-one
		711	17b-Acetoxy-5-androstene
		784	17b-Acetoxy-5a-androst-9(11)-en-12-one

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
352	17b-Acetoxy-5a-androstan-16-one; 5a-Androstan-17-b-ol-16-one acetate	1001	17b-Hydroxy-17a-propyl-4-androsten-3-one; 17a-Propyltestosterone
344	17b-Acetoxy-5a-androstan-2-one; 5a-Androstan-17-b-ol-2-one acetate	104	17b-Hydroxy-2,2,17a-trimethyl-4,9,11-estratrien-3-one
347	17b-Acetoxy-5a-androstan-3-one; 5a-Androstan-17-b-ol-3-one acetate	1013	17b-Hydroxy-2-hydroxymethylene-17a-methyl-5a-androstan-3-one
348	17b-Acetoxy-5b-androstan-3-one; 5b-Androstan-17-b-ol-3-one acetate	1055	17b-Hydroxy-3,5-seco-A-norandrostan-5-one-3-oic acid
1025	17b-Acetoxy-6-methylene-4-androsten-3-one; 6-Methylenetestosterone acetate	2019	17b-Hydroxy-3-oxo-17a-pregn-4-ene-21-carboxylic acid $\gamma$ -lactone
143	17b-Acetoxy-7,7-difluoro-5a-androstane	906	17b-Hydroxy-4,6-androstadien-3-one
464	17b-Acetoxy-9a-Hydroxy-5a-androstan-3-one; 9a,17b-Dihydroxy-5a-androstan-3-one 17-acetate	798	17b-Hydroxy-4-androsten-3-one 17-glucuronide, potassium salt; Testosterone b-D-glucuronide, potassium salt
1073	17b-Acetoxy-des-A-androst-9-en-5-one	773	17b-Hydroxy-4-androsten-3-one-[19-d3]; Testosterone-[19-d3]
395	17b-Benzoyloxy-2a,4a-dibromo-5a-androstan-3-one	13	17b-Hydroxy-4-estren-3-one; 19-Nortestosterone
912	17b-Benzoyloxy-4,6-androstadien-3-one	320	17b-Hydroxy-5a-androstan-16-one; 5a-Androstan-17-b-ol-16-one
795	17b-Benzoyloxy-4-androsten-3-one; 4-Androsten-17-b-ol-3-one benzoate; Testosterone benzoate	303	17b-Hydroxy-5a-androstan-2-one; 5a-Androstan-17-b-ol-2-one
361	17b-Benzoyloxy-5a-androstan-3-one; 5a-Androstan-17-b-ol-3-one benzoate	332	17b-Hydroxy-5a-androstan-3-one 17-sulfate, pyridinium salt
1067	17b-Benzoyloxy-6a-bromo-des-A-10b(H)-androstan-5-one	599	17b-Hydroxy-5a-androstane-2,7-dione
1074	17b-Benzoyloxy-des-A-androst-9-en-5-one	610	17b-Hydroxy-5a-androstane-3,11-dione
150	17b-Chloro-16a,17a-epoxy-5a-androstane	615	17b-Hydroxy-5a-androstane-3,16-dione
38	17b-Estradiol 17-(b-D-glucuronide), sodium salt	313	17b-Hydroxy-5b-androstan-3-one; 5b-Androstan-17b-ol-3-one
30	17b-Estradiol 17-acetate; 17b-Acetoxy-1,3,5(10)-estratrien-3-ol	4	17b-Hydroxy-5b-estran-3-one; 17b-Hydroxy-19-nor-5b-androstan-3-one
41	17b-Estradiol 17-cypionate; 17b-(Cyclopentanepropionyloxy)-1,3,5(10)-estratrien-3-ol	1021	17b-Hydroxy-6-methylene-4-androsten-3-one; 6-Methylenetestosterone
39	17b-Estradiol 17-glucuronide 3-sulfate, dipotassium salt + approx. 30% Trizma base	774	17b-Hydroxy-9b,10a-androst-4-en-3-one; Retrotestosterone
34	17b-Estradiol 17-hemisuccinate	1072	17b-Hydroxy-des-A-androst-9-en-5-one
32	17b-Estradiol 17-propionate; 17b-Propionyloxy-1,3,5(10)-estratrien-3-ol	997	17b-Hydroxymethyl-4-androsten-3-one
35	17b-Estradiol 17-valerate; 17b-Pentanoyloxy-1,3,5(10)-estratrien-3-ol	977	17b-Hydroxymethyl-5-androsten-3b-ol
36	17b-Estradiol 3-(b-D-glucuronide), sodium salt	297	17b-Iodo-5a-androstan-16-one
37	17b-Estradiol 3-glucuronide 17-sulfate, dipotassium salt	244	17b-Methoxy-3b,12b,15a-triacetoxy-5a-androstane
29	17b-Estradiol 3-methyl ether; 3-Methoxy-1,3,5(10)-estratrien-17b-ol	985	17b-Methyl-4-androsten-3-one
28	17b-Estradiol 3-sulfate, sodium salt + approx. 67% N-methyl-D-glucamine	105	17b-Nitro-5a-androstane
31	17b-Estradiol diacetate; 3,17b-Diacetoxy-1,3,5(10)-estratriene	17	17b-Phenylpropionyloxy-4-estren-3-one
33	17b-Estradiol dipropionate; 3,17b-Dipropionyloxy-1,3,5(10)-estratriene	910	17b-Propionyloxy-1,4-androstadien-3-one
26	17b-Estradiol; 1,3,5(10)-Estratriene-3,17b-diol	357	17b-Propionyloxy-5a-androstan-3-one; 5a-Androstan-17-b-ol-3-one propionate
131	17b-Ethoxy-5a-androstane	1697	18,20b-Diacetoxy-5-pregnen-3-one ethylene ketal
276	17b-Fluoro-5a-androstan-3-one	1478	18,21-Diacetoxy-11b,18-Epoxy-3a-hydroxy-5b-pregnan-20-one 3-glucuronide triacetyl methyl ester
905	17b-Hydroxy-1,4-androstadien-3-one	1887	18,21-Diacetoxy-4-pregnene-3,20-dione
1012	17b-Hydroxy-17a-hydroxymethyl-4-androsten-3-one; 17a-Hydroxymethyltestosterone	1879	18,21-Dihydroxy-4-pregnene-3,11,20-trione
976	17b-Hydroxy-17a-methyl-2'H-5a-androst-2-eno[3,2-c]pyrazole	1788	18,21-Dihydroxy-4-pregnene-3,20-dione; 18-Hydroxydeoxycorticosterone
943	17b-Hydroxy-17a-methyl-2-oxa-5a-androstan-3-one; Oxandrolone	1941	18,21-Dihydroxypregn-1,4-diene-3,20-dione 18->20-hemiketal
996	17b-Hydroxy-17a-methyl-4-androsten-3-one; 17a-Methyltestosterone	1163	18-Hydroxy-3-oxo-eti-4-en-20-oic acid 20,18-lactone
1018	17b-Hydroxy-17a-methyl-4-androstene-3,11-dione	1745	18-Hydroxy-4-pregnene-3,20-dione
941	17b-Hydroxy-17a-methyl-5a-androstan-3-one	1837	19,21-Diacetoxy-14-hydroxy-14b-pregn-4-ene-3,20-dione
960	17b-Hydroxy-17a-methyl-5a-androstane-3,11-dione	1804	19,21-Diacetoxy-4-pregnene-3,20-dione
942	17b-Hydroxy-17a-methyl-5b-androstan-3-one	1790	19,21-Dihydroxy-14b,17a-pregn-4-ene-3,20-dione
961	17b-Hydroxy-17a-methyl-5b-androstane-3,11-dione	1789	19,21-Dihydroxy-4-pregnene-3,20-dione
931	17b-Hydroxy-17a-methylandrostan-3,2-c]pyrazole; Stanozolol	1799	19-Acetoxy-14-hydroxy-14b-pregn-4-ene-3,20-dione
		1801	19-Acetoxy-21-hydroxy-4-pregnene-3,20-dione
		1180	19-Acetoxy-3,5,14-etiatrienic acid ethyl ester
		1166	19-Acetoxy-3-oxo-14b,17a-eti-4-enic acid
		1165	19-Acetoxy-3-oxoeti-4-enic acid
		1753	19-Acetoxy-4-pregnene-3,20-dione
		1746	19-Hydroxy-14b,17a-pregn-4-ene-3,20-dione
		1179	19-Hydroxy-3,5,14-etiatrienic acid ethyl ester

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1162	19-Hydroxy-3-oxo-14b,17a-eti-4-enic acid	83	2,3-Dihydroxy-1,3,5(10)-estratrien-17-one;
1164	19-Hydroxy-3-oxo-14b,17a-eti-4-enic acid ethyl ester		2-Hydroxyestrone
879	19-Hydroxy-4-androstene-3,17-dione	56	2,3-Dimethoxy-1,3,5(10)-estratrien-17b-ol
2064	19-Nor-14b,17a-pregn-4-ene-3,20-dione	2642	2,3-Seco-2,3-dicyano-5a-cholestane
2063	19-Nor-4-pregnene-3,20-dione	1054	2,3-Seco-4,4-dimethyl-5a-androstane-2,3-dioic acid
2936	19-Nor-ergosta-3,5,7,9-tetraene	1053	2,3-Seco-5a-androstane-2,3-dioic acid 3-amide
2929	19-Nor-ergosta-5,7,9-triene	1052	2,3-Seco-5a-androstane-2,3-dioic anhydride
2937	19-Norergosta-5,7,9,22-tetraen-3b-ol; Neoergosterol	2649	2,3-Seco-5a-cholestane-2,3-dioic acid
2931	19-Norergosta-5,7,9-trien-3-one	2655	2,3-Seco-5a-cholestane-2,3-dioic acid 2-methyl ester 3-amide
2930	19-Norergosta-5,7,9-triene-3b,4x-diol		
15	19-Nortestosterone 17-propionate; 17b-Propionyloxy-4-estren-3-one	2652	2,3-Seco-5a-cholestane-2,3-dioic acid anhydride
1051	19-Nortestosterone 3-(O-carboxymethyl)oxime	2653	2,3-Seco-5a-cholestane-2,3-dioic acid dimethyl ester
16	19-Nortestosterone cypionate; 17b-Cyclopentanepropionyloxy-4-estren-3-one	2644	2,3-Seco-5a-cholestane-2,3-diol
1519	1a,12b-Dihydroxy-5a-pregnane-7,20-dione	2650	2,3-Seco-5b-cholestane-2,3-dioic acid
2550	1a,24R,25-Trihydroxycholesterol; (24R)-Cholest-5-ene-1a,3b,24,25-tetrol	2943	2,3-Seco-5b-lumist-7-en-2,3-dioic acid
2551	1a,25S,26-Trihydroxycholesterol; (25S)-Cholest-5-ene-1a,3b,25,26-tetrol	2944	2,3-Seco-5b-lumista-7,22-diene-2,3-dioic acid
2980	1a,2a-Diacetoxylupane; Lupane-1a,2a-diol diacetate	2939	2,3-Seco-5b-lumistane-2,3-dioic acid
396	1a,2a-Epoxy-5a-androstan-3-one	2986	2,3-Seco-lupane-2,3-dioic acid anhydride
243	1a,3a,17b-Triacetoxo-5a-androstan-9a-ol; 5a-Androstane-1a,3a,9a,17b-tetrol 1,3,17-triacetate	2461	2,4-Dihydroximino-5a-cholestan-3-one; 5a-Cholestane-2,3,4-trione 2,4-di(oxime)
479	1a,3b-Diacetoxo-5a-androstan-17-one	2643	2,4-Seco-2,4-dinitro-5a-cholestane
431	1a,6a-Dihydroxy-5a-androstan-11-one	2942	2,4-Seco-A-nor-5a-ergost-7-ene-2,4-dioic acid
2526	1a-Benzoyloxy-5a-cholest-2-ene; 5a-Cholest-2-en-1a-ol benzoate	2918	2-(Hydroxymethylene)-5a-ergost-7-en-3-one
2394	1a-Iodo-2b,3b-di-(p-chlorobenzoyloxy)-5a-cholestane; 1a-Iodo-5a-cholestane-2b,3b-diol bis(4-chlorobenzoate)	2832	2-(Hydroxymethylene)-5b-ergosta-7,22-dien-3-one
542	1b,3a,15a-Triacetoxo-5a-androstan-17-one(crude)	1006	2-Anisylidene-17b-hydroxy-5a-androstan-3-one
541	1b,3a,6a-Triacetoxo-5a-androstan-17-one	2014	2-Anisylidene-20,20-ethylenedioxy-5a-pregnan-3-one; 2-Anisylidene-5a-pregnane-3,20-dione 20-ethylene ketal
1487	1b,3b,11a,19-Tetraacetoxo-5,14,21-trihydroxy-5b,14b-pregnan-20-one	995	2-Anisylidene-4,4-dimethyl-5a-androstan-3-one
1715	1b,3b,11a-Triacetoxo-5-pregnen-20-one	991	2-Anisylidene-5a-androstan-3-one
1705	1b,3b,11a-Trihydroxy-5-pregnen-20-one	2624	2-Anisylidene-5a-cholestan-3b-ol
1384	1b,3b,11a-Trihydroxy-5a-pregnan-20-one	2991	2-Anisylidenelup-20(29)-en-3-one
534	1b,3b,17b-Triacetoxo-5a-androstan-6-one	990	2-Benzylidene-5a-androstan-3-one
2110	1b,3b,5,11a,14,19-Hexahydroxy-5b,14b-card-20(22)-enolide 1,3-acetonide; Oubagenin 1,3-acetonide	2626	2-Benzylidene-5a-cholestan-3-one
1694	1b,3b-Diacetoxo-5-pregnen-20-one	2923	2-Benzylidene-5b-lumista-7,22-dien-3-one
446	1b,3b-Dihydroxy-5a-androstan-17-one	2990	2-Benzylidenelup-20(29)-en-3-one
497	1b,6a-Diacetoxo-3,3-difluoro-5a-androstan-17-one	984	2-Cyano-4,17a-dimethyl-4a,5-epoxy-5a-androst-2-ene-3,17b-diol
525	1b,6a-Dihydroxy-3a-methoxy-5a-androstan-17-one	983	2-Cyano-4a,5-epoxy-5a-androst-2-ene-3,17b-diol
447	1b,6a-Dihydroxy-5a-androstan-17-one	585	2-Hydroximino-5a-androstan-3-one; 5a-Androstane-2,3-dione 2-oxime
538	1b,7b,17b-Triacetoxo-5a-androstan-11-one	2453	2-Hydroximino-5a-cholestan-3-one; 5a-Cholestane-2,3-dione 2-oxime
513	1b,7b,17b-Trihydroxy-5a-androstan-11-one	2985	2-Hydroximinolupane-3-one; Lupane-2,3-dione 2-oxime
850	1b,7b-Dihydroxy-3b-fluoro-5-androsten-17-one	919	2-Hydroxy-1,4-androstadiene-3,17-dione
620	1b-Hydroxy-5b-androstane-6,17-dione	1000	2-Hydroxymethylene-4,4-dimethyl-5a-androstan-3-one
2056	2'-(4-Fluorophenyl)-2'H-pregna-2,4-dien-20-yno[3,2-c]pyrazol-17a-ol; Nivazol	998	2-Hydroxymethylene-5a-androstan-3-one
1026	2,16-Dianisylidene-5a-androstane-3,17-dione	2924	2-Hydroxymethylene-5b-lumista-7,22-dien-3-one
2933	2,2,4,4-Tetramethyl-19-norergosta-5,7,9-trien-3-one	72	2-Methoxy-1,3,5(10)-estratriene-3,16a,17b-triol; 2-Methoxyestriol
2625	2,2,4,4-Tetramethyl-5-cholesten-3-one	52	2-Methoxyestradiol; 2-Methoxy-1,3,5(10)-estratriene-3,17b-diol
2415	2,2-Dibromo-5a-cholestan-3-one	2992	2-Veratrylidenelup-20(29)-en-3-one
572	2,2-Ethylenedioxy-5a-androstan-17-one	1439	20,20-Ethylenedioxy-3a,11b,21-trihydroxy-5a-pregnane; 3a,11b,21-Trihydroxy-5a-pregnan-20-one ethylene ketal
379	2,2-Ethylenedioxy-5a-androstan-17a-ol; 17a-Hydroxy-5a-androstan-2-one, ethylene ketal	1494	20,20-Ethylenedioxy-5a-pregnan-11-one; 5a-Pregnane-11,20-dione 20-ethylene ketal (crude)
3009	2,3-Bis(3'-hydroxybenzyl)butyrolactone, trans-; Enterolactone	1324	20,20-Ethylenedioxy-5a-pregnan-11b-ol; 11b-Hydroxy-5a-pregnan-20-one ethylene ketal (crude)
2454	2,3-Dihydroximino-5a-cholestane; 5a-Cholestane-2,3-dione di(oxime)	1492	20,20-Ethylenedioxy-5a-pregnan-3-one
		1323	20,20-Ethylenedioxy-5a-pregnan-3b-ol
		1373	20,20-Ethylenedioxy-5b-pregnane-3a,6a-diol
		1634	20-Bromo-5a-pregn-17(20)-ene-3b,21-diol

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2024	20-Cyano-3b,21-diacetoxypregna-5,17(20)-diene; 20-Cyanopregna-5,17(20)-diene-3b,21-diol diacetate	2043	21-Acetoxy-11b,17a-dihydroxy-9a-fluoro-16a-methyl-1,4-pregnadiene-3,20-dione; Dexamethasone 21-acetate
1440	20a,21-(Isopropylidenedioxy)-3a-hydroxy-5b-pregnan-11-one	2044	21-Acetoxy-11b,17a-dihydroxy-9a-fluoro-16b-methyl-1,4-pregnadiene-3,20-dione; Betamethasone 21-acetate
1712	20a,21-Diacetoxy-17b-hydroxy-17a-pregn-4-en-3-one	1885	21-Acetoxy-11b-hydroxy-4-pregnen-18-al-3,20-dione; Aldosterone 21-acetate
1427	20a,21-Diacetoxy-3a-hydroxy-5b-pregnan-11-one 3-glucuronide triacetyl methyl ester	1812	21-Acetoxy-11b-hydroxy-4-pregnene-3,20-dione 11-glucuronide triacetyl methyl ester
1260	20a,21-Diacetoxy-5b-pregnan-17a-ol; 5b-Pregnane-17a,20a,21-triol 20,21-diacetate	1798	21-Acetoxy-11b-hydroxy-4-pregnene-3,20-dione; Corticosterone 21-acetate
1289	20a,21-Diacetoxy-5b-pregnane-3a,11b,17-triol 3-glucuronide triacetyl methyl ester	1549	21-Acetoxy-14a,15a-epoxy-5b-pregnane-3,20-dione
1235	20a,21-Diacetoxy-5b-pregnane; 5b-Pregnane-20a,21-diol diacetate	1832	21-Acetoxy-14a,17a-dihydroxy-4-pregnene-3,20-dione
1674	20a,21-Dihydroxy-4-pregnen-3-one	1719	21-Acetoxy-16,17-epoxypregnenolone; 21-Acetoxy-16,17-epoxy-3-hydroxy-5-pregnen-20-one
1977	20a-Acetoxy-16b-methyl-5a-pregnane; 16b-Methyl-5a-pregnan-20a-ol acetate	1965	21-Acetoxy-17a-hydroxy-1,4-pregnadiene-3,11,20-trione; Prednisone 21-acetate
1199	20a-Acetoxy-5a,17a-pregnane; 5a,17a-Pregnan-20a-ol acetate	1966	21-Acetoxy-17a-hydroxy-4,6-pregnadiene-3,11,20-trione
1224	20a-Acetoxy-5a-pregnan-17a-ol; 5a-Pregnane-17a,20a-diol 20-acetate	1800	21-Acetoxy-17a-hydroxy-4-pregnene-3,20-dione
1198	20a-Acetoxy-5a-pregnane; 5a-Pregnan-20a-ol acetate	1610	21-Acetoxy-17a-hydroxy-5a-pregnane-3,11,20-trione
1225	20a-Acetoxy-5b-pregnan-17a-ol; 5b-Pregnane-17a,20a-diol 20-acetate	1611	21-Acetoxy-17a-hydroxy-5b-pregnane-3,11,20-trione
1221	20a-Acetoxy-5b-pregnan-3a-ol; 5b-Pregnane-3a,20a-diol 20-acetate	1833	21-Acetoxy-19,21-dihydroxy-14b,17a-pregn-4-ene-3,20-dione
1202	20a-Benzoyloxy-5a-pregnane; 5a-Pregnan-20a-ol benzoate	1802	21-Acetoxy-19-hydroxy-4-pregnene-3,20-dione
1922	20a-Hydroxy-1,4-pregnadien-3-one	1470	21-Acetoxy-3a,11b,17a-trihydroxy-5b-pregnan-20-one 3-glucuronide triacetyl methyl ester
1657	20a-Hydroxy-4-pregnen-3-one	1406	21-Acetoxy-3a,12a-dihydroxy-5b-pregnan-20-one
1441	20b,21-(Isopropylidenedioxy)-3a-hydroxy-5b-pregnan-11-one	1471	21-Acetoxy-3a,17a,20b-trihydroxy-5b-pregnan-11-one 3-glucuronide triacetyl methyl ester
1834	20b,21-Diacetoxy-17a-hydroxy-4-pregnene-3,11-dione	1407	21-Acetoxy-3a,17a-dihydroxy-5b-pregnan-20-one
1428	20b,21-Diacetoxy-3a-hydroxy-5b-pregnan-11-one 3-glucuronide triacetyl methyl ester	1436	21-Acetoxy-3a,17a-dihydroxy-5b-pregnan-20-one 3-glucuronide triacetyl methyl ester
1261	20b,21-Diacetoxy-5a-pregnan-17a-ol; 5a-Pregnane-17a,20b,21-triol 20,21-diacetate	1563	21-Acetoxy-3a,17a-dihydroxy-5b-pregnane-11,20-dione
1236	20b,21-Diacetoxy-5b-pregnane; 5b-Pregnane-20b,21-diol diacetate	1575	21-Acetoxy-3a,17a-dihydroxy-5b-pregnane-11,20-dione 3-glucuronide triacetyl methyl ester
1675	20b,21-Dihydroxy-4-pregnen-3-one	1566	21-Acetoxy-3a,17a-dihydroxy-5b-pregnane-11,20-dione 3-sulfate, triethylammonium salt
1925	20b-Acetoxy-1,4-pregnadien-3-one	1528	21-Acetoxy-3a-hydroxy-5a-pregnane-11,20-dione; Alphadolane acetate
2089	20b-Acetoxy-5-oxo-3,5-seco-A-norpregnan-3-oic acid	1359	21-Acetoxy-3a-hydroxy-5b-pregnan-20-one
2090	20b-Acetoxy-5-oxo-3,5-seco-A-norpregnan-3-oic acid, sodium salt	1539	21-Acetoxy-3a-hydroxy-5b-pregnane-11,20-dione 3-glucuronide triacetyl methyl ester
1200	20b-Acetoxy-5a,17a-pregnane; 5a,17a-Pregnan-20b-ol acetate	1709	21-Acetoxy-3b,11a-dihydroxy-5-pregnen-20-one
1666	20b-Benzoyloxy-4-pregnen-3-one	1476	21-Acetoxy-3b,17a-dihydroxy-5,6a-epoxy-5a-pregnan-20-one
1203	20b-Benzoyloxy-5a-pregnane; 5a-Pregnan-20b-ol benzoate	1564	21-Acetoxy-3b,17a-dihydroxy-5a-pregnane-11,20-dione
1923	20b-Hydroxy-1,4-pregnadien-3-one	1565	21-Acetoxy-3b,17a-dihydroxy-5b-pregnane-11,20-dione
2069	20b-Hydroxy-18-methyl-19-nor-5-pregnene-3,18-dione 3-ethylene ketal	1485	21-Acetoxy-3b,5-dihydroxy-8b,19-epoxy-19-methoxy-5b,14b,17a-pregnan-20-one
1658	20b-Hydroxy-4-pregnen-3-one	1690	21-Acetoxy-3b-hydroxy-5-pregnen-20-one
2087	20b-Hydroxy-5-oxo-3,5-seco-A-norpregnan-3-oic acid	1360	21-Acetoxy-3b-hydroxy-5a-pregnan-20-one
2088	20b-Hydroxy-5-oxo-3,5-seco-A-norpregnan-3-oic acid methyl ester	1929	21-Acetoxy-3b-hydroxypregna-5,16-dien-20-one
1300	20b-Hydroxy-5a-pregnan-3-one	1872	21-Acetoxy-4-pregnene-3,11,20-trione
2076	20S-(Dimethylaminoethyl)amino-19-norpregna-1,3,5(10)-trien-3-ol	1873	21-Acetoxy-4-pregnene-3,18,20-trione
2075	20S-Amino-19-norpregna-1,3,5(10)-trien-3-ol	1874	21-Acetoxy-4-pregnene-3,19,20-trione
2046	21-(Adamantane-1-carboxy)-11b,17a-dihydroxy-9a-fluoro-16b-methyl-1,4-pregnadiene-3,20-dione; Betamethasone 21-adamantanoate	1446	21-Acetoxy-5,6a-epoxy-3b-hydroxy-5a-pregnan-20-one
1526	21-Acetoxy-11a-hydroxy-5a-pregnane-3,20-dione	1612	21-Acetoxy-5a,6b,17a-trihydroxy-5a-pregnane-3,11,20-trione
1723	21-Acetoxy-11b,17a,20a-trihydroxy-4-pregnen-3-one	1318	21-Acetoxy-5a-pregnan-20-one
1932	21-Acetoxy-11b,17a,20b-trihydroxy-1,4-pregnadien-3-one	1226	21-Acetoxy-5a-pregnan-20b-ol; 5a-Pregnane-20b,21-diol 21-acetate
1831	21-Acetoxy-11b,17a-dihydroxy-4-pregnene-3,20-dione; Cortisol 21-acetate	1253	21-Acetoxy-5a-pregnane-17a,20b-diol; 5a-Pregnane-17a,20b,21-triol 21-acetate
1562	21-Acetoxy-11b,17a-dihydroxy-5b-pregnane-3,20-dione	1509	21-Acetoxy-5a-pregnane-3,20-dione
		1287	21-Acetoxy-5a-pregnane-3b,11b,17a,20b-tetrol; 5a-Pregnane-3b,11b,17a,20b,21-pentol 21-acetate

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1254	21-Acetoxy-5b-pregnane-17a,20b-diol; 5b-Pregnane-17a,20b,21-triol 21-acetate	2583	22R-Acetoxy-2b,3b,14a,20b,25-pentahydroxy-7-cholesten-6-one; b-Ecdysterone 22-acetate
1604	21-Acetoxy-5b-pregnane-3,11,20-trione	2134	24-Hydroxy-4-cholen-3-one
1605	21-Acetoxy-5b-pregnane-3,12,20-trione	476	2a,12b-Diacetoxy-5a-androstan-15-one
1862	21-Acetoxy-6a,11b,17a-trihydroxy-4-pregnene-3,20-dione	2384	2a,3-(Bis-3',5'-dinitrobenzoyloxy)-5a-cholestane; 5a-Cholestane-2a,3-diol bis(3,5-dinitrobenzoate)
1960	21-Acetoxy-6a,9a-Difluoro-11b,16a,17a-trihydroxy-1,4-pregnadiene-3,20-dione 16,17-acetonide; Fluocinolide	2395	2a,3a-Epoxy-5a-cholestan-5-ol
2047	21-Acetoxy-6a,9a-difluoro-11b,17a-dihydroxy-16a-methyl-1,4-pregnadiene-3,20-dione; Flumethasone acetate	2376	2a,3b-Diacetoxy-5a-cholestane; 5a-Cholestane-2a,3b-diol diacetate
2037	21-Acetoxy-6a-methyl-11b,17a-dihydroxy-1,4-pregnadiene-3,20-dione; 6a-Methylprednisolone acetate	1942	2a,3b-Diacetoxypregna-5,16-diene-12,20-dione
1863	21-Acetoxy-6b,11b,17a-trihydroxy-4-pregnene-3,20-dione	597	2a,4a-Dibromo-5a-androstan-3,17-dione
1901	21-Acetoxy-6b,17a-dihydroxy-4-pregnene-3,11,20-trione	2416	2a,4a-Dibromo-5a-cholestan-3-one
1797	21-Acetoxy-6b-hydroxy-4-pregnene-3,20-dione	228	2a,6a,16a-Triacetoxy-5a-androstan-3-one; 5a-Androstan-2a,6a,16a-triol triacetate
1855	21-Acetoxy-8b,19-epoxy-19-hydroxy-14b,17a-pregn-4-ene-3,20-dione	448	2a,7a(+7b)-Dihydroxy-5a-androstan-17-one (7a:7b=3:1)
1854	21-Acetoxy-8b,19-epoxy-19-hydroxy-4-pregnene-3,20-dione	442	2a,7a-Dihydroxy-5a-androstan-16-one
1857	21-Acetoxy-8b,19-epoxy-19-methoxy-14b,17a-pregn-4-ene-3,20-dione	948	2a-Acetoxy-4,4-dimethyl-5a-androstan-3-one; 4,4-Dimethyl-5a-androstan-2a-ol-3-one acetate
1856	21-Acetoxy-8b,19-epoxy-19-methoxy-4-pregnene-3,20-dione	2572	2a-Acetoxy-4-cholesten-3-one; 2a-Hydroxy-4-cholesten-3-one acetate
1847	21-Acetoxy-9a-Fluoro-11b,17a-dihydroxy-4-pregnene-3,20-dione	345	2a-Acetoxy-5a-androstan-3-one; 5a-Androstan-2a-ol-3-one acetate
2002	21-Benzylidene-3b-acetoxy-5a-pregnan-20-one	177	2a-Acetoxy-5a-androstan-3b-ol
2001	21-Benzylidene-3b-hydroxy-5b-pregnan-20-one	122	2a-Acetoxy-5a-androstan-3-one
1544	21-Bromo-3a,17a-dihydroxy-5b-pregnane-11,20-dione	2983	2a-Acetoxy-5a-androstan-3-one
1326	21-Chloro-3b-hydroxy-5a-pregnan-20-one	2417	2a-Bromo-2b-nitro-5a-cholestan-3-one
2035	21-Chloro-9a-fluoro-11b-hydroxy-16b-methyl-17a-propionyloxy-1,4-pregnadiene-3,20-dione; Clobetasol propionate	291	2a-Bromo-5a-androstan-1-one
1878	21-Diazo-3,20-dioxo-17a-pregn-4-eno-19,8-lactone	292	2a-Bromo-5a-androstan-3-one
1669	21-Diazo-3b-acetoxy-5-pregnen-20-one	146	2a-Bromo-5a-androstan-3a-ol (impure)
1668	21-Diazo-3b-hydroxy-5-pregnen-20-one	595	2a-Bromo-5a-androstan-3,17-dione
1775	21-Diazo-8b,19-epoxy-4-pregnene-3,20-dione	873	2a-Hydroxy-4-androstene-3,17-dione
1938	21-Hydroxy-1,4-pregnadiene-3,20-dione	2023	2a-Methyl-4-pregnene-17a,21-diol-3,11,20-trione; 17a,21-Dihydroxy-2a-methyl-4-pregnene-3,11,20-trione
2068	21-Hydroxy-19-nor-14b,17a-pregn-4-ene-3,20-dione	432	2b,16b-Dihydroxy-5a-androstan-11-one
1903	21-Hydroxy-3,20-dioxo-14b,17a-pregn-4-eno-19,8-lactone	842	2b,17b-Diacetoxy-4-androsten-3-one; 2b,17b-Dihydroxy-4-androsten-3-one 2,17-diacetate
1902	21-Hydroxy-3,20-dioxo-4-pregneno-19,8-lactone	20	2b,17b-Diacetoxy-4-estren-3-one
1875	21-Hydroxy-4-pregnene-3,11,20-trione 21-glucuronide	424	2b,17b-Dihydroxy-5a-androstan-7-one
1876	21-Hydroxy-4-pregnene-3,11,20-trione glucuronide triacetyl methyl ester	2383	2b,3a-Di-(p-chlorobenzoyloxy)-5a-cholestane; 5a-Cholestane-2b,3a-diol bis(4-chlorobenzoate)
1877	21-Hydroxy-4-pregnene-3,11,20-trione glucuronide triacetyl methyl ester 3,20-disemicarbazone	2377	2b,3a-Diacetoxy-5a-cholestane; 5a-Cholestane-2b,3a-diol diacetate
1871	21-Hydroxy-4-pregnene-3,19,20-trione	2382	2b,3a-Dibenzoyloxy-5a-cholestane; 5a-Cholestane-2b,3a-diol dibenzoate
1747	21-Hydroxy-4-pregnene-3,20-dione	2767	2b,3b,14,20,22R,25-Hexahydroxy-5b-ergost-7-en-6-one; Makisterone A
1765	21-Hydroxy-4-pregnene-3,20-dione 21-glucuronide	2582	2b,3b,14a,20b,22R,25-Hexahydroxy-7-cholesten-6-one; b-Ecdysterone; b-Ecdysone (approx.95%)
1766	21-Hydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester	2580	2b,3b,14a,22R,25-Pentahydroxy-7-cholesten-6-one; a-Ecdysterone; a-Ecdysone (95-98%)
1767	21-Hydroxy-4-pregnene-3,20-dione 21-glucuronide triacetyl methyl ester 3,20-disemicarbazone	2584	2b,3b,5a,14,20,22R,25-Heptahydroxy-5b-cholest-7-en-6-one; Polygodine B
1761	21-Hydroxy-4-pregnene-3,20-dione 3-phenylpropionate	2431	2b-Acetoxy-5a-cholestan-1-one; 2b-Hydroxy-5a-cholestan-1-one acetate
1599	21-Hydroxy-5a-pregnane-3,6,20-trione	2369	2b-Acetoxy-5a-cholestan-1a-ol; 5a-Cholestane-1a,2b-diol 2-acetate
1502	21-Hydroxy-5b-pregnane-3,20-dione	2433	2b-Acetoxy-5a-cholestan-3-one; 2b-Hydroxy-5a-cholestan-3-one acetate
1312	21-Methoxy-5a-pregnan-20-one	2434	2b-Acetoxy-5a-cholestan-3-one; 2b-Hydroxy-5a-cholestan-3-one acetate
2053	21-Methylene-4,17(20)-pregnadien-3-one	1196	2b-Acetoxy-5a-pregnane; 5a-Pregnane-2b-ol acetate
1950	21-Pivaloyloxy-11b,17a-dihydroxy-1,4-pregnadiene-3,20-dione; Prednisolone 21-pivalate	2413	2b-Chloro-5a-cholestan-6-one
1760	21-Trimethylacetoxo-4-pregnene-3,20-dione; 21-Hydroxy-4-pregnene-3,20-dione pivalate	2423	2b-Hydroxy-5a-cholestan-3-one
2853	22,23-Dibromo-4-lumisten-3-one		
2941	22,23-Dibromo-5-oxo-A-nor-3,5-seco-22-lumisten-3-oic acid		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2639	2x-Hydroxy-3,3-dimethyl-A-nor-5-cholestene-2b-carboxylic acid	1595	3,3:20,20-Bis(ethylenedioxy)-5a-pregnan-11-one; 5a-Pregnane-3,11,20-trione 3,20-bis(ethylene ketal)
2640	2x-Hydroxy-3,3-dimethyl-A-nor-5-cholestene-2b-carboxylic acid methyl ester	1597	3,3:20,20-Bis(ethylenedioxy)-9a-bromo-5a-pregnan-11-one; 9a-Bromo-5a-pregnane-3,11,20-trione 3,20-bis(ethylene ketal)
2235	3,11-Dioxo-5b-cholan-24-oic acid methyl ester	685	3,3:7,7-Bis(ethylenedioxy)-5a-androstan-11-one; 5a-Androstane-3,7,11-trione 3,7-bis(ethylene ketal)
2317	3,12-Dioxo-5b-bisnorchol-9(11)-en-22-oic acid methyl ester	640	3,3:7,7-Bisethylenedioxy-1b,11a-epoxy-5a-androstane; 1b,11a-Epoxy-5a-androstane-3,7-dione bis(ethylene ketal)
2306	3,12-Dioxo-5b-bisnorcholan-22-oic acid methyl ester	635	3,3:7,7-Bisethylenedioxy-5a-androstan-11a-ol; 11a-Hydroxy-5a-androstane-3,7-dione bis(ethylene ketal)
2268	3,12-Dioxo-5b-chol-9(11)-en-24-oic acid methyl ester	636	3,3:7,7-Bisethylenedioxy-5a-androstan-11b-ol; 11b-Hydroxy-5a-androstane-3,7-dione bis(ethylene ketal)
2233	3,12-Dioxo-5b-cholan-24-oic acid	582	3,3:7,7-Bisethylenedioxy-5a-androstane; 5a-Androstane-3,7-dione bis(ethylene ketal)
2236	3,12-Dioxo-5b-cholan-24-oic acid methyl ester	2661	3,4-Seco-5-cholestene-3,4-dioic acid; Diel's acid
1138	3,15-Dioxo-14-hydroxy-5b,14b-etianic acid methyl ester	2654	3,4-Seco-5a-cholestane-3,4-dioic acid dimethyl ester; Dihydro Diel's ester
90	3,16a,17b-Trihydroxy-1,3,5(10)-estratrien-6-one; 6-Ketoestriol	2651	3,4-Seco-5b-cholestane-3,4-dioic acid
85	3,16a-Dihydroxy-1,3,5(10)-estratrien-17-one; 16a-Hydroxysterone	2657	3,4-Seco-A-bishomo-5a-cholestane-3,4-dioic acid
88	3,16b-Diacetoxy-1,3,5(10)-estratrien-17-one; 16b-Hydroxysterone diacetate	2656	3,4-Seco-A-homo-5a-cholestane-3,4-dioic acid
82	3,17b-Dihydroxy-1,3,5(10)-estratrien-16-one; 16-Ketoestradiol	925	3,5,8(9)-Androstatriene-7,17-dione (crude)
1172	3,19-Dioxoeti-4-enic acid ethyl ester	904	3,5-Androstadien-17-one (95-97%)
390	3,3-Difluoro-11a-hydroxy-5a-androstan-17-one	922	3,5-Androstadiene-1,7,17-trione
592	3,3-Difluoro-5a-androstane-11,17-dione	914	3,5-Androstadiene-2,7-dione
581	3,3-Dimethoxy-5a-androstan-17-one	900	3,5-Androstadiene-3,17b-diol dipropionate
267	3,3-Dimethoxy-5a-androstane	2645	3,5a-Dihydroxy-3,5-seco-A-norcholestane
268	3,3-Dimethoxy-5b-androstane	2646	3,5a-Seco-A-norcholestan-3-oic acid
2637	3,3-Dimethyl-A-nor-5-cholestene-2b-carboxylic acid	2455	3,6-Dihydroximino-5a-cholestane
2638	3,3-Dimethyl-A-nor-5-cholestene-2b-carboxylic acid methyl ester	2234	3,6-Dioxo-5a-cholan-24-oic acid methyl ester
548	3,3-Ethylenedioxy-11a,12b,15a-triacetoxy-5a-androstane	59	3,6a,17b-Triacetoxy-1,3,5(10)-estratriene; 6a-Hydroxyestradiol triacetate
958	3,3-Ethylenedioxy-12a-methyl-5b-androstan-11-one; 12a-Methyl-5b-androstane-3,11-dione 3-ethylene ketal	2246	3,7,12-Trioxo-5b-cholan-24-oic acid methyl ester
651	3,3-Ethylenedioxy-12b,15a-dihydroxy-5a-androstan-17-one	2240	3,7-Dioxo-5-hydroxy-5a-cholan-24-oic acid methyl ester
650	3,3-Ethylenedioxy-1b,6a-dihydroxy-5a-androstan-17-one; 1b,6a-Dihydroxy-5a-androstane-3,17-dione 3-ethylene ketal	2628	3-Acetoxy-4-methylcholesta-3,5-diene; 4-Methylcholesta-3,5-dien-3-ol acetate
759	3,3-Ethylenedioxy-5-androstene	2538	3-Acetoxy-5-cholesten-19-ol; 19-Hydroxycholesterol 3-acetate
575	3,3-Ethylenedioxy-5a-androstan-11-one; 5a-Androstane-3,11-dione 3-ethylene ketal	2593	3-Acetoxycholesta-3,5-diene; Cholesta-3,5-dien-3-ol acetate
483	3,3-Ethylenedioxy-5a-androstan-12b,15a-diol; 12b,15a-Dihydroxy-5a-androstan-3-one ethylene ketal	274	3-Aza-5a-androstan-2-one
577	3,3-Ethylenedioxy-5a-androstan-17-one; 5a-Androstane-3,17-dione 3-ethylene ketal	2598	3-Benzoyloxycholesta-4,6-diene; Cholesta-4,6-dien-3-ol benzoate
380	3,3-Ethylenedioxy-5a-androstan-17b-ol; 5a-Androstan-17b-ol-3-one ethylene ketal (crude)	2675	3-Benzylidene-A-homo-5a-cholestan-4-one
573	3,3-Ethylenedioxy-5a-androstan-6-one; 5a-Androstene-3,6-dione 3-ethylene ketal	73	3-Deoxysterone; 1,3,5(10)-Estratrien-17-one
574	3,3-Ethylenedioxy-5a-androstan-7-one; 5a-Androstane-3,7-dione 3-ethylene ketal (80% pure)	147	3-Doxyl-5a-androstan-17b-ol; 3-Doxyl-17b-hydroxy-5a-androstane
269	3,3-Ethylenedioxy-5a-androstane	2674	3-Furfurylidene-A-homo-5a-cholestan-4-one
482	3,3-Ethylenedioxy-5a-androstane-6a,12b-diol; 6a,12b-Dihydroxy-5a-androstan-3-one ethylene ketal	2562	3-Hydroximino-4-cholestene; 4-Cholesten-3-one oxime
2421	3,3-Ethylenedioxy-5a-cholestane; 5a-Cholestan-3-one ethylene ketal	264	3-Hydroximino-5a-androstane; 5a-Androstan-3-one oxime
2729	3,3-Ethylenedioxy-5a-ergostan-11-one; 5a-Ergostane-3,11-dione 3-ethylene ketal	2418	3-Hydroximino-5a-cholestane; 5a-Cholestan-3-one oxime
576	3,3-Ethylenedioxy-5b-androstan-11-one; 5b-Androstane-3,11-dione 3-ethylene ketal	2844	3-Hydroximino-5b-lumistane; 5b-Lumistan-3-one oxime
547	3,3-Ethylenedioxy-6a,11a,12b-triacetoxy-5a-androstane	79	3-Hydroxy-1,3,5(10)-estratrien-17-one b-D-glucuronide, sodium salt; Estrone b-D-glucuronide, sodium salt
484	3,3-Ethylenedioxy-6a,11a-diacetoxy-5a-androstane (crude)	1182	3-Hydroxy-19-nor-17(aH)-etia-1,3,5(10)-trienic acid ethyl ester
272	3,3-Ethylenedithio-5a-androstane	87	3-Hydroxy-2-methoxy-1,3,5(10)-estratrien-17-one; 2-Methoxysterone
896	3,3:17,17-Bis(ethylenedioxy)-5-androsten-11-one; 5-Androstene-3,11,17-trione 3,17-bis(ethylene ketal)	2676	3-Hydroxymethylene-A-homo-5a-cholestan-4-one
583	3,3:17,17-Bis(ethylenedioxy)-5a-androstane; 5a-Androstane-3,17-dione bis(ethylene ketal)	3003	3-Isopropylidene-4,4',4"-trisorlupane
		102	3-Methoxy-1,3,5(10),6,8(9)-estrapentaen-17-one
		907	3-Methoxy-3,5-androstadien-17-one
		898	3-Methoxy-3,5-androstadiene
		2604	3-Methoxycholesta-3,5-dien-7-one

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2592	3-Methoxycholesta-3,5-diene	1100	3a,12a-Dihydroxy-5b-etianic acid methyl ester 3-hemisuccinate
988	3-Methylene-5a-androstan-17-one	2288	3a,12a-Dihydroxy-5b-norcholan-23-oic acid methyl ester
972	3-Methylene-5a-androstan-17b-ol	1342	3a,12a-Dihydroxy-5b-pregnan-20-one
969	3-Methylene-5a-androstane (crude)	1119	3a,12a-Dihydroxy-7-oxo-5b-etianic acid methyl ester 3-hemisuccinate
2264	3-Oxo-5b-chol-11-en-24-oic acid methyl ester	2227	3a,12b-Diacetoxy-11-oxo-5b-cholan-24-oic acid methyl ester
2197	3-Oxo-5b-cholan-24-oic acid methyl ester; Methyl dehydrolithocholate	2293	3a,12b-Diacetoxy-5b,20R-bisnorcholan-22-oic acid methyl ester
2270	3-Oxochola-4,6-dien-24-oic acid	2292	3a,12b-Diacetoxy-5b-bisnorcholan-22-oic acid methyl ester
1919	3a,11,20b-Triacetoxy-5b-pregna-7,9(11)-diene; 5b-Pregna-7,9(11)-diene-3a,11,20b-triol triacetate	1105	3a,12b-Diacetoxy-5b-etianic acid methyl ester
510	3a,11a,17b-Trihydroxy-5a-androstan-6-one	2221	3a,12b-Dihydroxy-11-oxo-5b-cholan-24-oic acid
1393	3a,11a,21-Trihydroxy-5b-pregnan-20-one	2222	3a,12b-Dihydroxy-11-oxo-5b-cholan-24-oic acid methyl ester
649	3a,11a-Diacetoxy-5a-androstane-7,17-dione; 3a,11a-Dihydroxy-5a-androstane-7,17-dione 3,11-diacetate	1091	3a,12b-Dihydroxy-5a-etianic acid
452	3a,11a-Dihydroxy-5a-androstan-17-one	2285	3a,12b-Dihydroxy-5b,20R-bisnorcholan-22-oic acid methyl ester
1520	3a,11a-Dihydroxy-5a-pregnane-7,20-dione	2286	3a,12b-Dihydroxy-5b-bisnorcholan-22-oic acid methyl ester
1336	3a,11a-Dihydroxy-5b-pregnan-20-one	2162	3a,12b-Dihydroxy-5b-cholan-24-oic acid
1467	3a,11b,17a,21-Tetrahydroxy-5b-pregnan-20-one 21-glucuronide	1094	3a,12b-Dihydroxy-5b-etianic acid methyl ester
1468	3a,11b,17a,21-Tetrahydroxy-5b-pregnan-20-one 21-glucuronide methyl ester	1343	3a,16a-Dihydroxy-5a-pregnan-20-one
1458	3a,11b,17a,21-Tetrahydroxy-5b-pregnan-20-one 21-sulfate, sodium salt	459	3a,16a-Dihydroxy-5b-androstan-17-one; 16a-Hydroxyetiocholanolone
1455	3a,11b,17a,21-Tetrahydroxy-5b-pregnan-20-one; Tetrahydrocortisol	472	3a,16b-Diacetoxy-5a-androstan-7-one
1392	3a,11b,17a-Trihydroxy-5b-pregnan-20-one	1064	3a,17(a)b-Dihydroxy-17(a)-methyl-D-homo-5a-androstane-11,17-dione
1472	3a,11b,20a-Trihydroxy-5b-pregnan-21-oic acid methyl ester	1449	3a,17a,20a,21-Tetrahydroxy-5b-pregnan-11-one; a-Cortolone
1473	3a,11b,20b-Trihydroxy-5b-pregnan-21-oic acid methyl ester	1382	3a,17a,20a-Trihydroxy-5b-pregnan-11-one; Pregnanetriolone
1394	3a,11b,21-Trihydroxy-5a-pregnan-20-one	1465	3a,17a,20b,21-Tetrahydroxy-5b-pregnan-11-one 3-glucuronide
1429	3a,11b,21-Trihydroxy-5b-pregnan-20-one 21-glucuronide (free acid)	1466	3a,17a,20b,21-Tetrahydroxy-5b-pregnan-11-one 3-glucuronide methyl ester
1437	3a,11b,21-Trihydroxy-5b-pregnan-20-one 3,21-diglucuronide	1450	3a,17a,20b,21-Tetrahydroxy-5b-pregnan-11-one; b-Cortolone
1438	3a,11b,21-Trihydroxy-5b-pregnan-20-one 3,21-diglucuronide methyl ester	1417	3a,17a,20b-Triacetoxy-5b-pregnan-11-one
1433	3a,11b,21-Trihydroxy-5b-pregnan-20-one 3-glucuronide methyl ester	1383	3a,17a,20b-Trihydroxy-5b-pregnan-11-one
1431	3a,11b,21-Trihydroxy-5b-pregnan-20-one 3-glucuronide trihydrate, thallium salt	1400	3a,17a,21-Trihydroxy-5b-pregnan-20-one
1432	3a,11b,21-Trihydroxy-5b-pregnan-20-one 3a-glucuronide, barium acetate salt (1:1:1 complex)	1434	3a,17a,21-Trihydroxy-5b-pregnan-20-one 21-glucuronide methyl ester
1337	3a,11b-Dihydroxy-5a-pregnan-20-one	1430	3a,17a,21-Trihydroxy-5b-pregnan-20-one 21-glucuronide, barium salt
481	3a,11b-Dihydroxy-5b-androstan-17-one 3-glucuronide; 11b-Hydroxyetiocholanolone 3-glucuronide	1569	3a,17a,21-Trihydroxy-5b-pregnane-11,20-dione 21-glucuronide
462	3a,11b-Dihydroxy-5b-androstan-17-one 3-sulfate sodium salt	1572	3a,17a,21-Trihydroxy-5b-pregnane-11,20-dione 21-glucuronide methyl ester
454	3a,11b-Dihydroxy-5b-androstan-17-one; 11b-Hydroxyetiocholanolone	1560	3a,17a,21-Trihydroxy-5b-pregnane-11,20-dione 21-sulfate, sodium salt; Tetrahydrocortisone 21-sulfate, sodium salt
2291	3a,12a-Diacetoxy-5b,20R-bisnorcholan-22-oic acid methyl ester	1568	3a,17a,21-Trihydroxy-5b-pregnane-11,20-dione 3-glucuronide
2290	3a,12a-Diacetoxy-5b-bisnorcholan-22-oic acid methyl ester	1571	3a,17a,21-Trihydroxy-5b-pregnane-11,20-dione 3-glucuronide methyl ester
1104	3a,12a-Diacetoxy-5b-etianic acid methyl ester	1573	3a,17a,21-Trihydroxy-5b-pregnane-11,20-dione 3-glucuronide methyl ester 20-semicarbazone
2294	3a,12a-Diacetoxy-5b-norcholan-23-oic acid methyl ester	1346	3a,17a-Dihydroxy-5a-pregnan-20-one
2226	3a,12a-Diacetoxy-7-oxo-5b-cholan-24-oic acid methyl ester	1521	3a,17a-Dihydroxy-5a-pregnane-11,20-dione
1093	3a,12a-Dihydroxy-5b,17a-etianic acid methyl ester	1347	3a,17a-Dihydroxy-5b-pregnane-20-one
2283	3a,12a-Dihydroxy-5b,20R-bisnorcholan-22-oic acid	437	3a,17b-Dihydroxy-5a-androstan-11-one
480	3a,12a-Dihydroxy-5b-androstan-17-one 3,12-diacetate; 3a,12a-Diacetoxy-5b-androstan-17-one	444	3a,17b-Dihydroxy-5a-androstan-16-one
2287	3a,12a-Dihydroxy-5b-bisnorcholan-22-oic acid amide	422	3a,17b-Dihydroxy-5a-androstan-6-one
2284	3a,12a-Dihydroxy-5b-bisnorcholan-22-oic acid methyl ester	427	3a,17b-Dihydroxy-5a-androstan-7-one
2262	3a,12a-Dihydroxy-5b-chol-9(11)-en-24-oic acid methyl ester	438	3a,17b-Dihydroxy-5b-androstan-11-one
2177	3a,12a-Dihydroxy-5b-cholan-24-oic acid methyl ester 3-hemisuccinate	1063	3a,17b-Dihydroxy-D-homo-5b-androstane-11,17(a)-dione
1092	3a,12a-Dihydroxy-5b-etianic acid methyl ester	1403	3a,18,21-Trihydroxy-5b-pregnan-20-one 18->20-hemiketal

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1692	3a,20-Diacetoxy-5b-pregn-17(20)-en-11-one	2272	3a,24-Dihydroxy-24,24-diphenyl-5b-cholan-11-one
1425	3a,20a,21-Trihydroxy-5b-pregnane-11-one 3-glucuronide triacetyl methyl ester	2378	3a,4a-Diacetoxy-5a-cholestane; 5a-Cholestane-3a,4a-diol diacetate
1423	3a,20a,21-Trihydroxy-5b-pregnan-11-one 3-glucuronide, barium salt	2379	3a,4b-Diacetoxy-5a-cholestane; 5a-Cholestane-3a,4b-diol diacetate
1577	3a,20a-Diacetoxy-11-oxo-5b-pregnan-21-oic acid methyl ester	1614	3a,5,14-Trihydroxy-20-oxo-5b,17a-pregnane-19,21-dioic acid 21,14a-lactone
1409	3a,20a-Diacetoxy-17a-hydroxy-5b-pregnan-11-one	1033	3a,5-Cyclo-1b,6b-dihydroxy-5a-androstan-17-one
1228	3a,20a-Diacetoxy-5a-pregnane; 5a-Pregnane-3a,20a-diol diacetate	1029	3a,5-Cyclo-5a-androstan-17-one
1264	3a,20a-Dibenzoyloxy-5b-pregnan-11b-ol; 5b-Pregnane-3a,11b,20a-triol 3,20-dibenzoate	1028	3a,5-Cyclo-5a-androstan-6-one
1422	3a,20b,21-Trihydroxy-5b-pregnan-11-one 3-glucuronide	1035	3a,5-Cyclo-5a-androstane-6,17-dione
1424	3a,20b,21-Trihydroxy-5b-pregnan-11-one 3-glucuronide methyl ester	1027	3a,5-Cyclo-5a-androstane-6a,17b-diol
1426	3a,20b,21-Trihydroxy-5b-pregnan-11-one 3-glucuronide triacetyl methyl ester	2058	3a,5-Cyclo-5a-pregnan-20-one
1578	3a,20b-Diacetoxy-11-oxo-5b-pregnan-21-oic acid methyl ester	2821	3a,5-Epoxy-5a-ergosta-7,9(11),22-triene
1474	3a,20b-Diacetoxy-11b-hydroxy-5b-pregnan-21-oic acid methyl ester	2782	3a,5a-Diacetoxy-5a-ergosta-7,22-diene; 5a-Ergosta-7,22-diene-3a,5a-diol diacetate
1410	3a,20b-Diacetoxy-17a-hydroxy-5b-pregnan-11-one	1587	3a,5b-Tetrahydroaldosterone g-lactone
1361	3a,20b-Diacetoxy-5b-pregnan-11-one	1464	3a,6a,17a,21-Tetraacetoxy-5b-pregnan-20-one
1256	3a,20b-Diacetoxy-5b-pregnan-17a-ol; 5b-Pregnane-3a,17a,20b-triol 3,20-diacetate	1453	3a,6a,17a,21-Tetrahydroxy-5b-pregnan-20-one
1279	3a,20b-Diacetoxy-5b-pregnane-11b,17a-diol; 5b-Pregnane-3a,11b,17a,20b-tetrol 3,20-diacetate	1386	3a,6a,17a-Trihydroxy-5b-pregnan-20-one
1230	3a,20b-Diacetoxy-5b-pregnane; 5b-Pregnane-3a,20b-diol diacetate	1387	3a,6a,21-Trihydroxy-5b-pregnan-20-one
1576	3a,20b-Dihydroxy-11-oxo-5b-pregnan-21-oic acid methyl ester	1553	3a,6a,21-Trihydroxy-5b-pregnane-11,20-dione
1331	3a,20b-Dihydroxy-5b-pregnan-11-one	2187	3a,6a,7a-Trihydroxy-5b-cholan-24-oic acid methyl ester; Methyl hyocholate
1460	3a,21-Diacetoxy-17a,20b-dihydroxy-5b-pregnan-11-one	2184	3a,6a,7a-Trihydroxy-5b-cholan-24-oic acid; Hyocholic acid
1415	3a,21-Diacetoxy-17a-hydroxy-5b-pregnan-20-one	1363	3a,6a-Diacetoxy-5b-pregnan-20-one
1365	3a,21-Diacetoxy-5b-pregnan-20-one	449	3a,6a-Dihydroxy-5a-androstan-17-one
1532	3a,21-Diacetoxy-5b-pregnane-11,20-dione	450	3a,6a-Dihydroxy-5b-androstan-17-one; 6a-Hydroxyetiocholanolone
1350	3a,21-Dihydroxy-5a-pregnan-20-one; Tetrahydro-doc	2166	3a,6a-Dihydroxy-5b-cholan-24-oic acid N-(carboxymethyl)amide; Glycohyodeoxycholic acid
1523	3a,21-Dihydroxy-5a-pregnane-11,20-dione	2170	3a,6a-Dihydroxy-5b-cholan-24-oic acid N-(sulfoethyl)amide; Taurohyodeoxycholic acid
1351	3a,21-Dihydroxy-5b-pregnan-20-one	2157	3a,6a-Dihydroxy-5b-cholan-24-oic acid; Hyodeoxycholic acid (approx.95%)
1368	3a,21-Dihydroxy-5b-pregnan-20-one 21-glucuronide	1333	3a,6b-Dihydroxy-5a-pregnan-20-one
1370	3a,21-Dihydroxy-5b-pregnan-20-one 21-glucuronide methyl ester	2158	3a,6b-Dihydroxy-5b-cholan-24-oic acid; Murocholic acid
1369	3a,21-Dihydroxy-5b-pregnan-20-one 3-glucuronide methyl ester	1334	3a,6b-Dihydroxy-5b-pregnan-20-one
1372	3a,21-Dihydroxy-5b-pregnan-20-one 3a, 21-di(glucuronide methyl ester)	2296	3a,7a,12a-Trihydroxy-24-nor-5b-cholan-23-oic acid ethyl ester; Ethyl norcholate
1367	3a,21-Dihydroxy-5b-pregnan-20-one 3a-glucuronide	2295	3a,7a,12a-Trihydroxy-24-nor-5b-cholan-23-oic acid methyl ester; Methyl norcholate
1535	3a,21-Dihydroxy-5b-pregnane-11,20-dione 21-glucuronide methyl ester	2188	3a,7a,12a-Trihydroxy-5b-cholan-24-oic acid methyl ester; Methyl cholate
1537	3a,21-Dihydroxy-5b-pregnane-11,20-dione 21-glucuronide triacetyl methyl ester 20-semicarbazone	2190	3a,7a,12a-Trihydroxy-5b-cholan-24-oic acid N-(carboxymethyl)amide, sodium salt; Glycocholic acid, sodium salt
1541	3a,21-Dihydroxy-5b-pregnane-11,20-dione 3, 21-diglucuronide	2186	3a,7a,12a-Trihydroxy-5b-cholan-24-oic acid, sodium salt; Cholic acid, sodium salt
1542	3a,21-Dihydroxy-5b-pregnane-11,20-dione 3, 21-diglucuronide, diammonium salt	2185	3a,7a,12a-Trihydroxy-5b-cholan-24-oic acid; Cholic acid
1533	3a,21-Dihydroxy-5b-pregnane-11,20-dione 3-glucuronide	1389	3a,7a,12a-Trihydroxypregnan-20-one (approx.95%)
1536	3a,21-Dihydroxy-5b-pregnane-11,20-dione 3-glucuronide triacetyl methyl ester	2228	3a,7a-Diacetoxy-12-oxo-5b-cholan-24-oic acid methyl ester
1534	3a,21-Dihydroxy-5b-pregnane-11,20-dione 3-glucuronide, ammonium salt	2192	3a,7a-Diacetoxy-12a-hydroxy-5b-cholan-24-oic acid
1543	3a,21-Dihydroxy-5b-pregnane-11,20-dione 3a,21-di(glucuronide triacetyl methyl ester) 20-semicarbazone	2159	3a,7a-Dihydroxy-5b-cholan-24-oic acid; Chenodeoxycholic acid
2133	3a,24-Diacetoxy-5a-chol-11-ene; 5a-Chol-11-ene-3a,24-diol diacetate	2229	3a,7b-Diacetoxy-12-oxo-5b-cholan-24-oic acid methyl ester
		2304	3a,9a-Epoxy-11-oxo-23,24-bisnor-5b-cholan-22-oic acid
		2305	3a,9a-Epoxy-11-oxo-23,24-bisnor-5b-cholan-22-oic acid methyl ester
		2219	3a,9a-Epoxy-11-oxo-5b-cholan-24-oic acid methyl ester
		2260	3a,9a-Epoxy-5b-chol-11-en-24-oic acid methyl ester
		370	3a-(p-Toluenesulfonyloxy)-5a-androstan-17-one; 5a-Androstan-3a-ol-17-one tosylate

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2154	3a-Acetoxy-11,12-dibromo-5b-cholan-24-oic acid methyl ester	1149	3a-Acetoxy-5b-eti-9(11)-enic acid methyl ester
2211	3a-Acetoxy-11-oxo-5b-cholan-24-oic acid methyl ester	1354	3a-Acetoxy-5b-pregnan-12a-ol-20-one
2181	3a-Acetoxy-11a,12a-epoxy-5b-cholan-24-oic acid methyl ester	1314	3a-Acetoxy-5b-pregnan-20-one
2217	3a-Acetoxy-11a-bromo-12-oxo-5b-cholan-24-oic acid methyl ester	1353	3a-Acetoxy-5b-pregnane-20b-ol-11-one; 3a,20b-Dihydroxy-5b-pregnan-11-one 3-acetate
2289	3a-Acetoxy-11a-hydroxy-5b-bisnorcholan-22-oic acid methyl ester	2210	3a-Acetoxy-6-oxo-5a-cholan-24-oic acid methyl ester
2182	3a-Acetoxy-11b,12b-epoxy-5b-cholan-24-oic acid methyl ester	2242	3a-Acetoxy-7,12-dioxo-5b-cholan-24-oic acid methyl ester
1469	3a-Acetoxy-11b,17a,21-trihydroxy-5b-pregnan-20-one 21-glucuronide triacetyl methyl ester	2899	3a-Acetoxy-5,7,22-triene (crude); Lumista-5,7,22-trien-3a-ol acetate
2218	3a-Acetoxy-11b-bromo-12-oxo-5b-cholan-24-oic acid methyl ester	1319	3a-Benzoyloxy-5a-pregnan-20-one; 5a-Pregnan-3a-ol-20-one 3-benzoate
2173	3a-Acetoxy-11b-hydroxy-5b-cholan-24-oic acid methyl ester	1511	3a-Benzoyloxy-5b-pregnane-11,20-dione
1096	3a-Acetoxy-11b-hydroxy-5b-etianic acid methyl ester	500	3a-Bromo-1b,6a-dihydroxy-5a-androstan-17-one
2302	3a-Acetoxy-12-oxo-5b-bisnorcholan-22-oic acid methyl ester	940	3a-Bromo-4,4-dimethyl-5a-androstan-2-one
2266	3a-Acetoxy-12-oxo-5b-chol-9(11)-en-24-oic acid methyl ester	293	3a-Bromo-5a-androstan-17-one
2271	3a-Acetoxy-12-oxo-5b-chola-7,9(11)-dien-24-oic acid methyl ester	241	3a-Chloro-1b,6a,17b-triacetoxy-5a-androstane
2213	3a-Acetoxy-12-oxo-5b-cholan-24-oic acid methyl ester	287	3a-Chloro-5a-androstan-17-one
1168	3a-Acetoxy-12-oxo-5b-eti-9(11)-enic acid methyl ester	286	3a-Chloro-5a-androstan-7-one
394	3a-Acetoxy-12a-bromo-5b-androstan-11-one	240	3a-Chloro-5a-androstane-1b,15a,17b-triol
2174	3a-Acetoxy-12a-hydroxy-5b-cholan-24-oic acid methyl ester; Methyl deoxycholate 3-acetate	2343	3a-Chloro-5a-cholestan-5-ol
1097	3a-Acetoxy-14a-hydroxy-5b-etianic acid methyl ester	498	3a-Chloro-6a,11a-dihydroxy-5a-androstan-17-one
1098	3a-Acetoxy-14b-hydroxy-5b,14b-etianic acid methyl ester	132	3a-Chloroacetoxy-5a-androstane
1435	3a-Acetoxy-17a,21-dihydroxy-5b-pregnan-20-one 21-glucuronide triacetyl methyl ester	2467	3a-Dimethylamino-5-cholestene
1574	3a-Acetoxy-17a,21-dihydroxy-5b-pregnane-11,20-dione 21-glucuronide triacetyl methyl ester 20-semicarbazone	388	3a-Fluoro-11a-hydroxy-5a-androstan-17-one
1355	3a-Acetoxy-17a-hydroxy-5a-pregnan-20-one	490	3a-Fluoro-12b,15a-dihydroxy-5a-androstan-17-one
837	3a-Acetoxy-1a-hydroxy-5-androsten-17-one; 1a,3a-Dihydroxy-5-androsten-17-one 3-acetate (crude)	486	3a-Fluoro-1b,15a-dihydroxy-5a-androstan-17-one
1358	3a-Acetoxy-21-hydroxy-5b-pregnan-20-one	493	3a-Fluoro-1b,6a-diacetoxy-5a-androstan-17-one
1371	3a-Acetoxy-21-hydroxy-5b-pregnan-20-one 21-glucuronide triacetyl methyl ester	278	3a-Fluoro-5a-androstan-17-one
1527	3a-Acetoxy-21-hydroxy-5b-pregnane-11,20-dione (1:1 adduct with benzene)	277	3a-Fluoro-5a-androstan-7-one
1538	3a-Acetoxy-21-hydroxy-5b-pregnane-11,20-dione 21-glucuronide triacetyl methyl ester	686	3a-Fluoro-5a-androstane-1,6,17-trione
1540	3a-Acetoxy-21-hydroxy-5b-pregnane-11,20-dione 21-glucuronide triacetyl methyl ester 20-semicarbazone	591	3a-Fluoro-5a-androstane-11,17-dione
2274	3a-Acetoxy-24,24-diphenyl-5b-chol-23-ene	687	3a-Fluoro-5a-androstane-6,11,16-trione
2278	3a-Acetoxy-24,24-diphenyl-5b-chola-20,23-dien-11-one	688	3a-Fluoro-5a-androstane-7,11,17-trione
2276	3a-Acetoxy-24,24-diphenyl-5b-chola-20,23-diene	487	3a-Fluoro-6a,11a-dihydroxy-5a-androstan-17-one
2273	3a-Acetoxy-24-phenyl-5b-cholane-11,24-dione	485	3a-Fluoro-7b,11a-dihydroxy-5a-androstan-16-one
1177	3a-Acetoxy-5,16-etiadiene 17-nitrile	489	3a-Fluoro-7b,11a-dihydroxy-5a-androstan-17-one
179	3a-Acetoxy-5a-androstan-17a-ol	2298	3a-Hydroxy-11-oxo-23,24-bisnor-5b-cholan-22-oic acid
2435	3a-Acetoxy-5a-cholestan-6-one; 3a-Hydroxy-5a-cholestan-6-one acetate	2299	3a-Hydroxy-11-oxo-23,24-bisnor-5b-cholan-22-oic acid methyl ester
2777	3a-Acetoxy-5a-ergosta-7,22-dien-5-ol; 5a-Ergosta-7,22-diene-3a,5-diol 3-acetate	2204	3a-Hydroxy-11-oxo-5b-cholan-24-oic acid
1088	3a-Acetoxy-5a-etianic acid	1106	3a-Hydroxy-11b,12b-epoxy-5b-etianic acid methyl ester
1313	3a-Acetoxy-5a-pregnan-20-one; 5a-Pregnan-3a-ol-20-one 3-acetate	2205	3a-Hydroxy-12-oxo-5b-cholan-24-oic acid
2961	3a-Acetoxy-5a-stigmast-8(14)-ene; 5a-Stigmast-8(14)-en-3a-ol acetate	2208	3a-Hydroxy-12-oxo-5b-cholan-24-oic acid methyl ester
2257	3a-Acetoxy-5b-chol-11-en-24-oic acid methyl ester	2275	3a-Hydroxy-24,24-diphenyl-5b-chol-23-en-11-one
2256	3a-Acetoxy-5b-chol-9(11)-en-24-oic acid methyl ester	776	3a-Hydroxy-5-androsten-17-one; 5-Androsten-3a-ol-17-one
2152	3a-Acetoxy-5b-cholan-24-oic acid methyl ester	315	3a-Hydroxy-5a-androstan-7-one; 5a-Androstan-3a-ol-7-one
2372	3a-Acetoxy-5b-cholestan-5-ol; 5b-Cholestane-3a,5-diol 3-acetate	630	3a-Hydroxy-5a-androstane-11,17-dione b-D-glucuronide; 11-Ketoandrosterone b-D-glucuronide
		625	3a-Hydroxy-5a-androstane-11,17-dione; 11-Ketoandrosterone
		621	3a-Hydroxy-5a-androstane-6,17-dione
		2249	3a-Hydroxy-5a-chol-11-en-24-oic acid
		2429	3a-Hydroxy-5a-cholestan-7-one
		1301	3a-Hydroxy-5a-pregnan-20-one; 5a-Pregnan-3a-ol-20-one
		1504	3a-Hydroxy-5a-pregnane-11,20-dione; Alphaxolane
		317	3a-Hydroxy-5b-androstan-11-one; 5b-Androstan-3a-ol-11-one
		322	3a-Hydroxy-5b-androstan-17-one; Etiocholanolone
		627	3a-Hydroxy-5b-androstane-11,17-dione 3-sulfate, triethylammonium salt
		631	3a-Hydroxy-5b-androstane-11,17-dione glucuronide; 11-Ketoetiocholanolone glucuronide

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2252	3a-Hydroxy-5b-chole-11-en-24-oic acid methyl ester	1108	3b,12b-Diacetoxy-14b-hydroxy-5b,14b-etianic acid methyl ester
2142	3a-Hydroxy-5b-cholan-24-oic acid; Lithocholic acid	2620	3b,12b-Diacetoxy-5a,25(R)-spirost-9(11)-ene; 5a,25(R)-Spirost-9(11)-ene-3b,12b-diol diacetate
1142	3a-Hydroxy-5b-eti-9(11)-enic acid methyl ester	2613	3b,12b-Diacetoxy-5a,25(R)-spirostane; 5a,25(R)-Spirostane-3b,12b-diol diacetate
1083	3a-Hydroxy-5b-etianic acid	1630	3b,12b-Diacetoxy-5b-pregn-20-ene; 5b-Pregn-20-ene-3b,12b-diol diacetate
1602	3a-Hydroxy-5b-pregnane-7,12,20-trione	822	3b,12b-Dihydroxy-5-androsten-7-one
1042	3a-Hydroxy-A-nor(5->6b)abeo-androstan-5-one	1044	3b,12b-Dihydroxy-A-nor(5->6b)abeo-androstan-5-one
1071	3a-Hydroxy-des-D-5b-androst-12-en-14-one	2581	3b,14,20,22R,25-Pentahydroxy-5b-cholest-7-en-6-one; 2-Deoxy-b-ecdysterone
299	3a-Iodo-5a-androstan-17-one	517	3b,14a,16b-Trihydroxy-5a-androstan-11-one
384	3a-Methoxy-17,17-ethylenedioxy-5a-androstane; 3a-Methoxy-5a-androstan-17-one ethylene ketal	828	3b,14a-Dihydroxy-5-androsten-17-one; 14a-Hydroxydehydroepiandrosterone
340	3a-Methoxy-5a-androstan-17-one	823	3b,14a-Dihydroxy-5-androsten-7-one
223	3a-Methoxy-5a-androstane-1b,6a-diol	2579	3b,14b,22R,25-Tetrahydroxy-5b-cholest-7-en-6-one; 2-Deoxy-a-ecdysterone
1644	3b,11,20b-Triacetoxy-5a-pregn-9(11)-ene; 5a-Pregn-9(11)-ene-3b,11,20b-triol triacetate	234	3b,15a,17b-Triacetoxy-5a-androstane; 5a-Androstane-3b,15a,17b-triol triacetate
1918	3b,11-Diacetoxy-20x-methoxy-5a-pregna-7,9(11)-diene; 5a-Pregna-7,9(11)-diene-3b,11,20x-triol 3,11-diacetate 20-methyl ether	1125	3b,15a-Diacetoxy-14b-hydroxy-11-oxo-5b,14b-etianic acid methyl ester
2826	3b,11-Diacetoxy-5a-ergosta-7,9(11),22-triene; 5a-Ergosta-7,9(11),22-triene-3b,11-diol diacetate	1412	3b,15a-Diacetoxy-9a-hydroxy-5a-pregnan-20-one
232	3b,11a,16b-Triacetoxy-5a-androstane; 5a-Androstane-3b,11a,16b-triol triacetate	435	3b,15a-Dihydroxy-5a-androstan-11-one
536	3b,11a,17b-Triacetoxy-5a-androstan-7-one	420	3b,15a-Dihydroxy-5a-androstan-6-one
233	3b,11a,17b-Triacetoxy-5a-androstane; 5a-Androstane-3b,11a,17b-triol triacetate	1713	3b,15b-Diacetoxy-14a-hydroxy-5a-pregnan-16-en-20-one (crude)
1113	3b,11a-Diacetoxy-14a,15a-epoxy-5b-etianic acid methyl ester	421	3b,15b-Dihydroxy-5a-androstan-6-one
1447	3b,11a-Diacetoxy-5,6b-epoxy-5b-pregnan-20-one	1930	3b,16a-Diacetoxy-17a-hydroxypregna-5,14-dien-20-one
2698	3b,11a-Diacetoxy-5a,9b-ergostane; 5a,9b-Ergostane-3b,11a-diol diacetate	845	3b,16a-Diacetoxy-5-androsten-17-one; 3b,16a-Dihydroxy-5-androsten-17-one 3,16-diacetate
2697	3b,11a-Diacetoxy-5a-ergostane; 5a-Ergostane-3b,11a-diol diacetate	473	3b,16a-Diacetoxy-5a-androstan-11-one (crude)
1531	3b,11a-Diacetoxy-5a-pregnane-7,20-dione	470	3b,16a-Diacetoxy-5a-androstan-6-one
1443	3b,11a-Dihydroxy-5,6b-epoxy-5b-pregnan-20-one	829	3b,16a-Dihydroxy-5-androsten-17-one
826	3b,11a-Dihydroxy-5-androsten-17-one	832	3b,16a-Dihydroxy-5-androsten-17-one 3-sulfate; 16a-Hydroxydehydroepiandrosterone sulfate
1676	3b,11a-Dihydroxy-5-pregnen-20-one	1792	3b,16a-Dihydroxy-5-pregnene-11,20-dione
443	3b,11a-Dihydroxy-5a-androstan-16-one	436	3b,16a-Dihydroxy-5a-androstan-11-one
455	3b,11a-Dihydroxy-5a-androstan-17-one	1344	3b,16a-Dihydroxy-5a-pregnan-20-one
425	3b,11a-Dihydroxy-5a-androstan-7-one	1345	3b,16a-Dihydroxy-5b-pregnan-20-one
1339	3b,11a-Dihydroxy-5a-pregnan-20-one (crude)	740	3b,16b,17b-Triacetoxy-5-androstene; 5-Androstene-3b,16b,17b-triol triacetate
1340	3b,11a-Dihydroxy-5b-pregnan-20-one	535	3b,16b,18-Triacetoxy-5a-androstan-6-one
1066	3b,11a-Dihydroxy-D-homo-5-androsten-17(a)-one	537	3b,16b,18-Triacetoxy-5a-androstan-7-one
2808	3b,11a-Dihydroxyergosta-5,8-dien-7-one	512	3b,16b,18-Trihydroxy-5a-androstan-7-one
1557	3b,11b,21-Trihydroxy-5a-pregnan-18-al-20-one; 3a,5a-Tetrahydroaldosterone	846	3b,16b-Diacetoxy-5-androsten-17-one; 3b,16b-Dihydroxy-5-androsten-17-one 3,16-diacetate
1396	3b,11b,21-Trihydroxy-5a-pregnan-20-one; Reichstein's compound R	528	3b,16b-Diacetoxy-5-hydroxy-5a-androstan-11-one
2785	3b,11b-Diacetoxy-5a-ergosta-7,22-diene; 5a-Ergosta-7,22-diene-3b,11b-diol diacetate	474	3b,16b-Diacetoxy-5a-androstan-11-one
2699	3b,11b-Diacetoxy-5a-ergostane; 5a-Ergostane-3b,11b-diol diacetate	830	3b,16b-Dihydroxy-5-androsten-17-one
827	3b,11b-Dihydroxy-5-androsten-17-one; 11b-Hydroxydehydroepiandrosterone	426	3b,16b-Dihydroxy-5a-androstan-7-one
1677	3b,11b-Dihydroxy-5-pregnen-20-one	1062	3b,17(a)x-Diacetoxy-17(a)-methyl-D-homo-5a-androstan-17-one
1341	3b,11b-Dihydroxy-5a-pregnan-20-one	1060	3b,17(a)x-Dihydroxy-17(a)-methyl-D-homo-5a-androstan-17-one
821	3b,12a-Dihydroxy-5-androsten-7-one	1708	3b,17a,21-Trihydroxy-5-pregnen-20-one
522	3b,12b,14a-Trihydroxy-5a-androstan-17-one	1828	3b,17a,21-Trihydroxy-5-pregnene-11,20-dione
545	3b,12b,15a-Triacetoxy-5a-androstan-17-one	183	3b,17a-Diacetoxy-5a-androstane; 5a-Androstane-3b,17a-diol diacetate
1421	3b,12b,15a-Triacetoxy-5a-pregnan-20-one	2004	3b,17a-Dihydroxy-16a-methyl-5-pregnen-20-one
523	3b,12b,15a-Trihydroxy-5a-androstan-17-one	1984	3b,17a-Dihydroxy-16b-methyl-5a-pregnan-20-one
1398	3b,12b,15a-Trihydroxy-5a-pregnan-20-one	1444	3b,17a-Dihydroxy-5,6a-epoxy-5a-pregnan-20-one
511	3b,12b,17b-Trihydroxy-5a-androstan-7-one(crude)		
1114	3b,12b-Diacetoxy-14b,15b-epoxy-5b,14b-etianic acid methyl ester		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1682	3b,17a-Dihydroxy-5-pregnen-20-one 17-sulfate, potassium salt	1160	3b,19-Diacetoxy-17b-eti-4-enic acid ethyl ester
1684	3b,17a-Dihydroxy-5-pregnen-20-one 3,17-disulfate, diammonium salt	1158	3b,19-Diacetoxy-17b-eti-4-enic acid methyl ester
1685	3b,17a-Dihydroxy-5-pregnen-20-one 3,17-disulfate, dipotassium salt	1442	3b,19-Diacetoxy-21-diazo-5-hydroxy-5b-pregnan-20-one
1681	3b,17a-Dihydroxy-5-pregnen-20-one 3-sulfate, triethylammonium salt	1462	3b,19-Diacetoxy-5,14-dihydroxy-5b,14b,17a-pregnan-20-one
1678	3b,17a-Dihydroxy-5-pregnen-20-one; 17a-Hydroxypregnenolone	1159	3b,19-Diacetoxy-5-eticnic acid ethyl ester
1791	3b,17a-Dihydroxy-5-pregnene-11,20-dione	1582	3b,19-Diacetoxy-5-hydroxy-20-oxo-5b,14b-pregnano-21,14-lactone
1349	3b,17a-Dihydroxy-5b-pregnan-20-one	1109	3b,19-Diacetoxy-5-hydroxy-5b,14b,17a-eticnic acid ethyl ester
1348	3b,17a-Hydroxy-5a-pregnan-20-one	1411	3b,19-Diacetoxy-5-hydroxy-5b-pregnan-20-one
854	3b,17b,18-Trihydroxy-5-androsten-16-one	1151	3b,19-Dihydroxy-17b-eti-4-enic acid
980	3b,17b-Diacetoxy-11-methylene-5a-androstane; 11-Methylene-5a-androstane-3b,17b-diol diacetate	1154	3b,19-Dihydroxy-17b-eti-4-enic acid ethyl ester
934	3b,17b-Diacetoxy-11a-methyl-5a-androstan-11b-ol; 11a-Methyl-5a-androstane-3b,11b,17b-triol 3,17-diacetate	1152	3b,19-Dihydroxy-17b-eti-4-enic acid methyl ester
982	3b,17b-Diacetoxy-16-benzylidene-5a-androstane; 16-Benzylidene-5a-androstane-3b,17b-diol diacetate	1178	3b,19-Dihydroxy-5,14-etiadienic acid ethyl ester
981	3b,17b-Diacetoxy-17-methyl-5-androstene; 17-Methyl-5-androstene-3b,17b-diol diacetate	831	3b,19-Dihydroxy-5-androsten-17-one
1157	3b,17b-Diacetoxy-17a-eti-5-enic acid	834	3b,19-Dihydroxy-5-androsten-17-one 19-sulfate, potassium salt
1695	3b,17b-Diacetoxy-17a-pregn-5-en-20-one	833	3b,19-Dihydroxy-5-androsten-17-one 3-sulfate, sodium salt
1971	3b,17b-Diacetoxy-17a-pregn-5-en-20-yne; 17a-Pregn-5-en-20-yne-3b,17b-diol diacetate	835	3b,19-Dihydroxy-5-androsten-17-one disulfate, dipotassium salt
1718	3b,17b-Diacetoxy-17a-pregn-5-en-21-oic acid methyl ester	1693	3b,20-Diacetoxy-5a-pregn-17(20)-en-11-one
2077	3b,17b-Diacetoxy-19-nor-4-pregnen-20-yne; 3b,17b-Diacetoxy-17a-ethynyl-4-estrene	1648	3b,20a,21-Triacetoxy-17a-pregn-5-en-17b-ol; 17a-Pregn-5-ene-3b,17b,20a,21-tetrol 3,20,21-triacetate
843	3b,17b-Diacetoxy-5-androsten-11-one; 3b,17b-Dihydroxy-5-androsten-11-one 3,17-diacetate	1281	3b,20a,21-Triacetoxy-5a,17a-pregnan-17b-ol; 5a,17a-Pregnane-3b,17b,20a,21-tetrol 3,20,21-triacetate
738	3b,17b-Diacetoxy-5-androsten-11b-ol; 5-Androstene-3b,11b,17b-triol 3,17-diacetate	1262	3b,20a,21-Triacetoxy-5a-pregnane; 5a-Pregnane-3b,20a,21-triol triacetate
844	3b,17b-Diacetoxy-5-androsten-16-one; 3b,17b-Dihydroxy-5-androsten-16-one 3,17-diacetate	1642	3b,20a-Diacetoxy-17a-pregn-5-en-17b-ol; 17a-Pregn-5-ene-3b,17b,20a-triol 3,20-diacetate
1102	3b,17b-Diacetoxy-5a,17a-eticnic acid	1641	3b,20a-Diacetoxy-5-pregnen-17a-ol; 5-Pregnene-3b,17a,20a-triol 3,20-diacetate
1629	3b,17b-Diacetoxy-5a,17a-pregn-20-ene; 5a,17a-Pregn-20-ene-3b,17b-diol diacetate	1257	3b,20a-Diacetoxy-5a-pregnan-17a-ol; 5a-Pregnane-3b,17a,20a-triol 3,20-diacetate
1227	3b,17b-Diacetoxy-5a,17a-pregnane; 5a,17a-Pregnane-3b,17b-diol diacetate	1915	3b,20a-Diacetoxypregna-5,16-diene; Pregna-5,16-diene-3b,20a-diol diacetate
475	3b,17b-Diacetoxy-5a-androstan-11-one	1632	3b,20a-Dibenzoyloxy-5-pregnene; 5-Pregnene-3b,20a-diol dibenzoate
477	3b,17b-Diacetoxy-5a-androstan-16-one	1265	3b,20a-Dibenzoyloxy-5a-pregnan-17a-ol; 5a-Pregnane-3b,17a,20a-triol 3,20-dibenzoate
226	3b,17b-Diacetoxy-5a-androstan-1a-ol; 5a-Androstane-1a,3b,17b-triol 3,17-diacetate	1282	3b,20b,21-Triacetoxy-5a,17a-pregnan-17b-ol; 5a,17a-Pregnane-3b,17b,20b,21-tetrol 3,20,21-triacetate
468	3b,17b-Diacetoxy-5a-androstan-2-one	1263	3b,20b,21-Triacetoxy-5a-pregnane; 5a-Pregnane-3b,20b,21-triol triacetate
471	3b,17b-Diacetoxy-5a-androstan-6-one	1978	3b,20b-Diacetoxy-11a-methyl-5a-pregnan-11b-ol; 11a-Methyl-5a-pregnane-3b,11b,20b-triol 3,20-diacetate
2310	3b,17b-Dihydroxy-17a,20x-bisnor-5-cholen-22-oic acid methyl ester	1376	3b,20b-Diacetoxy-12a-bromo-5a-pregnan-11-one
1153	3b,17b-Dihydroxy-17a-eti-5-enic acid methyl ester	1280	3b,20b-Diacetoxy-21-methoxy-5a-pregnan-17a-ol; 5a-Pregnane-3b,17a,20b,21-tetrol 3,20-diacetate 21-methyl ether
1717	3b,17b-Dihydroxy-17a-pregn-5-en-21-oic acid methyl ester	1259	3b,20b-Diacetoxy-5a,17a-pregnan-17b-ol; 5a,17a-Pregnane-3b,17b,20b-triol 3,20-diacetate
901	3b,17b-Dihydroxy-3,5-androstadiene bis(trimethylsilyl ether)	1362	3b,20b-Diacetoxy-5a-pregnan-11-one
824	3b,17b-Dihydroxy-5-androsten-16-one	1255	3b,20b-Diacetoxy-5a-pregnan-11b-ol; 5a-Pregnane-3b,11b,20b-triol 3,20-diacetate
1095	3b,17b-Dihydroxy-5a,17a-eticnic acid methyl ester	1258	3b,20b-Diacetoxy-5a-pregnan-17a-ol; 5a-Pregnane-3b,17a,20b-triol 3,20-diacetate
423	3b,17b-Dihydroxy-5a-androstan-6-one	1231	3b,20b-Diacetoxy-5b-pregnane; 5b-Pregnane-3b,20b-diol diacetate
847	3b,18-Diacetoxy-5-androsten-17-one; 3b,18-Dihydroxy-5-androsten-17-one 3,18-diacetate	1986	3b,20b-Diacetoxy-9a,12b-dimethyl-5a-pregnan-11-one
739	3b,18-Diacetoxy-5-androsten-17b-ol; 5-Androstene-3b,17b,18-triol 3,18-diacetate	1375	3b,20b-Diacetoxy-9a-bromo-5a-pregnan-11-one
478	3b,18-Diacetoxy-5a-androstan-16-one	1408	3b,20b-Diacetoxy-9a-hydroxy-5a-pregnan-11-one
1679	3b,18-Dihydroxy-5-pregnen-20-one 18->20 hemiketal	1985	3b,20b-Diacetoxy-9a-methyl-5a-pregnan-11-one
428	3b,18-Dihydroxy-5a-androstan-7-one	1633	3b,20b-Dibenzoyloxy-5-pregnene; 5-Pregnene-3b,20b-diol dibenzoate
1724	3b,19,20b/F-trihydroxy-5-pregnen-18-oic acid 3b-acetate 18,20b/F-lactone		
2101	3b,19-Diacetoxy-14,21-epoxy-5-hydroxy-5b,14b-cardanolide		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1237	3b,20b-Dibenzoyloxy-5a-pregnane; 5a-Pregnane-3b,20b-diol dibenzoate	2811	3b,5-Diacetoxy-5a-ergosta-7,22-dien-11-one
2021	3b,20b/F-Dihydroxy-16a-ethyl-5-pregnen-18-oic acid 3b-acetate 18,20b/F-lactone	2810	3b,5-Diacetoxy-5a-ergosta-9(11),22-dien-7-one (crude)
1916	3b,20x-Diacetoxypregna-5,17(20)-diene; Pregna-5,17(20)-diene-3b,20x-diol diacetate	1132	3b,5-Diacetoxy-5b,14b,17a-etiane-19,20-dioic acid 20-ethyl ester
2092	3b,21(a)-Diacetoxy-21-homo-5a,17a-pregnan-17b-ol	1133	3b,5-Diacetoxy-5b,14b,17a-etiane-19,20-dioic acid diethyl ester
1414	3b,21-Diacetoxy-14b-hydroxy-5b-pregnan-20-one	2446	3b,5-Diacetoxy-5b-cholestan-6-one; 3b,5-Dihydroxy-5b-cholestan-6-one diacetate
1416	3b,21-Diacetoxy-17a-pregn-5-en-17b-ol; 17a-Pregn-5-ene-3b,17b,21-triol 3,21-diacetate	2870	3b,5-Diacetoxy-5b-lumista-7,22-diene; 5b-Lumista-7,22-diene-3b,5-diol diacetate
1643	3b,21-Diacetoxy-17a-pregn-5-en-17b-ol; 17a-Pregn-5-ene-3b,17b,21-triol 3,21-diacetate	529	3b,5-Diacetoxy-6b-hydroxy-5a-androstan-17-one
1988	3b,21-Diacetoxy-20a-methyl-5-pregnen-20b-ol; 20a-Methyl-5-pregnene-3b,20b,21-triol 3,21-diacetate	647	3b,5-Dihydroxy-5a-androstane-6,17-dione 3-acetate; 3b-Acetoxy-5-hydroxy-5a-androstane-6,17-dione
1477	3b,21-Diacetoxy-5,6a-epoxy-17a-hydroxy-5b-pregnan-20-one	1126	3b,5-Dihydroxy-5b,14b,17a-etiane-19,20-dioic acid
1583	3b,21-Diacetoxy-5-hydroxy-20-oxo-5b,14b,17a-pregnano-19,8-lactone	1129	3b,5-Dihydroxy-5b,14b,17a-etiane-19,20-dioic acid 19-methyl ester 20-ethyl ester
1696	3b,21-Diacetoxy-5-pregnen-20-one	1128	3b,5-Dihydroxy-5b,14b,17a-etiane-19,20-dioic acid 20-ethyl ester
1805	3b,21-Diacetoxy-5-pregnene-11,20-dione	1127	3b,5-Dihydroxy-5b,14b,17a-etiane-19,20-dioic acid dimethyl ester
1806	3b,21-Diacetoxy-5a-pregn-16-ene-11,20-dione	2125	3b,5-Dihydroxy-5b,17a-carda-14,20(22)-dienolid-19-oic acid methyl ester
1366	3b,21-Diacetoxy-5a-pregnan-20-one	2124	3b,5-Dihydroxy-5b,17a-carda-14,20(22)-dienolid-19-oic acid; 14-Anhydro-17a-strophanthidinic acid
1486	3b,21-Diacetoxy-8b,19-epoxy-5-hydroxy-19-methoxy-5b,14b,17a-pregnan-20-one	1174	3b,5-Dihydroxy-5b,17a-eti-14-ene-19,20-dioic acid
1683	3b,21-Dihydroxy-5-pregnen-20-one 21-sulfate, potassium salt	1175	3b,5-Dihydroxy-5b,17a-eti-14-ene-19,20-dioic acid 20-ethyl ester
1686	3b,21-Dihydroxy-5-pregnen-20-one disulfate, triethylammonium salt	2443	3b,5-Dihydroxy-5b-cholestan-6-one
1680	3b,21-Dihydroxy-5-pregnen-20-one; 21-Hydroxypregnenolone	2444	3b,5-Dihydroxy-5b-cholestan-7-one
1793	3b,21-Dihydroxy-5-pregnene-11,20-dione	1475	3b,5-Dihydroxy-8b,19-epoxy-19-methoxy-5b-pregnan-20-one
2313	3b,22-Diacetoxy-11-oxobisnorchol-20-ene	1110	3b,5-Dihydroxy-8b,19-epoxy-5b-etianic acid
2132	3b,24-Diacetoxy-5-cholelene; 5-Cholelene-3b,24-diol diacetate	1111	3b,5-Dihydroxy-8b,19-epoxy-5b-etianic acid 3-acetate; 3b-Acetoxy-8b,19-epoxy-5-hydroxy-5b-etianic acid
2135	3b,24-Dihydroxy-5-chole-7-one	1445	3b,5-Dihydroxy-8b,19-epoxy-5b-pregnan-20-one
2578	3b,24x-Diacetoxy-5-cholesten-7-one; 3b,24x-Dihydroxy-5-cholesten-7-one diacetate	2783	3b,5a-Diacetoxy-5a-ergosta-7,22-diene; 5a-Ergosta-7,22-diene-3b,5a-diol diacetate
2540	3b,4a-Diacetoxy-5-cholestene; 5-Cholestene-3b,4a-diol diacetate	2824	3b,5a-Diacetoxy-5a-ergosta-7,9(11),22-triene; 5a-Ergosta-7,9(11),22-triene-3b,5a-diol diacetate
2403	3b,4b-(Isopropylidenedioxy)-5a-cholestan-5-ol; 5a-Cholestane-3b,4b,5-triol 3,4-acetonide	2781	3b,5a-Diacetoxy-5a-ergosta-7,9(11)-diene; 5a-Ergosta-7,9(11)-diene-3b,5a-diol diacetate
2346	3b,4b-Epoxy-5b-cholestane	2825	3b,5a-Diacetoxy-5a-ergosta-8,14,22-triene; 5a-Ergosta-8,14,22-triene-3b,5a-diol diacetate
1479	3b,5,14,19,21-Pentahydroxy-5b,14b-pregnan-20-one	2784	3b,5a-Diacetoxy-5a-ergosta-8,14-diene; 5a-Ergosta-8,14-diene-3b,5a-diol diacetate
1586	3b,5,14,19-Tetrahydroxy-20-oxo-5b,14b,17a-pregnan-21-oic acid	2316	3b,5b-Diacetoxy-9b,11b-epoxybisnorchol-7-en-22-al
1452	3b,5,14,19-Tetrahydroxy-5b,14,17a-pregnan-20-one	2895	3b,5b-Dihydroxy-7,22-lumistadien-6-one
2115	3b,5,14-Trihydroxy-19-oxo-5b,14b,17a-card-20(22)-enolide; 17a-Strophanthidin	1559	3b,5b-Tetrahydroaldosterone; 11b,18-Epoxy-3b,18,21-trihydroxy-5b-pregnan-20-one
2114	3b,5,14-Trihydroxy-19-oxo-5b,14b-card-20(22)-enolide; Strophanthidin	229	3b,6a,16a-Triacetoxo-5a-androstane; 5a-Androstane-3b,6a,16a-triol triacetate
1613	3b,5,14-Trihydroxy-20-oxo-5b,14b,17a-pregnane-19,21-dioic acid	2882	3b,6a-Diacetoxy-5b,9a-lumista-7,22-dien-5-ol; 5b,9a-Lumista-7,22-diene-3b,5,6a-triol 3,6-diacetate
1135	3b,5,14-Trihydroxy-5b,14b,17a-etiane-19,20-dioic acid	433	3b,6a-Dihydroxy-5a-androstan-11-one
2121	3b,5,14-Trihydroxy-5b,14b-card-20(22)-enolide-19-oic acid; Strophanthidinic acid	451	3b,6a-Dihydroxy-5a-androstan-17-one
514	3b,5,16b-Trihydroxy-5a-androstan-11-one	1461	3b,6b-Diacetoxy-5,17a-dihydroxy-5a-pregnan-20-one
1418	3b,5,19-Triacetoxo-5b-pregnan-20-one	2387	3b,6b-Diacetoxy-5-fluoro-5a-cholestane; 5-Fluoro-5a-cholestane-3b,6b-diol diacetate
1385	3b,5,19-Trihydroxy-5b-pregnan-20-one	530	3b,6b-Diacetoxy-5-hydroxy-5a-androstan-17-one
1581	3b,5,21-Trihydroxy-20-oxo-5b,14b,17a-pregnano-19,8-lactone	2404	3b,6b-Diacetoxy-5a-cholestan-5-ol; 5a-Cholestane-3b,5,6b-triol 3,6-diacetate
1580	3b,5,21-Trihydroxy-20-oxo-5b-pregnano-19,8-lactone	2380	3b,6b-Diacetoxy-5a-cholestane; 5a-Cholestane-3b,6b-diol diacetate
519	3b,5,6a-Trihydroxy-5a-androstan-17-one	1103	3b,6b-Diacetoxy-5a-etianic acid methyl ester
1451	3b,5,6b,17a-Tetrahydroxy-5a-pregnan-20-one	2851	3b,6b-Diacetoxy-5b-lumist-22-ene
543	3b,5,6b-Triacetoxo-5a-androstan-17-one		
1137	3b,5,8-Trihydroxy-5b,17a-etiane-19,20-dioic acid 19,8-lactone		
1136	3b,5,8-Trihydroxy-5b-etiane-19,20-dioic acid 19,8-lactone		
2765	3b,5-Diacetoxy-5a-ergost-22-en-7-one		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2883	3b,6b-Diacetoxy-5b-lumista-7,22-dien-5-ol; 5b-Lumista-7,22-diene-3b,5,6b-triol 3,6-diacetate	1115	3b-Acetoxy-12-oxo-5b-etianic acid methyl ester
2889	3b,6b-Dibenzoyloxy-5b-lumista-7,22-dien-5-ol; 5b-Lumista-7,22-diene-3b,5,6b-triol 3,6-dibenzoate	2723	3b-Acetoxy-12a-bromo-5a-ergostan-11-one
2874	3b,6b-Epoxy-9a-lumista-7,22-dien-5a-ol	1120	3b-Acetoxy-14b-hydroxy-12-oxo-5b,14b-etianic acid methyl ester
1707	3b,7a,11a-Trihydroxy-5-pregnen-20-one	1116	3b-Acetoxy-15-oxo-5b,14a-etianic acid methyl ester
521	3b,7a,14a-Trihydroxy-5a-androstan-17-one	1117	3b-Acetoxy-15-oxo-5b,14b-etianic acid methyl ester
539	3b,7a,16b-Triacetoxy-5a-androstan-11-one	2025	3b-Acetoxy-16-methyl-5a-pregna-9(11),16-dien-20-one
230	3b,7a,16b-Triacetoxy-5a-androstane; 5a-Androstane-3b,7a,16b-triol triacetate	1818	3b-Acetoxy-16a,17a-epoxy-5-pregnene-11,20-dione
515	3b,7a,16b-Trihydroxy-5a-androstan-11-one	1998	3b-Acetoxy-16b-methyl-5a-pregna-9(11)-en-20-one
231	3b,7a,17b-Triacetoxy-5a-androstane; 5a-Androstane-3b,7a,17b-triol triacetate	2095	3b-Acetoxy-17,21(a)-epoxy-21-homo-5a-pregna-20-ene
509	3b,7a,17b-Trihydroxy-5a-androstan-6-one	1710	3b-Acetoxy-17a,21-dihydroxy-5-pregnen-20-one
434	3b,7a-Dihydroxy-5a-androstan-11-one	2000	3b-Acetoxy-17a,21-dimethyl-5-pregnen-20-one (crude)
1335	3b,7a-Dihydroxy-5a-pregnan-20-one (crude)	1687	3b-Acetoxy-17a-hydroxy-5-pregnen-20-one
1419	3b,7b,11-Triacetoxy-5a-pregnan-20-one (crude)	1691	3b-Acetoxy-17a-hydroxy-5-pregnen-20-one 17-sulfate, triethylammonium salt
1716	3b,7b,11a-Triacetoxy-5-pregnen-20-one	1356	3b-Acetoxy-17a-hydroxy-5a-pregnan-20-one
544	3b,7b,11a-Triacetoxy-5a-androstan-17-one	1970	3b-Acetoxy-17a-pregna-5-en-20-yn-17b-ol; 17a-Pregna-5-en-20-yne-3b,17b-diol 3-acetate
1706	3b,7b,11a-Trihydroxy-5-pregnen-20-one	1912	3b-Acetoxy-17a-pregna-5,20-dien-17b-ol; 17a-Pregna-5,20-diene-3b,17b-diol 3-acetate
520	3b,7b,11a-Trihydroxy-5a-androstan-17-one	1972	3b-Acetoxy-17b-benzoyloxy-17a-pregna-5-en-20-yne; 17a-Pregna-5-en-20-yne-3b,17b-diol 3-acetate 17-benzoate
1388	3b,7b,11a-Trihydroxy-5a-pregnan-20-one	1156	3b-Acetoxy-17b-hydroxy-17a-eti-5-enic acid methyl ester
1420	3b,7b,12b-Triacetoxy-5a-pregnan-20-one	1689	3b-Acetoxy-18-hydroxy-5-pregnen-20-one 18->20 hemiketal
1390	3b,7b,12b-Trihydroxy-5a-pregnan-20-one (crude)	840	3b-Acetoxy-19-hydroxy-5-androsten-17-one 19-sulfate, pyridinium salt
1391	3b,7b,15a-Trihydroxy-5a-pregnan-20-one	2935	3b-Acetoxy-19-nor-5,7,9-ergostatrien-11-one
516	3b,7b,16b-Trihydroxy-5a-androstan-11-one	838	3b-Acetoxy-1a-hydroxy-5-androsten-17-one; 1a,3b-Dihydroxy-5-androsten-17-one 3-acetate
540	3b,7b,17a-Triacetoxy-5a-androstan-11-one	1404	3b-Acetoxy-1b,11a-dihydroxy-5a-pregnan-20-one
2543	3b,7b-Dibenzoyloxy-5-cholestene; 5-Cholestene-3b,7b-diol dibenzoate	1205	3b-Acetoxy-20,20-difluoro-5a-pregnane
825	3b,7b-Dihydroxy-5-androsten-17-one	1513	3b-Acetoxy-20,20-ethylenedioxy-5a-pregnan-2-one
1043	3b,7b-Dihydroxy-A-nor(5->6b)abeo-androstan-5-one	1377	3b-Acetoxy-20,21-dibromo-17b-hydroxy-5a,17a-pregnan-20-ene
1101	3b,7x-Diacetoxy-5a-etianic acid	1928	3b-Acetoxy-20-hydroximinopregna-5,16-diene; 16a-Dehydropregnenolone acetate oxime
1381	3b,9a,20b-Trihydroxy-5a-pregnan-11-one	1267	3b-Acetoxy-20a,21-epoxy-5a,17a-pregnan-17b-ol; 20a,21-Epoxy-5a,17a-pregnane-3b,17b-diol 3-acetate
2894	3b,9b-Epoxy-5a-lumista-7,22-dien-6-one	1268	3b-Acetoxy-20b,21-epoxy-5a,17a-pregnan-17b-ol; 20b,21-Epoxy-5a,17a-pregnane-3b,17b-diol 3-acetate
2875	3b,9b-Epoxy-5a-lumista-7,22-dien-6a-ol	1448	3b-Acetoxy-21-diazo-8b,19-epoxy-5-hydroxy-5b-pregnan-20-one
913	3b-(3,5-Dinitrobenzoyloxy)-5,7-androstadien-17-one	2094	3b-Acetoxy-21-homo-5a,17a-pregna-21-en-17b-ol
372	3b-(p-Fluorobenzenesulfonyloxy)-5a-androstan-17-one; 3b-Hydroxy-5a-androstan-17-one p-fluorobenzenesulfonate	2280	3b-Acetoxy-23,24-bisnor-5a-cholane; 23,24-Bisnor-5a-cholan-3-ol acetate; 3b-Acetoxy-20-methyl-5a-pregnane
2448	3b-(p-Toluenesulfonyloxy)-5-hydroxy-5a-cholestan-6-one	2748	3b-Acetoxy-23,24-dibromo-5a-ergost-6-ene 5,8a-maleic anhydride adduct (crude)
371	3b-(p-Toluenesulfonyloxy)-5a-androstan-17-one; 3b-Hydroxy-5a-androstan-17-one tosylate	947	3b-Acetoxy-4,4-dimethyl-5a-androstan-2-one; 4,4-Dimethyl-5a-androstan-3b-ol-2-one acetate
632	3b-(p-Toluenesulfonyloxy)-5a-androstane-11,17-dione	2574	3b-Acetoxy-4-cholesten-6-one; 3b-Hydroxy-4-cholesten-6-one acetate
2385	3b-(p-Toluenesulfonyloxy)-5a-cholestan-1a-ol; 5a-Cholestane-1a,3b-diol 3-tosylate	2735	3b-Acetoxy-4-ergostene; 4-Ergosten-3b-ol acetate
2386	3b-(p-Toluenesulfonyloxy)-5a-cholestan-5-ol; 5a-Cholestane-3b,5-diol 3-tosylate	2117	3b-Acetoxy-5,14-dihydroxy-19-oxo-5b,14b,17a-card-20(22)-enolide; 17a-Strophanthidin 3-acetate
2341	3b-(p-Toluenesulfonyloxy)-5a-cholestan-5-ol; 5a-Cholestan-3b-ol tosylate	899	3b-Acetoxy-5,16-androstadiene
1512	3b-(p-Toluenesulfonyloxy)-5a-pregnane-11,20-dione	2397	3b-Acetoxy-5,6b-epoxy-5b-cholestan-5,6b-Epoxy-5b-cholestan-3b-ol acetate
2571	3b-(Tetrahydropyranyloxy)-5-cholesten-7-one; 3b-Hydroxy-5-cholesten-7-one tetrahydropyranyl ether	2873	3b-Acetoxy-5,6b-epoxy-5b-lumista-7,22-diene; 5,6b-Epoxy-5b-lumista-7,22-dien-3b-ol acetate
2324	3b-Acetamido-5a-cholestane	2766	3b-Acetoxy-5,8a-peroxido-5a-ergost-9(11)-en-12-one
2616	3b-Acetoxy-11,11-dibromo-5a,25(R)-spirostan-12-one; Hecogenin acetate 11,11-dibromide	2712	3b-Acetoxy-5,8a-peroxido-5a-ergostane-9a,11a-diol; 5,8a-Peroxido-5a-ergostane-3b,9,11a-triol 3-acetate
2212	3b-Acetoxy-11-oxo-5b-cholan-24-oic acid methyl ester	2726	3b-Acetoxy-5,9a,11b-trihydroxy-5a-ergostan-7-one
2705	3b-Acetoxy-11a,12a-epoxy-5a-ergostane; 11a,12a-Epoxy-5a-ergostan-3b-ol acetate	802	3b-Acetoxy-5-androsten-17-one semicarbazone
2183	3b-Acetoxy-11b,12b-epoxy-5b-cholan-24-oic acid methyl ester		
2224	3b-Acetoxy-11b-hydroxy-12-oxo-5b-cholan-24-oic acid methyl ester		
2303	3b-Acetoxy-12-oxo-5b-bisnorcholan-22-oic acid methyl ester		
2267	3b-Acetoxy-12-oxo-5b-chole-9(11)-en-24-oic acid methyl ester		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
884	3b-Acetoxy-5-androstene-11,17-dione	628	3b-Acetoxy-5a-androstane-11,17-dione
2393	3b-Acetoxy-5-bromo-5a-cholestan-6b-ol; 5-Bromo-5a-cholestan-3b,6b-diol 3-acetate	178	3b-Acetoxy-5a-androstane-6a-ol
2255	3b-Acetoxy-5-cholesten-24-oic acid	2307	3b-Acetoxy-5a-bisnorchol-20-ene; 23,24-Bisnor-5a-chole-20-en-3b-ol acetate
2576	3b-Acetoxy-5-cholesten-22-one; 3b-Hydroxy-5-cholesten-22-one acetate	2282	3b-Acetoxy-5a-bisnorchol-20-ol; 23,24-Bisnor-5a-chole-3b,20-diol 3-acetate
2573	3b-Acetoxy-5-cholesten-7-one; 3b-Hydroxy-5-cholesten-7-one acetate	1290	3b-Acetoxy-5a-bromo-6b,19:18,20b-diepoxy-18-methoxypregnane
2736	3b-Acetoxy-5-ergostene; 5-Ergosten-3b-ol acetate	2910	3b-Acetoxy-5a-campest-8(14)-ene; 5a-Campest-8(14)-en-3b-ol acetate
1146	3b-Acetoxy-5-etienic acid	2476	3b-Acetoxy-5a-cholest-1-ene; 5a-Cholest-1-en-3b-ol acetate
1147	3b-Acetoxy-5-etienic acid methyl ester	2537	3b-Acetoxy-5a-cholest-6-en-5-ol; 5a-Cholest-6-ene-3b,5-diol 3-acetate
2438	3b-Acetoxy-5-fluoro-5a-cholestan-6-one; 5-Fluoro-3b-hydroxy-5a-cholestan-6-one acetate	2575	3b-Acetoxy-5a-cholest-8-en-7-one; 3b-Hydroxy-5a-cholest-8-en-7-one acetate
2607	3b-Acetoxy-5-hydroxy-5a-cholesta-7,22-dien-6-one	2596	3b-Acetoxy-5a-cholesta-6,8(9)-diene; 5a-Cholesta-6,8(9)-dien-3b-ol acetate
2445	3b-Acetoxy-5-hydroxy-5a-cholestan-6-one	2597	3b-Acetoxy-5a-cholesta-7,14-diene; 5a-Cholesta-7,14-dien-3b-ol acetate
2809	3b-Acetoxy-5-hydroxy-5a-ergosta-7,22-dien-11-one	2432	3b-Acetoxy-5a-cholestan-2-one; 3b-Hydroxy-5a-cholestan-2-one acetate
1130	3b-Acetoxy-5-hydroxy-5b,14b,17a-etiane-19,20-dioic acid 20-ethyl ester	2436	3b-Acetoxy-5a-cholestan-7-one; 3b-Hydroxy-5a-cholestan-7-one acetate
1131	3b-Acetoxy-5-hydroxy-5b,14b,17a-etiane-19,20-dioic acid diethyl ester	2374	3b-Acetoxy-5a-cholestan-7a-ol; 5a-Cholestane-3b,7a-diol 3-acetate
1628	3b-Acetoxy-5-pregnen-17a-ol; 5-Pregnene-3b,17a-diol 3-acetate	2375	3b-Acetoxy-5a-cholestan-7b-ol; 5a-Cholestane-3b,7b-diol 3-acetate
1662	3b-Acetoxy-5-pregnen-20-one	2401	3b-Acetoxy-5a-cholestan-5,6b-diol; 5a-Cholestane-3b,5,6b-triol 3-acetate
1621	3b-Acetoxy-5-pregnene; 5-Pregnen-3b-ol acetate	2740	3b-Acetoxy-5a-ergost-11-ene; 5a-Ergost-11-en-3b-ol acetate
2818	3b-Acetoxy-5a,14b-ergosta-7,9(11),22-triene; 5a,14b-Ergosta-7,9(11),22-trien-3b-ol acetate	2753	3b-Acetoxy-5a-ergost-22-en-11b-ol; 5a-Ergost-22-ene-3b,11b-diol 3-acetate
2806	3b-Acetoxy-5a,14b-ergosta-8,22-dien-11-one	2758	3b-Acetoxy-5a-ergost-22-ene-5,8a-diol; 5a-Ergost-22-ene-3b,5,8a-triol 3-acetate
1914	3b-Acetoxy-5a,17a-pregn-20-yn-17b-ol; 5a,17a-Pregn-20-yne-3b,17b-diol 3-acetate	2768	3b-Acetoxy-5a-ergost-22-ene-7,11-dione
1220	3b-Acetoxy-5a,17a-pregnan-17b-ol; 5a,17a-Pregnane-3b,17b-diol 3-acetate	2742	3b-Acetoxy-5a-ergost-22-ene; 5a-Ergost-22-en-3b-ol acetate
2612	3b-Acetoxy-5a,25(R)-spirostan-12-ol; 5a,25(R)-Spirostan-3b,12b-diol 3-acetate	2737	3b-Acetoxy-5a-ergost-6-ene; 5a-Ergost-6-en-3b-ol acetate
505	3b-Acetoxy-5a,6a-epoxyandrostan-17-one	2764	3b-Acetoxy-5a-ergost-8(14)-en-11-one
806	3b-Acetoxy-5a,6b-dichloro-17,17-ethylenedioxyandrost-14-ene	2752	3b-Acetoxy-5a-ergost-8(14)-en-11b-ol; 5a-Ergost-8(14)-ene-3b,11b-diol 3-acetate
2720	3b-Acetoxy-5a,8a-ergostan-11-one	2751	3b-Acetoxy-5a-ergost-8(14)-en-5-ol; 5a-Ergost-8(14)-ene-3b,5-diol 3-acetate
2695	3b-Acetoxy-5a,8a-ergostan-11b-ol; 5a,8b-Ergostane-3b,11b-diol 3-acetate	2747	3b-Acetoxy-5a-ergost-8(14)-ene 7a,15a-maleic anhydride adduct
2711	3b-Acetoxy-5a,8a-peroxido-22x,23x-epoxyergostane; 5a,8a-Peroxido-22x,23x-epoxyergostan-3b-ol acetate	2738	3b-Acetoxy-5a-ergost-8(14)-ene; 5a-Ergost-8(14)-en-3b-ol acetate
2755	3b-Acetoxy-5a,8a-peroxido-5a-ergost-22-ene; 5a,8a-Peroxido-5a-ergost-22-en-3b-ol acetate	2757	3b-Acetoxy-5a-ergost-9(11)-ene-5,8a-diol; 5a-Ergost-9(11)-ene-3b,5,8a-triol 3-acetate (crude)
2314	3b-Acetoxy-5a,8a-peroxidobisnorchol-9(11)-en-22-ol	2739	3b-Acetoxy-5a-ergost-9(11)-ene; 5a-Ergost-9(11)-en-3b-ol acetate
2707	3b-Acetoxy-5a,8a-peroxidoergostane; 5a,8a-Peroxidoergostan-3b-ol acetate	2792	3b-Acetoxy-5a-ergosta-6,22-diene-5,8a-diol; 5a-Ergosta-6,22-diene-3b,5,8a-triol 3-acetate (crude)
2850	3b-Acetoxy-5a,9a-lumist-8(14)-ene	2780	3b-Acetoxy-5a-ergosta-7,22-dien-11b-ol; 5a-Ergosta-7,22-diene-3b,11b-diol 3-acetate
2741	3b-Acetoxy-5a,9b-ergost-11-ene; 5a,9b-Ergost-11-en-3b-ol acetate	2778	3b-Acetoxy-5a-ergosta-7,22-dien-5-ol; 5a-Ergosta-7,22-diene-3b,5-diol 3-acetate
2807	3b-Acetoxy-5a,9b-ergosta-7,22-dien-11-one	2793	3b-Acetoxy-5a-ergosta-7,22-diene-5,11b-diol; 5a-Ergosta-7,22-diene-3b,5,11b-triol 3-acetate (crude)
2721	3b-Acetoxy-5a,9b-ergostan-11-one	2791	3b-Acetoxy-5a-ergosta-7,22-diene-5,6a-diol; 5a-Ergosta-7,22-diene-3b,5,6a-triol 3-acetate
2688	3b-Acetoxy-5a,9b-ergostane; 5a,9b-Ergostan-3b-ol acetate	2822	3b-Acetoxy-5a-ergosta-7,9(11),22-trien-5-ol; 5a-Ergosta-7,9(11),22-triene-3b,5-diol 3-acetate
710	3b-Acetoxy-5a-androst-16-ene	2775	3b-Acetoxy-5a-ergosta-7,9(11)-dien-5-ol; 5a-Ergosta-7,9(11)-diene-3b,5-diol 3-acetate
709	3b-Acetoxy-5a-androst-9(11)-ene		
350	3b-Acetoxy-5a-androstan-12-one; 5a-Androstan-3b-ol-12-one acetate		
351	3b-Acetoxy-5a-androstan-16-one; 5a-Androstan-3b-ol-16-one acetate		
353	3b-Acetoxy-5a-androstan-17-one; 5a-Androstan-3b-ol-17-one acetate		
343	3b-Acetoxy-5a-androstan-2-one; 5a-Androstan-3b-ol-2-one acetate		
349	3b-Acetoxy-5a-androstan-7-one; 5a-Androstan-3b-ol-7-one acetate		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2823	3b-Acetoxy-5a-ergosta-8,14,22-trien-5-ol; 5a-Ergosta-8,14,22-triene-3b,5-diol 3-acetate	716	3b-Acetoxy-6-fluoro-5-androstene
2819	3b-Acetoxy-5a-ergosta-8,14,22-triene; 5a-Ergosta-8,14,22-trien-3b-ol acetate	2796	3b-Acetoxy-6a-benzoyloxy-5a,9b-ergosta-7,22-dien-5-ol; 5a,9b-Ergosta-7,22-dien-3b,5,6a-triol 3-acetate 6-benzoate
2805	3b-Acetoxy-5a-ergosta-8,22-dien-11-one	2439	3b-Acetoxy-6a-bromo-5a-cholestan-7-one; 6a-Bromo-3b-hydroxy-5a-cholestan-7-one acetate
2804	3b-Acetoxy-5a-ergosta-8,22-dien-7-one	2440	3b-Acetoxy-6b-bromo-5a-cholestan-7-one; 6b-Bromo-3b-hydroxy-5a-cholestan-7-one acetate
2719	3b-Acetoxy-5a-ergostan-11-one	2388	3b-Acetoxy-6b-chloro-5a-cholestan-5-ol; 6b-Chloro-5a-cholestane-3b,5-diol 3-acetate
2694	3b-Acetoxy-5a-ergostan-11b-ol; 5a-Ergostane-3b,11b-diol 3-acetate	491	3b-Acetoxy-6b-fluoro-5-hydroxy-5a-androstan-17-one
2687	3b-Acetoxy-5a-ergostane; 5a-Ergostan-3b-ol acetate	2402	3b-Acetoxy-6b-methoxy-5a-cholestan-5-ol; 6b-Methoxy-5a-cholestane-3b,5-diol 3-acetate
1148	3b-Acetoxy-5a-eti-16-enic acid methyl ester	2265	3b-Acetoxy-7-oxo-5-cholen-24-oic acid methyl ester
1089	3b-Acetoxy-5a-etianic acid methyl ester	1167	3b-Acetoxy-7-oxoeti-5-enic acid methyl ester
2318	3b-Acetoxy-5a-hydroxy-23,24-bisnorchola-7,9(11)-dien-22-oic acid methyl ester	2542	3b-Acetoxy-7a-benzoyloxy-5-cholestene; 5-Cholestene-3b,7a-diol 3-acetate 7-benzoate
466	3b-Acetoxy-5a-hydroxy-5a-androstan-17-one	2759	3b-Acetoxy-7b,8b:9b,11b-diepoxy-5a-ergost-22-en-5-ol; 7b,8b:9b,11b-Diepoxy-5a-ergost-22-ene-3b,5-diol 3-acetate
2848	3b-Acetoxy-5a-lumist-7-ene	2713	3b-Acetoxy-7b,8b:9b,11b-diepoxy-5b-ergostan-5b-ol; 7b,8b:9b,11b-Diepoxy-5b-ergostane-3b,5b-diol 3-acetate
2953	3b-Acetoxy-5a-poriferastane; 5a-Poriferastan-3b-ol acetate	2530	3b-Acetoxy-7b-anilino-5-cholestene; 7b-Anilinocholesteryl acetate
1664	3b-Acetoxy-5a-pregn-16-en-20-one	2714	3b-Acetoxy-7x,8x,9x,11x-diepoxy-5a-ergostan-5-ol; 7x,8x,9x,11x-Diepoxy-5a-ergostane-3b,5-diol 3-acetate
1762	3b-Acetoxy-5a-pregn-16-ene-11,20-dione 20-oxime	2450	3b-Acetoxy-8a,14a-epoxy-5a-cholestan-7-one; 8a,14a-Epoxy-3b-hydroxy-5a-cholestan-7-one acetate
1763	3b-Acetoxy-5a-pregn-16-ene-11,20-dione 3-tosylate	2701	3b-Acetoxy-8b,14b-epoxy-5a-ergostane; 8b,14b-Epoxy-5a-ergostan-3b-ol acetate
1756	3b-Acetoxy-5a-pregn-16-ene-12,20-dione	1112	3b-Acetoxy-8b,19-epoxy-5-hydroxy-5b-etianic acid ethyl ester
1622	3b-Acetoxy-5a-pregn-17(20)-ene; 5a-Pregn-17(20)-en-3b-ol acetate	2689	3b-Acetoxy-9a,11a-dichloro-5a-ergostane; 9a,11a-Dichloro-5a-ergostan-3b-ol acetate
1663	3b-Acetoxy-5a-pregn-9(11)-en-20-one	2725	3b-Acetoxy-9a,11a-epoxy-5-hydroxy-5a-ergostan-7-one
1510	3b-Acetoxy-5a-pregnane-11,20-dione	2703	3b-Acetoxy-9a,11a-epoxy-5a,8a-ergostane; 9a,11a-Epoxy-5a,8a-ergostan-3b-ol acetate
2960	3b-Acetoxy-5a-stigmast-7-ene; 5a-Stigmast-7-en-3b-ol acetate	2797	3b-Acetoxy-9a,11a-epoxy-5a-ergosta-7,22-dien-5-ol; 9a,11a-Epoxy-5a-ergosta-7,22-diene-3b,5-diol 3-acetate (crude)
2974	3b-Acetoxy-5a-stigmasta-7,22(E)-diene	2702	3b-Acetoxy-9a,11a-epoxy-5a-ergostane; 9a,11a-Epoxy-5a-ergostan-3b-ol acetate
2952	3b-Acetoxy-5a-stigmastane; Stigmastanol acetate	2704	3b-Acetoxy-9a,11b-epoxy-5a-ergostane; 9a,11b-Epoxy-5a-ergostan-3b-ol acetate
2311	3b-Acetoxy-5b,8-peroxido-23,24-bisnorchol-9(11)-en-22-oic acid methyl ester	1546	3b-Acetoxy-9a,21-dibromo-17a-hydroxy-5a-pregnane-11,20-dione
2709	3b-Acetoxy-5b,8b:9a,11a-diepoxy-5b-ergostane; 5b,8b:9a,11a-Diepoxy-5b-ergostan-3b-ol acetate	2722	3b-Acetoxy-9a-bromo-5a-ergostan-11-one
2710	3b-Acetoxy-5b,8b:9b,11b-diepoxy-5b-ergostane; 5b,8b:9b,11b-Diepoxy-5b-ergostan-3b-ol acetate	2700	3b-Acetoxy-9a-bromo-5a-ergostan-11b-ol; 9a-Bromo-5a-ergostane-3b,11b-diol 3-acetate (crude)
629	3b-Acetoxy-5b-androstane-11,17-dione	467	3b-Acetoxy-9a-hydroxy-5a-androstan-17-one; 3b,9a-Dihydroxy-5a-androstan-17-one 3-acetate
2258	3b-Acetoxy-5b-chol-11-en-24-oic acid methyl ester	465	3b-Acetoxy-9a-hydroxy-5a-androstan-7-one; 3b,9a-Dihydroxy-5a-androstan-7-one 3-acetate
2153	3b-Acetoxy-5b-cholan-24-oic acid methyl ester	2724	3b-Acetoxy-9a-hydroxy-5a-ergostan-11-one
2605	3b-Acetoxy-5b-cholesta-8(14),15-dien-7-one; 3b-Hydroxy-5b-cholesta-8(14),15-dien-7-one acetate	2857	3b-Acetoxy-9a-lumista-5,22-diene
2817	3b-Acetoxy-5b-ergosta-7,9(11),22-triene; 5b-Ergosta-7,9(11),22-trien-3b-ol acetate	946	3b-Acetoxy-9a-methyl-5a-androstan-11-one; 9a-Methyl-5a-androstan-3b-ol-11-one acetate
2776	3b-Acetoxy-5b-ergosta-7,9(11)-dien-5-ol; 5b-Ergosta-7,9(11)-diene-3b,5-diol 3-acetate	930	3b-Acetoxy-9a-methyl-5a-androstane
1150	3b-Acetoxy-5b-eti-9(11)-enic acid methyl ester	2915	3b-Acetoxy-9a-methyl-5a-ergost-8(14)en-11-one
1090	3b-Acetoxy-5b-etianic acid methyl ester	2920	3b-Acetoxy-9a-methyl-5a-ergosta-7,22-dien-11-one
2849	3b-Acetoxy-5b-lumist-7-ene	2798	3b-Acetoxy-9b,11b-epoxy-5a-ergosta-7,22-dien-5-ol; 9b,11b-Epoxy-5a-ergosta-7,22-diene-3b,5-diol 3-acetate
2864	3b-Acetoxy-5b-lumista-6,22-diene 5b,8b-maleic anhydride adduct	1107	3b-Acetoxy-9b,11b-epoxy-5b-etianic acid methyl ester
2869	3b-Acetoxy-5b-lumista-7,22-dien-5-ol; 5b-Lumista-7,22-diene-3b,5-diol 3-acetate	1059	3b-Acetoxy-D-homo-5a-androstan-17(a)-one
2893	3b-Acetoxy-5b-lumista-7,22-dien-6-one	2308	3b-Acetoxybisnor-5-cholen-22-oic acid methyl ester
2858	3b-Acetoxy-5b-lumista-7,22-diene	2309	3b-Acetoxybisnorchol-11-en-22-oic acid methyl ester
2881	3b-Acetoxy-5b-lumista-7,22-diene-5,6b-diol; 5b-Lumista-7,22-diene-3b,5,6b-triol 3-acetate	2909	3b-Acetoxycampest-5-ene; Campest-5-en-3b-ol acetate
2900	3b-Acetoxy-5b-lumista-8,14,22-triene; 5b-Lumista-8,14,22-trien-3b-ol acetate		
2845	3b-Acetoxy-5b-lumistan-6-one		
1357	3b-Acetoxy-5b-pregnan-17a-ol-20-one; 3b,17a-Dihydroxy-5b-pregnan-20-one 3-acetate		
1316	3b-Acetoxy-5b-pregnan-20-one		
2459	3b-Acetoxy-5x-cholestane-6,7-dione; 3b-Hydroxy-5x-cholestane-6,7-dione acetate		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2912	3b-Acetoxycampesta-5,7-diene; Campesta-5,7-dien-3b-ol acetate	1655	3b-Chloro-5-pregnen-20-one
2594	3b-Acetoxycholesta-5,22(E)-diene; Cholesta-5,22(E)-dien-3b-ol acetate	288	3b-Chloro-5a-androstan-17-one
2595	3b-Acetoxycholesta-5,24-diene; Cholesta-5,24-dien-3b-ol acetate	2414	3b-Chloro-5a-cholestan-6-one
2833	3b-Acetoxyergosta-5,7,9(11),22-tetraene; Ergosta-5,7,9(11),22-tetraen-3b-ol acetate	853	3b-Chloro-7b,11a-dihydroxy-5-androsten-17-one
2907	3b-Acetoxyllumista-5,7,9(11),22-tetraene; Lumista-5,7,9(11),22-tetraen-3b-ol acetate	786	3b-Chloroacetoxy-5-androsten-17-one; 5-Androsten-3b-ol-17-one chloroacetate
2982	3b-Acetoxyilupan-2-one	713	3b-Chloroacetoxy-5-androstene
2636	3b-Acetoxy-nor-5-cholesten-24-one	133	3b-Chloroacetoxy-5a-androstane; 5a-Androstan-3b-ol chloroacetate
2959	3b-Acetoxy-poriferast-5-ene; Poriferast-5-en-3b-ol acetate	714	3b-Dichloroacetoxy-5-androstene
2972	3b-Acetoxy-poriferasta-5,22(E)-diene; Poriferastrerol acetate	385	3b-Ethoxy-17,17-ethylenedioxy-5a-androstane; 3b-Ethoxy-5a-androstan-17-one ethylene ketal
1926	3b-Acetoxy-pregna-5,16-dien-20-one	936	3b-Ethyl-5a-androstan-1-one
1911	3b-Acetoxy-pregna-5,17(20)-dien-21-ol; Pregna-5,17(20)-diene-3b,21-diol 3-acetate	805	3b-Fluoro-11a-hydroxy-5-androsten-17-one; 3b-Fluoro-5-androsten-11a-ol-17-one
2971	3b-Acetoxy-stigmasta-5,22(E)-diene; Stigmasteryl acetate	389	3b-Fluoro-11a-hydroxy-5a-androstan-17-one
2973	3b-Acetoxy-stigmasta-5,24(28)(E)-diene; Fucosteryl acetate	765	3b-Fluoro-5-androsten-17-one
2976	3b-Acetoxy-stigmasta-5,7,22(E)-triene; Stigmasta-5,7,22(E)-trien-3b-ol acetate	897	3b-Fluoro-5-androstene-7,11,17-trione
2970	3b-Acetoxy-stigmasta-5,7-diene; Stigmasta-5,7-dien-3b-ol acetate	279	3b-Fluoro-5a-androstan-17-one
2478	3b-Acetylthio-5-cholestene	590	3b-Fluoro-5a-androstane-11,17-dione
2706	3b-Benzoyloxy-14a,15a-epoxy-5a-ergostane; 14a,15a-Epoxy-5a-ergostan-3b-ol benzoate	488	3b-Fluoro-6a,11a-dihydroxy-5a-androstan-17-one
2743	3b-Benzoyloxy-4-ergostene; 4-Ergosten-3b-ol benzoate	804	3b-Fluoro-7b-hydroxy-5-androsten-17-one; 3b-Fluoro-5-androsten-7b-ol-17-one
2122	3b-Benzoyloxy-5,14-dihydroxy-5b,14b,17a-card-20(22)-enolide-19-oic acid Me ester; 17a-Strophanthidinic acid Me ester 3-benzoate	2958	3b-Formyloxystigmast-5-ene; b-Sitosteryl formate
2577	3b-Benzoyloxy-5-cholesten-7-one oxime; 7-Hydroximino-5-cholesten-3b-ol 3-benzoate	2207	3b-Hydroxy-11-oxo-5b-cholan-24-oic acid methyl ester
2744	3b-Benzoyloxy-5-ergostene; 5-Ergosten-3b-ol benzoate	1994	3b-Hydroxy-16-methyl-5a-pregn-16-en-20-one
2315	3b-Benzoyloxy-5a,8a-peroxidobisnorchol-9(11)-en-22-al dimethyl acetal	1991	3b-Hydroxy-16a-methyl-5-pregnen-20-one; 16a-Methylpregnenolone
362	3b-Benzoyloxy-5a-androstan-17-one; 5a-Androstan-3b-ol-17-one benzoate	1002	3b-Hydroxy-16b-isopropyl-5-androsten-17-one
2441	3b-Benzoyloxy-5a-cholestan-6-one oxime; 6-Hydroximino-5a-cholestan-3b-ol 3-benzoate	1982	3b-Hydroxy-16b-methyl-5a-pregnan-20-one
2745	3b-Benzoyloxy-5a-ergost-14-ene; 5a-Ergost-14-en-3b-ol benzoate	2061	3b-Hydroxy-19-nor-14b,17a-pregn-5(10)-en-20-one
2820	3b-Benzoyloxy-5a-ergosta-7,14,22-triene; 5a-Ergosta-7,14,22-trien-3b-ol benzoate	2062	3b-Hydroxy-19-nor-14b,17a-pregn-5(10)-en-20-one 3',5'-dinitrobenzoate
2319	3b-Benzoyloxy-5a-hydroxy-23,24-bisnorchola-7,9(11)-dien-22-oic acid methyl ester	1711	3b-Hydroxy-1b,11a-(isopropylidenedioxy)-5-pregnen-20-one
1320	3b-Benzoyloxy-5a-pregnan-20-one	2635	3b-Hydroxy-27-nor-5-cholesten-24-one
2407	3b-Benzoyloxy-5b,6b-epoxy-5b-cholestan-7b-ol; 5b,6b-Epoxy-5b-cholestane-3b,7b-diol 3-benzoate	944	3b-Hydroxy-4,4,14a-trimethyl-5a-androstan-17-one
2312	3b-Benzoyloxy-5b,8-peroxido-23,24-bisnorchol-9(11)-en-22-oic acid methyl ester	2564	3b-Hydroxy-4-cholesten-6-one
2623	3b-Benzoyloxy-7a-methyl-5-cholestene; 7a-Methylcholesteryl benzoate	1176	3b-Hydroxy-5,16-etiadienic acid
2921	3b-Benzoyloxy-9a-methyl-5a-ergosta-7,22-dien-11-one	781	3b-Hydroxy-5-androsten-17-one 3b-sulfate, triethylammonium salt
2859	3b-Benzoyloxy-5a-lumista-7,22-diene	800	3b-Hydroxy-5-androsten-17-one b-D-glucuronide; 5-Androsten-3b-ol-17-one b-D-glucuronide
787	3b-Bromoacetoxy-5-androsten-17-one; 5-Androsten-3b-ol-17-one bromoacetate	881	3b-Hydroxy-5-androsten-19-al-17-one 3-sulfate, pyridinium salt
2482	3b-Bromoacetoxy-5-cholestene; 5-Cholesten-3b-ol bromoacetate; Cholesteryl bromoacetate	880	3b-Hydroxy-5-androstene-11,17-dione; 11-Ketodehydroepiandrosterone
2487	3b-Butoxy-5-cholestene; Cholesteryl butyl ether	2248	3b-Hydroxy-5-cholesten-24-oic acid
769	3b-Chloro-17,17-ethylenedioxy-5-androstene	2566	3b-Hydroxy-5-cholesten-22-one; 22-Ketocholesterol
499	3b-Chloro-1b,15a-diacetoxy-5a-androstan-17-one	1145	3b-Hydroxy-5-etienic acid hydrazide
2622	3b-Chloro-3a-methyl-5-cholestene	1141	3b-Hydroxy-5-etienic acid methyl ester
2327	3b-Chloro-3a-nitro-5a-cholestane	1144	3b-Hydroxy-5-etienic acid methyl ester 3-hemisuccinate
768	3b-Chloro-5-androsten-17-one	1661	3b-Hydroxy-5-pregnen-20-one 3-sulfate
		1667	3b-Hydroxy-5-pregnen-20-one glucuronide, sodium salt
		1748	3b-Hydroxy-5-pregnene-11,20-dione
		779	3b-Hydroxy-5a-androst-9(11)-en-17-one
		700	3b-Hydroxy-5a-androst-9(11)-ene; 5a-Androst-9(11)-en-3b-ol
		378	3b-Hydroxy-5a-androstan-11-hydrazone; 5a-Androstan-3b-ol-11-one hydrazone
		365	3b-Hydroxy-5a-androstan-17-one 2'-furoate
		336	3b-Hydroxy-5a-androstan-17-one 3-sulfate, potassium salt
		366	3b-Hydroxy-5a-androstan-17-one 5'-phenyl-2'-furoate
		373	3b-Hydroxy-5a-androstan-17-one b-D-glucoside

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
376	3b-Hydroxy-5a-androstan-17-one glucuronide; Epiandrosterone glucuronide	2670	3x-Bromo-A-homo-5a-cholestan-4-one
377	3b-Hydroxy-5a-androstan-17-one tetraacetyl-b-D-glucoside	2206	3x-Methoxy-12-oxo-5b-cholan-24-oic acid
367	3b-Hydroxy-5a-androstan-17-one thiophene-2-carboxylate	917	4,11-Androstadiene-3,17-dione
316	3b-Hydroxy-5a-androstan-7-one; 5a-Androstan-3b-ol-7-one	903	4,16-Androstadien-3-one
624	3b-Hydroxy-5a-androstane-11,16-dione	2890	4,22-Lumistadien-3-one
622	3b-Hydroxy-5a-androstane-6,17-dione	3004	4,4',4''-Trisnor-3-deoxy-b-amyrin
2143	3b-Hydroxy-5a-cholan-24-oic acid	3005	4,4',4''-Trisnor-a-amyrone
2426	3b-Hydroxy-5a-cholestan-6-one	3006	4,4',4''-Trisnor-b-amyrone
2422	3b-Hydroxy-5a-cholestan-6-one; 6-Ketocholestanol	2627	4,4,14a-Trimethyl-5a-cholesta-8(9),24-dien-3b-ol; Lanosterol
2430	3b-Hydroxy-5a-cholestan-7-one	963	4,4-Dimethyl-11a,16b-dihydroxy-5a-androstane-3,6-dione
2803	3b-Hydroxy-5a-ergosta-8,22-dien-7-one (crude)	2662	4,4-Dimethyl-2,3-seco-5-cholestene-2,3-dioic acid
1084	3b-Hydroxy-5a-etianic acid methyl ester	2658	4,4-Dimethyl-2,3-seco-5-cholestene-2,3-diol
1660	3b-Hydroxy-5a-pregn-16-en-20-one	989	4,4-Dimethyl-5-androsten-3-one
1749	3b-Hydroxy-5a-pregn-16-ene-11,20-dione	935	4,4-Dimethyl-5a-androstan-3-one
1750	3b-Hydroxy-5a-pregn-16-ene-12,20-dione	965	4,4-Dimethyl-5a-androstane-3,11,17-trione
1303	3b-Hydroxy-5a-pregnan-20-one; 5a-Pregnan-3b-ol-20-one	957	4,4-Dimethyl-5a-androstane-3,17-dione
1506	3b-Hydroxy-5a-pregnane-11,20-dione	968	4,4-Dimethyl-5a-androstane-3,6,11,16-tetraone
1507	3b-Hydroxy-5a-pregnane-12,20-dione	956	4,4-Dimethyl-5a-androstane-3,6-dione
1086	3b-Hydroxy-5b,17a-etianic acid methyl ester	964	4,4-Dimethyl-5a-androstane-3,7,17-trione
337	3b-Hydroxy-5b-androstan-17-one 3-sulfate, sodium salt	955	4,4-Dimethyl-7b,11a,17b-trihydroxy-5a-androstan-3-one
338	3b-Hydroxy-5b-androstan-17-one 3-sulfate, triethylammonium salt	2659	4,5-Secocholest-3-en-5-one semicarbazone (crude)
2253	3b-Hydroxy-5b-chol-11-en-24-oic acid methyl ester	2904	4,6,22-Lumistatrien-3-one
2146	3b-Hydroxy-5b-cholan-24-oic acid methyl ester	902	4,6-Androstadien-3-one
1143	3b-Hydroxy-5b-eti-11-enic acid methyl ester	916	4,6-Androstadiene-3,17-dione
1081	3b-Hydroxy-5b-etianic acid	2667	4-Acetoxy-3,4-epoxy-A-homo-5a-cholestane; 3,4-Epoxy-A-homo-5a-cholestan-4-ol acetate
1085	3b-Hydroxy-5b-etianic acid methyl ester	2672	4-Acetoxy-A-homo-5a-cholest-3-ene; A-Homo-5a-cholest-3-en-4-ol acetate
1990	3b-Hydroxy-6-methyl-5-pregnen-20-one; 6-Methylpregnenolone	705	4-Androsten-17b-ol
1155	3b-Hydroxy-7a-methoxy-5-etienic acid	748	4-Androsten-3-one
2263	3b-Hydroxy-7b-methoxy-5-cholen-24-oic acid	763	4-Androsten-3-one semicarbazone
2250	3b-Hydroxychol-5-en-24-oic acid 3-sulfate, sodium salt	891	4-Androstene-3,11,17-trione; Adrenosterone
2251	3b-Hydroxychol-5-en-24-oic acid methyl ester	859	4-Androstene-3,11-dione
1924	3b-Hydroxypregna-5,16-dien-20-one	860	4-Androstene-3,16-dione
298	3b-Iodo-5a-androstan-17-one (crude)	863	4-Androstene-3,17-dione
788	3b-Iodoacetoxo-5-androsten-17-one; 5-Androsten-3b-ol-17-one iodoacetate	858	4-Androstene-3,6-dione
708	3b-Methoxy-5-androstene	718	4-Androstene-3b,17b-diol
2570	3b-Methoxy-5-cholesten-7-one	2938	4-Bromo-19-norergosta-1,3,5,7,9-pentaen-3-ol
224	3b-Methoxy-5a-androstane-5,6b-diol	2552	4-Cholesten-3-one
2368	3b-Methoxy-5a-cholestan-5-ol; 5a-Cholestane-3b,5-diol 3-methyl ether	2557	4-Cholesten-6-one
2335	3b-Methoxy-5a-cholestane	2585	4-Cholestene-3,6-dione
2774	3b-Methoxy-5a-ergosta-7,22-dien-11a-ol; 5a-Ergosta-7,22-diene-3b,11a-diol 3-methyl ether	2731	4-Ergosten-3b-ol
1195	3b-Methoxy-5a-pregnane	12	4-Estren-17-one
2913	3b-Methoxy-9a-methyl-5a-ergost-7-en-11a-ol	11	4-Estren-3-one
2919	3b-Methoxy-9a-methyl-5a-ergosta-7,22-dien-11-one	23	4-Estrene-3,17-dione (crude)
2917	3b-Methoxy-9a-methyl-5a-ergosta-7,22-dien-11a-ol	44	4-Hydroxyestradiol; 1,3,5(10)-Estratriene-3,4,17b-triol
107	3b-N,N-Diacetylamino-5a-androstane	807	4-Oxa-5-androsten-17b-ol-3-one acetate; 17-Acetoxy-4-oxa-5-androsten-3-one
358	3b-Propionoxy-5a-androstan-17-one; 3b-Hydroxy-5a-androstan-17-one propionate	764	4-Oxa-5-androsten-3-one
2484	3b-Propoxy-5-cholestene; Cholesteryl propyl ether	2561	4-Oxa-5-cholesten-3-one
242	3b-Propoxy-5a-androstane-7b,12b,15a-triol; 5a-Androstane-3b,7b,12b,15a-tetrol 3-propyl ether	1065	4-Oxa-A-homo-5a-androstane-3,11,17-trione
792	3b-Succinoyloxy-5-androsten-17-one; 5-Androsten-3b-ol-17-one hemisuccinate	2669	4-Oxa-A-homo-5a-cholestan-3-one
2447	3x,7x-Dibenzoyloxy-5a-cholestane-6-one; 3x,7x-Dihydroxy-5a-cholestan-6-one dibenzoate	1869	4-Pregnen-3,18,20-trione; 18-Oxoprogesterone
2671	3x-Acetoxy-A-homo-5a-cholestan-4-one; 3x-Hydroxy-A-homo-5a-cholestan-4-one acetate	1649	4-Pregnen-3-one
		1783	4-Pregnene-11b,17a-diol-3,20-dione; 11b,17a-Dihydroxy-4-pregnene-3,20-dione
		1819	4-Pregnene-17a,20a,21-triol-3,11-dione; 17a,20a,21-Trihydroxy-4-pregnene-3,11-dione
		1820	4-Pregnene-17a,20b,21-triol-3,11-dione; 17a,20b,21-Trihydroxy-4-pregnene-3,11-dione

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2013	4-Pregnene-20b-carboxaldehyde-3-one; 3-Oxo-4-pregnene-20b-carboxaldehyde	2827	5-Acetoxy-3b-(p-toluenesulfonyloxy)-5a-ergosta-7,9(11),22-triene; 5a-Ergosta-7,9(11),22-triene-3b,5a-diol 5-acetate 3-tosylate
1867	4-Pregnene-3,11,20-trione; 11-Ketoprogesterone	2373	5-Acetoxy-5a-cholestan-3b-ol; 5a-Cholestane-3b,5-diol 5-acetate
1905	4-Pregnene-3,12,15,20-tetrone	2337	5-Acetoxy-5a-cholestane; 5a-Cholestan-5-ol acetate
1726	4-Pregnene-3,20-dione-[19-d3]	2779	5-Acetoxy-5a-ergosta-7,22-dien-3b-ol; 5a-Ergosta-7,22-diene-3b,5-diol 5-acetate
1725	4-Pregnene-3,20-dione; Progesterone	2678	5-Acetoxy-des-A-cholest-5(10)-ene; Des-A-cholest-5(10)-en-5-ol acetate
1904	4-Pregnene-3,6,11,20-tetrone	2677	5-Acetoxy-des-A-cholest-5-ene; Des-A-cholest-5-en-5-ol acetate
1865	4-Pregnene-3,6,20-trione; 6-Ketoprogesterone	753	5-Androsten-17-one
1722	4-Pregnene-3b,11b,17a,21-tetrol-20-one; 3b,11b,17a,21-Tetrahydroxy-4-pregnen-20-one	715	5-Androsten-3b-ol 3,5-dinitrobenzoate
1624	4-Pregnene-3b,20b-diol	750	5-Androsten-7-one
429	4a,17b-Dihydroxy-5a-androstan-7-one	696	5-Androstene
501	4a,5-Epoxy-17b-hydroxy-5a-androstan-3-one	892	5-Androstene-3,11,17-trione
397	4a,5-Epoxy-5a-androstan-3-one	719	5-Androstene-3a,17b-diol
638	4a,5-Epoxy-5a-androstane-3,17-dione	743	5-Androstene-3b,16a,17b,18-tetrol
2336	4a-Acetoxy-5a-cholestane; 5a-Cholestan-4a-ol acetate	736	5-Androstene-3b,16a,17b-triol
2371	4a-Acetoxy-5b-cholestan-3b-ol; 5a-Cholestane-3b,4a-diol 4-acetate	744	5-Androstene-3b,16b,17b,18-tetrol
967	4a-Hydroxymethyl-4b-methyl-5a-androstane-3,6,16-trione	720	5-Androstene-3b,17a-diol
2381	4b,5-Diacetoxy-5a-cholestane; 5a-Cholestane-4b,5-diol diacetate	727	5-Androstene-3b,17a-diol disulfate, dipotassium salt
502	4b,5-Epoxy-17b-hydroxy-5b-androstan-3-one	721	5-Androstene-3b,17b-diol
398	4b,5-Epoxy-5b-androstan-3-one	726	5-Androstene-3b,17b-diol 17-sulfate, potassium salt
639	4b,5-Epoxy-5b-androstane-3,17-dione	728	5-Androstene-3b,17b-diol 3-acetate; 3b-Acetoxy-5-androsten-17b-ol
2406	4b,5-Epoxy-5b-cholestane-3b,6b-diol	724	5-Androstene-3b,17b-diol 3-sulfate, ammonium salt
123	4b-Acetoxy-5a-androstane	725	5-Androstene-3b,17b-diol 3-sulfate, sodium salt
2370	4b-Acetoxy-5a-cholestan-3a-ol; 5a-Cholestane-3a,4b-diol 4-acetate	729	5-Androstene-3b,17b-diol diacetate; 3b,17b-Diacetoxy-5-androstene
2338	4b-Benzoyloxy-5a-cholestane; 5a-Cholestan-4b-ol benzoate	730	5-Androstene-3b,17b-diol dipropionate; 3b,17b-Dipropionyloxy-5-androstene
2325	4b-Nitro-5b-cholestane; 4b-Nitrocoprostan	866	5-Androstene-7,12-dione
2934	4x-Acetoxy-19-nor-5,7,9-ergostatrien-3-one	2682	5-Benzoyloxy-des-A-cholesta-5,7,9,11-tetraene; Des-A-cholesta-5,7,9,11-tetraen-5-ol benzoate
1124	5,14,19-Trihydroxy-3-oxo-5b,14b-etianic acid methyl ester	2681	5-Benzoyloxy-des-A-cholesta-5,7,9-triene; Des-A-cholesta-5,7,9-trien-5-ol benzoate
1550	5,14,19-Trihydroxy-5b,14b,17a-pregnane-3,20-dione	2344	5-Bromo-5a-cholestan-6b-ol
1482	5,14-Dihydroxy-3b,19,21-triacetoxy-5b,14b,17a-pregnan-20-one	506	5-Bromo-6b,19-epoxy-3b-hydroxy-5a-androstan-17-one 3-sulfate, pyridinium salt
1481	5,14-Dihydroxy-3b,19,21-triacetoxy-5b,14b-pregnan-20-one	1325	5-Chloro-3b-hydroxy-5a-pregnan-20-one
1483	5,17a-Dihydroxy-3b,6b,21-triacetoxy-5a-pregnan-20-one	2130	5-Cholene-3b,24-diol
1069	5,17b-Diacetoxy-des-A-10b(H)-androst-5-ene; Des-A-10b(H)-androst-5-ene-5,17b-diol diacetate	2525	5-Cholest-3b-ol digitonide; Cholesteryl digitonide
439	5,17b-Dihydroxy-5a-androstan-11-one	2553	5-Cholesten-3-one (approx.97%)
1118	5,19-Dihydroxy-3-oxo-5b-etianic acid ethyl ester	2528	5-Cholesten-3b-ol 3',5'-dinitrobenzoate; Cholesteryl 3',5'-dinitrobenzoate
1514	5,19-Dihydroxy-5b-pregnane-3,20-dione	2481	5-Cholesten-3b-ol dichloroacetate; Cholesteryl dichloroacetate
2091	5,20-Dioxo-3,5-seco-A-norpregnan-3-oic acid	2483	5-Cholesten-3b-ol diiodoacetate; Cholesteryl diiodoacetate
2968	5,24(28)-Stigmastadien-3b-ol; Fucosterol	2491	5-Cholesten-3b-ol hemisuccinate; Cholesteryl hydrogen succinate
1380	5,6-Epoxy-3b-acetoxy-5a-pregnan-20-one	2498	5-Cholesten-3b-ol heptanoate; Cholesteryl enanthate
195	5,6-Epoxy-3b-methoxy-5a-androstan	2512	5-Cholesten-3b-ol octadecanoate; Cholesteryl stearate
2812	5,6a-(Isopropylidenedioxy)-5a-ergosta-7,22-dien-3-one	2485	5-Cholesten-3b-ol propionate; Cholesteryl propionate
1378	5,6a-Epoxy-3b-hydroxy-5a-pregnan-20-one	2565	5-Cholesten-3b-ol-7-one; 3b-Hydroxy-5-cholesten-7-one
194	5,6a-Epoxy-5a-androstan-3b-ol	2558	5-Cholesten-7-one
2884	5,6b-(Isopropylidenedioxy)-5b-lumista-7,22-dien-3b-ol; 5b-Lumista-7,22-diene-3b,5,6b-triol 5,6-acetonide	2463	5-Cholestene
1374	5,6b-Dichloro-3b,17a-dihydroxy-5a-pregnan-20-one	2534	5-Cholestene-3b,19-diol; 19-Hydroxycholesterol
733	5,6b-Dichloro-5a-androst-15-ene-3b,17b-diol diacetate	2535	5-Cholestene-3b,22R-diol; 22R-Hydroxycholesterol
646	5,6b-Dihydroxy-5a-androstane-3,17-dione 6-acetate; 6b-Acetoxy-5-hydroxy-5a-androstane-3,17-dione	2536	5-Cholestene-3b,22S-diol; 22S-Hydroxycholesterol
504	5,6b-Epoxy-3b-hydroxy-5b-androstan-17-one	2531	5-Cholestene-3b,4a-diol; 4a-Hydroxycholesterol
2872	5,6b-Epoxy-5b-lumista-7,22-dien-3b-ol; Lumisterol b-epoxide	2532	5-Cholestene-3b,4b-diol; 4b-Hydroxycholesterol
1699	5,8a-Epidioxy-3b-hydroxy-5a-pregn-6-en-20-one		
2828	5,8a-Peroxydo-5a-ergosta-6,9(11),22-trien-3b-ol (impure)		
1070	5-Acetoxy-17b-benzoyloxy-des-A-androst-5(10)-ene; Des-A-androst-5(10)-ene-5,17b-diol 5-acetate 17-benzoate		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
10	5-Estrene-3b,17b-diol 3-acetate 17-benzoate	503	5a,6a-Epoxy-3b-hydroxyandrostan-17-one
9	5-Estrene-3b,17b-diol; 19-Nor-5-androstene-3b,17b-diol	2691	5a,8a-Ergostane-3b,11a-diol
2814	5-Hydroxy-5a,9b-ergosta-7,22-diene-3,6-dione	2754	5a,8a-Peroxido-5a-ergost-22-en-3b-ol
691	5-Hydroxy-5a-androstane-3,6,17-trione	2891	5a,9a-Lumista-7,22-dien-3-one
2425	5-Hydroxy-5a-cholestan-4-one	2842	5a,9a-Lumistan-3-one
2427	5-Hydroxy-5a-cholestan-6-one	2834	5a,9a-Lumistan-3a-ol
2458	5-Hydroxy-5a-cholestane-3,6-dione	2836	5a,9a-Lumistan-3b-ol acetate; 3b-Acetoxy-5a,9a-lumistane
1598	5-Hydroxy-5a-pregnane-3,6,20-trione	2787	5a,9b-Ergosta-1,22-diene-3b,5,6a-triol
314	5-Hydroxy-5b-androstan-6-one; 5b-Androstan-5-ol-6-one	2715	5a,9b-Ergostan-3-one
2424	5-Hydroxy-5b-cholestan-3-one	2685	5a,9b-Ergostan-3b-ol
2428	5-Hydroxy-5b-cholestan-6-one	2728	5a,9b-Ergostane-3,11-dione
2641	5-Methyl-19-nor-5b-cholest-9-ene-3,6-dione	2692	5a,9b-Ergostane-3b,11a-diol
2940	5-Oxo-A-nor-3,5-seco-22-lumisten-3-oic acid	747	5a-Androst-1-en-3-one
2647	5-Oxo-A-nor-3,5-secocholestan-3-oic acid	698	5a-Androst-1-en-3a-ol
1288	5-Pregnane-3a,11b,17,20,21-pentol 3-gluconuride; 3a,11b,17,20,21-Pentahydroxy-5b-pregnane 3-gluconuride	699	5a-Androst-1-en-3b-ol
1659	5-Pregnen-3b-ol-20-one; 3b-Hydroxy-5-pregnen-20-one	695	5a-Androst-1-ene
1728	5-Pregnene-3,20-dione	862	5a-Androst-1-ene-3,17-dione
1635	5-Pregnene-3b,16a,20a-triol	890	5a-Androst-1-ene-3,6,17-trione
1647	5-Pregnene-3b,17a,20a,21-tetrol	723	5a-Androst-1-ene-6a,17b-diol
1636	5-Pregnene-3b,17a,20a-triol	755	5a-Androst-14-en-17-one
1703	5-Pregnene-3b,17a,20a-triol-11-one; 3b,17a,20a-Trihydroxy-4-pregnen-11-one	756	5a-Androst-15-en-17-one
1637	5-Pregnene-3b,17a,20b-triol	707	5a-Androst-15-en-17b-ol
1704	5-Pregnene-3b,17a,20b-triol-11-one; 3b,17a,20b-Trihydroxy-4-pregnen-11-one	749	5a-Androst-16-en-3-one
1623	5-Pregnene-3b,17b-diol; 17-Ethyl-5-androstene-3b,17b-diol	701	5a-Androst-16-en-3a-ol
1639	5-Pregnene-3b,20a,21-triol	702	5a-Androst-16-en-3b-ol
1625	5-Pregnene-3b,20a-diol	746	5a-Androst-2-en-1-one
1631	5-Pregnene-3b,20a-diol diacetate; 3b,20-Diacetoxy-5-pregnene	752	5a-Androst-2-en-16-one
1640	5-Pregnene-3b,20b,21-triol	754	5a-Androst-2-en-17-one
1626	5-Pregnene-3b,20b-diol	697	5a-Androst-2-en-1b-ol
722	5?-Androst-6-ene-3b,17b-diol	751	5a-Androst-2-en-7-one
527	5a,11a,17b-Trihydroxy-5a-androstan-6-one 11,17-diacetate	856	5a-Androst-2-ene-1,17-dione
757	5a,13a-Androst-15-en-17-one	717	5a-Androst-2-ene-1a,17b-diol
704	5a,13a-Androst-16-en-15a-ol	706	5a-Androst-9(11)-en-17b-ol
257	5a,13a-Androstan-15-one	745	5a-Androst-9(11)-en-3-one
115	5a,13a-Androstan-15a-ol	857	5a-Androst-9(11)-ene-2,17-dione
893	5a,14b-Androst-1-ene-3,12,15-trione	861	5a-Androst-9(11)-ene-3,16-dione
758	5a,14b-Androst-15-en-17-one (crude)	864	5a-Androst-9(11)-ene-3,17-dione
259	5a,14b-Androstan-15-one	245	5a-Androstan-1-one
117	5a,14b-Androstan-15b-ol	254	5a-Androstan-11-one
262	5a,14b-Androstan-17-one	256	5a-Androstan-12-one
570	5a,14b-Androstane-12,15-dione	113	5a-Androstan-12a-ol
663	5a,14b-Androstane-2,12,15-trione	258	5a-Androstan-15-one
682	5a,14b-Androstane-6,12,15-trione (deuterated)	114	5a-Androstan-15a-ol
564	5a,14b-Androstane-7,12-dione	116	5a-Androstan-15b-ol
167	5a,14b-Androstane-7b,12b-diol	260	5a-Androstan-16-one
2716	5a,14b-Ergostan-3-one	118	5a-Androstan-16b-ol
2686	5a,14b-Ergostan-3b-ol	261	5a-Androstan-17-one
2835	5a,14b-Lumistan-3a-ol	266	5a-Androstan-17-one tosylhydrazone
2838	5a,14b-Lumistan-3a-ol 3,5-dinitrobenzoate	119	5a-Androstan-17a-ol
1079	5a,17a-Etanic acid methyl ester	120	5a-Androstan-17b-ol
1734	5a,17a-Pregn-2-ene-11,20-dione (crude)	141	5a-Androstan-17b-ol tosylate; 17b-(p-Toluenesulfonyloxy)-5a-androstane
1296	5a,17a-Pregnan-20-one	312	5a-Androstan-17b-ol-3-one; 17b-Hydroxy-5a-androstan-3-one
1204	5a,17a-Pregnan-20a-ol phenylglyoxylate	247	5a-Androstan-2-one
1194	5a,17a-Pregnan-20b-ol	248	5a-Androstan-3-one
1249	5a,17a-Pregnane-3b,17b,20x-triol	135	5a-Androstan-3a-ol 2'-hydroxyethyl ether; 3a-(2-Hydroxyethoxy)-5a-androstane
1206	5a,17a-Pregnane-3b,17b-diol	136	5a-Androstan-3a-ol carboxymethyl ether
403	5a,17b-Dihydroxy-5a-androstan-3-one	108	5a-Androstan-3b-ol

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
138	5a-Androstan-3b-ol 2'-acetoxyethyl ether; 3b-(2-Acetoxyethoxy)-5a-androstane	667	5a-Androstane-3,7,11,-trione
140	5a-Androstan-3b-ol benzenesulfonate; 3b-Benzenesulfonyloxy-5a-androstane	693	5a-Androstane-3,7,11,16-tetraone
137	5a-Androstan-3b-ol methoxycarbonylmethyl ether	694	5a-Androstane-3,7,11,17-tetraone
139	5a-Androstan-3b-ol methyl succinate	668	5a-Androstane-3,7,16-dione
319	5a-Androstan-3b-ol-16-one; 3b-Hydroxy-5a-androstan-16-one (90-95%)	669	5a-Androstane-3,7,17-trione
1082	5a-Androstan-3b-ol-17b-carboxylic acid; 3b-Hydroxy-5a-androstan-17b-carboxylic acid; 3b-Hydroxy-5a-etianic acid	210	5a-Androstane-3a,11a,17b-triol
250	5a-Androstan-4-one (crude)	212	5a-Androstane-3a,11b,17b-triol (approx.95%)
109	5a-Androstan-4a-ol	216	5a-Androstane-3a,12b,15a-triol (crude)
252	5a-Androstan-6-one	157	5a-Androstane-3a,17b-diol
111	5a-Androstan-6a-ol	193	5a-Androstane-3a,17b-diol 17-b-D-glucuronide, potassium salt
112	5a-Androstan-6b-ol (crude)	208	5a-Androstane-3b,11a,16b-triol
550	5a-Androstane-1,17-dione	214	5a-Androstane-3b,11a,17b-triol
655	5a-Androstane-1,3,17-trione	209	5a-Androstane-3b,11b,16b-triol
656	5a-Androstane-1,6,16-trione	215	5a-Androstane-3b,11b,17b-triol
549	5a-Androstane-1,6-dione (crude)	456	5a-Androstane-3b,11b-diol-17-one; 3b,11b-Dihydroxy-5a-androstan-17-one
692	5a-Androstane-1,7,11,17-tetraone	218	5a-Androstane-3b,12a,16a-triol
658	5a-Androstane-1,7,12-trione	217	5a-Androstane-3b,12b,15a-triol
567	5a-Androstane-11,16-dione	219	5a-Androstane-3b,16a,17b-triol
568	5a-Androstane-11,17-dione	160	5a-Androstane-3b,17b-diol
170	5a-Androstane-11a,17b-diol	176	5a-Androstane-3b,17b-diol ditrityl ether
169	5a-Androstane-11b,16b-diol	198	5a-Androstane-3b,4a,12a-triol
569	5a-Androstane-12,15-dione	152	5a-Androstane-3b,4a-diol
571	5a-Androstane-12,17-dione	199	5a-Androstane-3b,5,6b-triol
188	5a-Androstane-15b,17b-diol 17-tosylate; 17b-(p-Toluenesulfonyloxy)-5a-androstan-15b-ol	200	5a-Androstane-3b,6a,16a-triol
171	5a-Androstane-16a,17a-diol	201	5a-Androstane-3b,6a,16b-triol
189	5a-Androstane-16a,17a-diol 16-tosylate; 16a-(p-Toluenesulfonyloxy)-5a-androstan-17a-ol	202	5a-Androstane-3b,6a,17a-triol
172	5a-Androstane-16a,17b-diol	154	5a-Androstane-3b,6a-diol
190	5a-Androstane-16a,17b-diol 16-tosylate	204	5a-Androstane-3b,7a,16b-triol
182	5a-Androstane-16a,17b-diol 17-mesylyate; 17b-Mesyloxy-5a-androstan-16a-ol	203	5a-Androstane-3b,7b,11a-triol
173	5a-Androstane-16b,17b-diol	205	5a-Androstane-3b,7b,17b-triol
191	5a-Androstane-16b,17b-diol 16-tosylate	155	5a-Androstane-3b,7b-diol
149	5a-Androstane-1a,6a-diol	206	5a-Androstane-3b,9a,16a-triol
196	5a-Androstane-1b,6a,17b-triol	207	5a-Androstane-3b,9a,17b-triol
151	5a-Androstane-1b,6a-diol	156	5a-Androstane-3b,9a-diol
662	5a-Androstane-2,11,17-trione	677	5a-Androstane-4,11,16-trione
552	5a-Androstane-2,17-dione	678	5a-Androstane-4,11,17-trione
659	5a-Androstane-2,6,17-trione	679	5a-Androstane-6,11,15-trione
660	5a-Androstane-2,7,16-trione	680	5a-Androstane-6,11,16-trione
661	5a-Androstane-2,7,17-trione	681	5a-Androstane-6,11,17-trione
551	5a-Androstane-2,7-dione	560	5a-Androstane-6,11-dione
197	5a-Androstane-2a,6a,16a-triol	683	5a-Androstane-6,12,16-trione
671	5a-Androstane-3,11,15-trione	561	5a-Androstane-6,12-dione
672	5a-Androstane-3,11,16-trione	562	5a-Androstane-6,17-dione
554	5a-Androstane-3,11-dione	220	5a-Androstane-6a,11a,17b-triol
674	5a-Androstane-3,12,15-trione	162	5a-Androstane-6a,11a-diol
675	5a-Androstane-3,12,16-trione	221	5a-Androstane-6a,12b,15a-triol
556	5a-Androstane-3,12-dione	163	5a-Androstane-6a,12b-diol
557	5a-Androstane-3,15-dione (crude)	164	5a-Androstane-6a,16a-diol
558	5a-Androstane-3,16-dione	165	5a-Androstane-6a,17b-diol
559	5a-Androstane-3,17-dione	684	5a-Androstane-7,11,17-trione
664	5a-Androstane-3,6,11-trione	563	5a-Androstane-7,12-dione
665	5a-Androstane-3,6,16-trione	565	5a-Androstane-7,16-dione
666	5a-Androstane-3,6,17-trione	566	5a-Androstane-7,17-dione
553	5a-Androstane-3,6-dione	168	5a-Androstane-7a,16b-diol
		222	5a-Androstane-7b,12b,15a-triol
		225	5a-Androstane-7b,12b,15a-triol 7,15-acetonide
		166	5a-Androstane-7b,12b-diol
		2908	5a-Campest-8(14)-en-3b-ol

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2361	5a-Chloestane-4b,5-diol	2786	5a-Ergosta-7,22-diene-3a,5-diol carbonate
2131	5a-Chol-11-ene-3a,24-diol	2771	5a-Ergosta-7,22-diene-3a,5a-diol
2137	5a-Cholan-24-oic acid methyl ester	2773	5a-Ergosta-7,22-diene-3b,11b-diol
2202	5a-Cholanic acid-3a-ol-6-one; 3a-Hydroxy-6-oxo-5a-cholan-24-oic acid	2790	5a-Ergosta-7,22-diene-3b,5,11b-triol
2462	5a-Cholest-1-ene	2788	5a-Ergosta-7,22-diene-3b,5,6a-triol
2554	5a-Cholest-6-en-3-one	2772	5a-Ergosta-7,22-diene-3b,5a-diol
2555	5a-Cholest-7-en-3-one	2684	5a-Ergostan-3b-ol
2471	5a-Cholest-7-en-3b-ol	2727	5a-Ergostane-3,11-dione
2560	5a-Cholest-8(14)-en-7-one	2690	5a-Ergostane-3b,11a-diol
2556	5a-Cholest-8(9)-3-one (crude)	2693	5a-Ergostane-3b,11b-diol
2559	5a-Cholest-8-en-7-one	2708	5a-Ergostane-3b,7a,8a-triol
2408	5a-Cholestan-1-one	1	5a-Estran-11a-ol
2409	5a-Cholestan-2-one	3	5a-Estran-17-one
2410	5a-Cholestan-3-one	8	5a-Estrane-3,11,16-trione
2334	5a-Cholestan-3b-ol sulfate, sodium salt; Dihydrocholesterol 3-sulfate, sodium salt	7	5a-Estrane-3,16-dione
2342	5a-Cholestan-3b-ol trityl ether	2	5a-Estrane-3b,11a,16b-triol
2330	5a-Cholestan-4b-ol	1075	5a-Etanic acid
2331	5a-Cholestan-5-ol	1080	5a-Etanic acid hydrazide
2332	5a-Cholestan-6a-ol	1076	5a-Etanic acid methyl ester
2339	5a-Cholestan-6a-ol 3',5'-dinitrobenzoate	304	5a-Hydroxyandrostan-3-one; Androstan-5a-ol-3-one
2333	5a-Cholestan-6b-ol (crude)	2951	5a-Poriferastan-3b-ol; 24a-Stigmastanol
2340	5a-Cholestan-6b-ol 3',5'-dinitrobenzoate	1619	5a-Preg-17(20)-ene
2322	5a-Cholestane	1866	5a-Pregn-1-ene-3,11,20-trione
2347	5a-Cholestane-1a,2a-diol	1653	5a-Pregn-16-en-20-one
2348	5a-Cholestane-1a,2b-diol	1868	5a-Pregn-16-ene-3,11,20-trione
2351	5a-Cholestane-1a,3b-diol	1731	5a-Pregn-16-ene-3,20-dione
2349	5a-Cholestane-1b,2a-diol	1651	5a-Pregn-2-en-20-one
2350	5a-Cholestane-1b,2b-diol	1733	5a-Pregn-2-ene-11,20-dione
2451	5a-Cholestane-2,3-dione	1652	5a-Pregn-3-en-20-one
2352	5a-Cholestane-2b,3a-diol	1650	5a-Pregn-9(11)-en-3-one
2452	5a-Cholestane-3,6-dione	1620	5a-Pregn-9(11)-en-3b-ol
2354	5a-Cholestane-3a,4b-diol	1936	5a-Pregna-2,16-diene-11,20-dione
2355	5a-Cholestane-3b,4a-diol	1191	5a-Pregnan-16a-ol
2398	5a-Cholestane-3b,5,6a-triol	1294	5a-Pregnan-20-one
2405	5a-Cholestane-3b,5,6b-triol 3-(3',5'-dinitrobenzoate)	1981	5a-Pregnan-20-one 16a,17a-pyrazoline
2358	5a-Cholestane-3b,6a-diol	1192	5a-Pregnan-20a-ol
2359	5a-Cholestane-4a,5-diol	1298	5a-Pregnan-20a-ol-3-one; 20a-Hydroxy-5a-pregnan-3-one
2364	5a-Cholestane-6a,7a-diol	1193	5a-Pregnan-20b-ol
2365	5a-Cholestane-6a,7b-diol	1501	5a-Pregnan-21-ol-3,20-dione; 21-Hydroxy-5a-pregnane-3,20-dione
2366	5a-Cholestane-6b,7a-diol	1186	5a-Pregnan-3a-ol
2367	5a-Cholestane-7a,8a-diol	1187	5a-Pregnan-3b-ol
2456	5a-Cholestane[2,3-c]furan	1315	5a-Pregnan-3b-ol-20-one acetate; 3b-Acetoxy-5a-pregnan-20-one
2457	5a-Cholestane[2,3-c]furoxan	1183	5a-Pregnane
1122	5a-Dihydroaldosterone g-etiolactone; 11b,18-Epoxy-18-hydroxy-3-oxo-5a-androstane-17b-carboxylic acid 17,18-lactone	1615	5a-Pregnane-1,7,12,20-tetrone
2734	5a-Ergost-14-en-3b-ol	1551	5a-Pregnane-11b,17a,21-triol-3,20-dione; 11b,17a,21-Trihydroxy-5a-pregnane-3,20-dione
2762	5a-Ergost-22-en-3-one	1515	5a-Pregnane-11b,21-diol-3,20-dione; 11b,21-Dihydroxy-5a-pregnane-3,20-dione
2756	5a-Ergost-22-ene-3a,5a,8a-triol	1214	5a-Pregnane-16b,20a-diol
2750	5a-Ergost-22-ene-3b,11a-diol	1215	5a-Pregnane-16b,20b-diol
2733	5a-Ergost-8(14)-en-3b-ol	1216	5a-Pregnane-17a,20a-diol
2761	5a-Ergost-8-en-3-one	1252	5a-Pregnane-17a,20b,21-triol
2732	5a-Ergost-8-en-3b-ol	1609	5a-Pregnane-17a,21-diol-3,11,20-trione 21-acetate; 21-Acetoxy-17a-hydroxy-5a-pregnane-3,11,20-one (approx.95%)
2746	5a-Ergost-8-en-3b-ol 3',5'-dinitrobenzoate	1517	5a-Pregnane-17a,21-diol-3,20-dione; 17a,21-Dihydroxy-5a-pregnane-3,20-dione
2730	5a-Ergost-8-ene	1488	5a-Pregnane-2,20-dione
2789	5a-Ergosta-6,22-diene-3a,5,8a-triol	1591	5a-Pregnane-3,11,20-trione
2801	5a-Ergosta-7,22-dien-3-one		
2769	5a-Ergosta-7,22-diene		
2813	5a-Ergosta-7,22-diene-3,6-dione		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1618	5a-Pregnane-3,12,15,20-tetrone	676	5b-Androstane-3,12,17-trione
1489	5a-Pregnane-3,20-dione	670	5b-Androstane-3,7,17-trione
1588	5a-Pregnane-3,6,20-trione (approx.96%)	211	5b-Androstane-3a,11a,17b-triol
1616	5a-Pregnane-3,7,11,20-tetrone (crude)	213	5b-Androstane-3a,11b,17b-triol
1590	5a-Pregnane-3,7,20-trione	158	5b-Androstane-3a,17b-diol
1454	5a-Pregnane-3a,11b,17a,21-tetrol-20-one; 3a,11b,17a,21-Tetrahydroxy-5a-pregnan-20-one (90-95%)	153	5b-Androstane-3a,6b-diol
1245	5a-Pregnane-3a,17a,20b-triol	457	5b-Androstane-3b,11b-diol-17-one; 3b,11b-Dihydroxy-5b-androstan-17-one
1554	5a-Pregnane-3a,17a,21-triol-11,20-dione; 3a,17a,21-Trihydroxy-5a-pregnane-11,20-dione	159	5b-Androstane-3b,17a-diol
1207	5a-Pregnane-3a,20a-diol	161	5b-Androstane-3b,17b-diol
1209	5a-Pregnane-3a,20b-diol	2390	5b-Bromo-5a-cholestane-3b,7a-diol
1286	5a-Pregnane-3b,11b,17a,20b,21-pentol	2247	5b-Chol-11-en-24-oic acid methyl ester
1456	5a-Pregnane-3b,11b,17a,21-tetrol-20-one; 3b,11b,17a,21-Tetrahydroxy-5a-pregnan-20-one	2269	5b-Chola-3,11-dien-24-oic acid methyl ester
1275	5a-Pregnane-3b,11b,20b,21-tetrol	2138	5b-Cholan-24-oic acid methyl ester
1277	5a-Pregnane-3b,17a,20a,21-tetrol	2136	5b-Cholan-24-oic acid; 5b-Cholanic acid
1247	5a-Pregnane-3b,17a,20a-triol; Reichstein's compound O	2126	5b-Cholan-24-ol
1278	5a-Pregnane-3b,17a,20b,21-tetrol	2128	5b-Cholane-3a,12a,24-triol
1248	5a-Pregnane-3b,17a,20b-triol; Reichstein's compound J	2127	5b-Cholane-3a,24-diol
1556	5a-Pregnane-3b,17a,21-triol-11,20-dione; 3b,17a,21-Trihydroxy-5a-pregnane-11,20-dione	2129	5b-Cholane-3a,7a,12a,24-tetrol
1401	5a-Pregnane-3b,17a,21-triol-20-one; 3b,17a,21-Trihydroxy-5a-pregnan-20-one	2232	5b-Cholanic acid-3,6-dione; 3,6-Dioxo-5b-cholan-24-oic acid
1211	5a-Pregnane-3b,20a-diol	2194	5b-Cholanic acid-3-one; 3-Oxo-5b-cholan-24-oic acid
1213	5a-Pregnane-3b,20b-diol	2203	5b-Cholanic acid-3a-ol-6-one; 3a-Hydroxy-6-oxo-5b-cholan-24-oic acid
1525	5a-Pregnane-3b,21-diol-11,20-dione; 3b,21-Dihydroxy-5a-pregnane-11,20-dione	2328	5b-Cholestan-3a-ol; Coprostan-3-ol
1491	5a-Pregnane-7,20-dione	2412	5b-Cholestan-6-one; Coprostan-6-one
2955	5a-Stigmast-7-en-3b-ol	2353	5b-Cholestane-2b,3a-diol
2956	5a-Stigmast-8(14)-en-3a-ol	2356	5b-Cholestane-3a,5-diol
2957	5a-Stigmast-8(14)-en-3b-ol	2400	5b-Cholestane-3a,7a,12a-triol
2969	5a-Stigmasta-7,22(E)-dien-3b-ol	2357	5b-Cholestane-3b,5-diol
2950	5a-Stigmastan-3b-ol; Stigmastanol	2360	5b-Cholestane-4a,5-diol
1078	5b,14b-Etanic acid methyl ester	2362	5b-Cholestane-4b,5-diol
1291	5b,14b-Pregnane-3b,5,14,19,20b,21-hexol	2363	5b-Cholestane-5,6a-diol
1594	5b,17a-Pregnane-3,12,20-trione	2802	5b-Ergosta-7,22-dien-3-one
2610	5b,25(R)-Spirostan-3b-ol; Smilagenin	2770	5b-Ergosta-7,22-dien-3b-ol trimethylsilyl ether
2896	5b,6b-(Isopropylidenedioxy)-9a-lumista-7,22-dien-3-one; 5b,6b-Dihydroxy-9a-lumista-7,22-dien-3-one 5,6-acetonide	1077	5b-Etanic acid methyl ester
2871	5b,6b-Epoxyllumista-7,22-dien-3a-ol	2897	5b-Hydroxy-9a-lumista-7,22-diene-3,6-dione
2862	5b,9a-Lumista-7,22-dien-3b-ol 3,5-dinitrobenzoate	2852	5b-Lumist-7-en-3-one
2878	5b,9a-Lumista-7,22-diene-3b,5,6a-triol	2846	5b-Lumist-7-en-3a-ol
2683	5b,9b-Ergostane	2847	5b-Lumist-7-en-3b-ol
889	5b-Androst-1-ene-3,6,11-trione	2865	5b-Lumista-6,22-dien-3b-ol 3,5-dinitrobenzoate 5b,8b-maleic anhydride adduct
865	5b-Androst-11-ene-3,17-dione	2863	5b-Lumista-6,22-dien-3b-ol 5b,8b-maleic anhydride adduct
888	5b-Androst-2-ene-1,6,16-trione	2892	5b-Lumista-7,22-dien-3-one
855	5b-Androst-2-ene-1,6-dione	2877	5b-Lumista-7,22-dien-3a,5,6b-triol
703	5b-Androst-3-en-5-ol	2854	5b-Lumista-7,22-dien-3a-ol
778	5b-Androst-9(11)-en-3a-ol-17-one; 3a-Hydroxy-5b-androst-9(11)-en-17-one	2860	5b-Lumista-7,22-dien-3a-ol 3,5-dinitrobenzoate
246	5b-Androstan-1-one	2866	5b-Lumista-7,22-dien-3b,5-diol
255	5b-Androstan-11-one	2868	5b-Lumista-7,22-dien-3b,5-diol 3,5-carbonate
263	5b-Androstan-17-one	2855	5b-Lumista-7,22-dien-3b-ol
121	5b-Androstan-17b-ol	2861	5b-Lumista-7,22-dien-3b-ol 3,5-dinitrobenzoate
249	5b-Androstan-3-one	2888	5b-Lumista-7,22-diene-3b,5,6a-triol 3-(3,5-dinitrobenzoate)
251	5b-Androstan-4-one	2906	5b-Lumista-7,9,(11),22-triene-3,6-dione
110	5b-Androstan-5-ol	2903	5b-Lumista-8,14,22-trien-3b-ol 3,5-dinitrobenzoate
253	5b-Androstan-6-one	2856	5b-Lumista-8,14-dien-3b-ol
657	5b-Androstane-1,6,17-trione	2843	5b-Lumistan-3-one
673	5b-Androstane-3,11,17-trione	2837	5b-Lumistan-3b-ol acetate; 3b-Acetoxy-5b-lumistane
555	5b-Androstane-3,11-dione	2839	5b-Lumistane-2a,3b-diol
		2840	5b-Lumistane-2b,3b-diol
		2841	5b-Lumistane-3b,7x,8x-triol
		1654	5b-Pregn-16-en-20-one
		1732	5b-Pregn-16-ene-3,20-dione

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
1729	5b-Pregn-6-ene-3,20-dione	1524	5b-Pregnane-3a,21-diol-11,20-dione; 3a,21-Dihydroxy-5b-pregnane-11,20-dione
1627	5b-Pregn-6-ene-3a,20a-diol	1269	5b-Pregnane-3a,6a,17a,20a-tetrol
1638	5b-Pregn-9(11)-ene-3a,20b,21-triol	1283	5b-Pregnane-3a,6a,17a,20b,21-pentol
1292	5b-Pregnan-11-one	1270	5b-Pregnane-3a,6a,17a,20b-tetrol
1293	5b-Pregnan-12-one	1240	5b-Pregnane-3a,6a,20a-triol
1190	5b-Pregnan-12a-ol	1271	5b-Pregnane-3a,6a,20b,21-tetrol
1500	5b-Pregnan-17a-ol-3,20-dione; 17a-Hydroxy-5b-pregnane-3,20-dione	1241	5b-Pregnane-3a,6a,20b-triol
1295	5b-Pregnan-20-one	1332	5b-Pregnane-3a,6a-diol-20-one; 3a,6a-Dihydroxy-5b-pregnan-20-one
1299	5b-Pregnan-20a-ol-3-one; 20a-Hydroxy-5b-pregnan-3-one	1457	5b-Pregnane-3b,11b,17a,21-tetrol-20-one; 3b,11b,17a,21-Tetrahydroxy-5b-pregnan-20-one
1601	5b-Pregnan-21-ol-3,11,20-trione; 21-Hydroxy-5b-pregnane-3,11,20-trione	1397	5b-Pregnane-3b,11b,21-triol-20-one; 3b,11b,21-Trihydroxy-5b-pregnan-20-one
1188	5b-Pregnan-3a-ol	1243	5b-Pregnane-3b,16b,20a-triol
1505	5b-Pregnan-3a-ol-11,20-dione; 3a-Hydroxy-5b-pregnane-11,20-dione (approx.90%)	1402	5b-Pregnane-3b,17a,21-triol-20-one; 3b,17a,21-Trihydroxy-5b-pregnan-20-one (approx.95%)
1322	5b-Pregnan-3a-ol-20-one 3-glucuronide, sodium salt	1212	5b-Pregnane-3b,20a-diol
1302	5b-Pregnan-3a-ol-20-one; 3a-Hydroxy-5b-pregnan-20-one	1139	5x-Etiochol-14-enic acid; 5x-Eti-14-enic acid
1189	5b-Pregnan-3b-ol	2567	5x-Hydroxy-4-oxacholest-1-en-3-one
1304	5b-Pregnan-3b-ol-20-one; 3b-Hydroxy-5b-pregnan-20-one	1670	5x-Pregn-17(20)-en-21-oic acid ethyl ester
1184	5b-Pregnane	2028	6,16-Dimethyl-16-dehydropregnenolone acetate; 3b-Acetoxy-6,10-dimethyl-5,16-pregnadien-20-one
1552	5b-Pregnane-11b,17a,21-triol-3,20-dione; 11b,17a,21-Trihydroxy-5b-pregnane-3,20-dione	1999	6,16a-Dimethylpregnenolone 3-acetate; 3b-Acetoxy-6,16a-dimethyl-5-pregnen-20-one
1516	5b-Pregnane-11b,21-diol-3,20-dione; 11b,21-Dihydroxy-5b-pregnane-3,20-dione	1993	6,16a-Dimethylpregnenolone; 3b-Hydroxy-6,16a-dimethyl-5-pregnen-20-one
1606	5b-Pregnane-17a,21-diol-3,11,20-trione; 17a,21-Dihydroxy-5b-pregnane-3,11,20-trione	283	6,6-Difluoro-5a-androstan-17-one
1518	5b-Pregnane-17a,21-diol-3,20-dione; 17a,21-Dihydroxy-5b-pregnane-3,20-dione	2442	6-Acetoximino-5a-cholestan-3b-ol acetate; 3b-Acetoxy-5a-cholestan-6-one oxime acetate
1592	5b-Pregnane-3,11,20-trione	1975	6-Chloro-11b,17a,21-trihydroxypregna-1,4,6-triene-3,20-dione; Cloprednol
1593	5b-Pregnane-3,12,20-trione	98	6-Dehydroestrone; 3-Hydroxy-1,3,5(10),6-estratetraen-17-one
1490	5b-Pregnane-3,20-dione	766	6-Fluoro-5-androsten-17-one
1495	5b-Pregnane-3,20-dione semicarbazone	2586	6-Hydroximino-4-cholesten-3-one; 4-Cholestene-3,6-dione 6-oxime
1589	5b-Pregnane-3,6,20-trione	2460	6-Hydroximino-5-hydroxy-5x-cholestan-3-one; 5-Hydroxy-5x-cholestane-3,6-dione 6-oxime
1617	5b-Pregnane-3,7,12,20-tetrone	2419	6-Hydroximino-5a-cholestane; 5a-Cholestan-6-one oxime
1284	5b-Pregnane-3a,11b,17a,20a,21-pentol; a-Cortol	2449	6-Hydroximino-5x-cholestane-3b,5-diol; 3b,5-Dihydroxy-5x-cholestan-6-one oxime
1272	5b-Pregnane-3a,11b,17a,20a-tetrol	89	6-Ketoestradiol 6-(O-carboxymethyl)oxime
1273	5b-Pregnane-3a,11b,17a,20b-tetrol	81	6-Ketoestradiol; 3,17b-Dihydroxy-1,3,5(10)-estratrien-6-one
1274	5b-Pregnane-3a,11b,20b,21-tetrol	91	6-Ketoestriol 6-(O-carboxymethyl)oxime
1242	5b-Pregnane-3a,11b,20b-triol	92	6-Ketoestriol 6-(O-carboxymethyl)oxime:BSA conjugate
1338	5b-Pregnane-3a,11b-diol-20-one; 3a,11b-Dihydroxy-5b-pregnan-20-one	2027	6-Methyl-16-dehydropregnenolone acetate; 3b-Acetoxy-6-methyl-5,16-pregnadien-20-one
1364	5b-Pregnane-3a,12a-diol-20-one diacetate; 3a,12a-Diacetoxy-5b-pregnan-20-one	2007	6-Methyl-16a,17a-epoxypregnenolone acetate; 3b-Acetoxy-16a,17a-epoxy-6-methyl-5-pregnen-20-one
1244	5b-Pregnane-3a,17a,20a-triol	2005	6-Methyl-16a,17a-epoxypregnenolone; 16a,17a-Epoxy-3b-hydroxy-6-methyl-5-pregnen-20-one
1399	5b-Pregnane-3a,17a,20a-triol-11-one; 3a,17a,20a-Trihydroxy-5b-pregnan-11-one	2003	6-Methyl-17a-hydroxypregnenolone diacetate; 3b,17a-Diacetoxy-6-methyl-5-pregnen-20-one
1276	5b-Pregnane-3a,17a,20b,21-tetrol	2619	6-Methyldiosgenin acetate; 3b-Acetoxy-6a-methyl-5-spirostene
1246	5b-Pregnane-3a,17a,20b-triol	1995	6-Methylpregnenolone acetate; 3b-Acetoxy-6-methyl-5-pregnen-20-one
1522	5b-Pregnane-3a,17a-diol-11,20-dione; 3a,17a-Dihydroxy-5b-pregnane-11,20-dione	235	6a,11a,17b-Triacetoxy-5a-androstane; 5a-Androstane-6a,11a,17b-triol triacetate
1250	5b-Pregnane-3a,20a,21-triol	653	6a,11a-Diacetoxy-3,3-ethylenedioxy-5a-androstan-17-one
1208	5b-Pregnane-3a,20a-diol	1034	6a,11a-Diacetoxy-3a,5-cyclo-5a-androstan-17-one
1218	5b-Pregnane-3a,20a-diol 20-sulfate, sodium salt	1777	6a,11a-Dihydroxy-17a-pregn-4-ene-3,20-dione
1239	5b-Pregnane-3a,20a-diol 3-glucuronide, sodium salt	812	6a,11a-Dihydroxy-5a-androst-1-en-3-one
1217	5b-Pregnane-3a,20a-diol 3-sulfate, potassium salt		
1229	5b-Pregnane-3a,20a-diol diacetate; 3a,20a-Diacetoxy-5b-pregnane		
1219	5b-Pregnane-3a,20a-diol disulfate, dipotassium salt		
1238	5b-Pregnane-3a,20a-diol glucuronide		
1251	5b-Pregnane-3a,20b,21-triol		
1210	5b-Pregnane-3a,20b-diol		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
445	6a,11a-Dihydroxy-5a-androstan-16-one	318	6a-Hydroxy-5a-androstan-11-one; 5a-Androstan-6a-ol-11-one
460	6a,11a-Dihydroxy-5a-androstan-17-one	598	6a-Hydroxy-5a-androstane-1,17-dione
402	6a,11a-Dihydroxy-5a-androstan-2-one	612	6a-Hydroxy-5a-androstane-3,16-dione
404	6a,11a-Dihydroxy-5a-androstan-3-one	1497	6a-Hydroxy-5a-pregnane-3,20-dione
533	6a,12b,15a-Triacetoxo-5a-androstan-3-one	341	6a-Methoxy-5a-androstan-17-one
507	6a,12b,15a-Trihydroxy-5a-androstan-3-one	2031	6a-Methyl-16-dehydroprogesterone; 6a-Methyl-5,16-pregnadiene-3,20-dione (approx.85%)
406	6a,12b-Dihydroxy-5a-androstan-3-one	2016	6a-Methyl-17a-hydroxyprogesterone acetate; 17a-Acetoxy-6a-methyl-4-pregnene-3,20-dione
2012	6a,16a-Dimethylprogesterone; 6a,16-Dimethyl-4-pregnene-3,20-dione	2036	6a-Methylprednisolone; 6a-Methyl-11b,17a,21-trihydroxy-1,4-pregnadiene-3,20-dione
811	6a,16b-Dihydroxy-5a-androst-2-en-1-one	531	6b,11a,16b-Triacetoxo-5a-androstan-3-one
407	6a,16b-Dihydroxy-5a-androstan-3-one	532	6b,11a,17b-Triacetoxo-5a-androstan-3-one
1821	6a,17,21-Trihydroxy-4-pregnene-3,20-dione	496	6b,11a-Diacetoxy-16,16-difluoro-5a-androstan-3-one
1779	6a,17a-Dihydroxy-4-pregnene-3,20-dione	1803	6b,11a-Diacetoxy-4-pregnene-3,20-dione
409	6a,17a-Dihydroxy-5a-androstan-3-one	841	6b,11a-Diacetoxy-5a-androst-1-en-3-one; 6b,11a-Dihydroxy-5a-androst-1-en-3-one 6,11-diacetate
1011	6a,17b-Dihydroxy-3-methyl-5a-androstan-2-en-1-one	1529	6b,11a-Diacetoxy-5a-pregnane-3,20-dione
817	6a,17b-Dihydroxy-4-androsten-3-one	524	6b,11a-Dihydroxy-17b-methoxy-5a-androstan-3-one
814	6a,17b-Dihydroxy-5a-androst-1-en-3-one	815	6b,11a-Dihydroxy-4-androsten-3-one
440	6a,17b-Dihydroxy-5a-androstan-11-one	1778	6b,11a-Dihydroxy-4-pregnene-3,20-dione
410	6a,17b-Dihydroxy-5a-androstan-3-one	813	6b,11a-Dihydroxy-5a-androst-1-en-3-one
654	6a,19-Diacetoxy-2x,3x-epoxy-5a-androstane-1,17-dione	405	6b,11a-Dihydroxy-5a-androstan-3-one
887	6a,19-Diacetoxy-5a-androst-2-ene-1,17-dione	1045	6b,11a-Dihydroxy-A-nor(5->6b)abeo-androstan-5-one
1835	6a,21-Diacetoxy-17a-hydroxy-4-pregnene-3,20-dione	1861	6b,11b,17a,21-Tetrahydroxy-4-pregnene-3,20-dione; 6b-Hydroxycortisol
1961	6a,9a-Difluoro-11b,16a,17a,21-tetrahydroxy-1,4-pregnadiene-3,20-dione 16,17-acetonide 21-propionate; Procinonide	886	6b,15a-Dihydroxy-4-androstene-3,17-dione
1959	6a,9a-Difluoro-11b,16a,17a,21-tetrahydroxy-1,4-pregnadiene-3,20-dione 16,17-acetonide; Fluocinolone acetonide	816	6b,15b-Dihydroxy-4-androsten-3-one
2040	6a,9a-Difluoro-16a-methyl-11b,17a,21-trihydroxy-1,4-pregnadiene-3,20-dione; Flumethasone	408	6b,16b-Dihydroxy-5a-androstan-3-one
2041	6a,9a-Difluoro-16b-methyl-11b,17a,21-trihydroxy-1,4-pregnadiene-3,20-dione; Flumethasone, 16b-methyl epimer	1899	6b,17a,21-Trihydroxy-4-pregnene-3,11,20-trione; 6b-Hydroxycortisone (approx.95%)
2876	6a-Acetoxy-3b,9b-epoxy-5a-lumista-7,22-diene; 3b,9b-Epoxy-5a-lumista-7,22-dien-6a-ol acetate	1822	6b,17a,21-Trihydroxy-4-pregnene-3,20-dione
882	6a-Acetoxy-4-androstene-3,17-dione	1780	6b,17a-Dihydroxy-4-pregnene-3,20-dione
124	6a-Acetoxy-5a-androstane	21	6b,17b-Diacetoxy-4-estren-3-one
2905	6a-Acetoxy-5b-lumista-7,9(11),22-trien-3-one; 6a-Hydroxy-5b-lumista-7,9(11),22-trien-3-one acetate	818	6b,17b-Dihydroxy-4-androsten-3-one
2795	6a-Benzoyloxy-5a,9b-ergosta-7,22-diene-3b,5-diol; 5a,9b-Ergosta-7,22-diene-3b,5,6a-triol 6-benzoate	1836	6b,21-Diacetoxy-17a-hydroxy-4-pregnene-3,20-dione
2794	6a-Benzoyloxy-5a-ergosta-7,22-diene-3b,5-diol; 5a-Ergosta-7,22-diene-3b,5,6a-triol 6-benzoate	1480	6b,21-Diacetoxy-3b,5,17a-trihydroxy-5a-pregnan-20-one
2885	6a-Benzoyloxy-5b-lumista-7,22-diene-3b,5-diol; 5b-Lumista-7,22-diene-3b,5,6a-triol 6-benzoate	1584	6b,21-Diacetoxy-5,17a-dihydroxy-5a-pregnane-3,20-dione
2231	6a-Bromo-3a,12a-diacetoxy-7-oxo-5b-cholan-24-oic acid ethyl ester	1781	6b,21-Dihydroxy-4-pregnene-3,20-dione
2230	6a-Bromo-3a,12a-diacetoxy-7-oxo-5b-cholan-24-oic acid methyl ester	1776	6b,9a-Dihydroxy-4-pregnene-3,20-dione (crude)
2391	6a-Bromo-5a-cholestane-3b,7a-diol	1795	6b-Acetoxy-17a-hydroxy-4-pregnene-3,20-dione
134	6a-Chloroacetoxy-5a-androstane	1796	6b-Acetoxy-21-hydroxy-4-pregnene-3,20-dione
1957	6a-Fluoro-11b,16a,17a,21-tetrahydroxy-1,4-pregnadiene-3,20-dione 16,17-acetonide; Flunisolid	1032	6b-Acetoxy-3a,5-cycloandrostan-17-one
1864	6a-Fluoro-11b,16a,17a,21-tetrahydroxy-4-pregnene-3,20-dione 16,17-acetonide; Flurandrenolide	1459	6b-Acetoxy-3b,5,17a-trihydroxy-5a-pregnan-20-one
1953	6a-Fluoro-11b,17a,21-trihydroxy-1,4-pregnadiene-3,20-dione; 6a-Fluoroprednisolone	1405	6b-Acetoxy-3b,5-dihydroxy-5a-pregnan-20-one
1030	6a-Hydroxy-3a,5-cyclo-5a-androstan-17-one (crude)	2389	6b-Acetoxy-3b-chloro-5a-cholestan-5-ol; 3b-Chloro-5a-cholestane-5,6b-diol 6-acetate
463	6a-Hydroxy-3b-methoxy-5a-androstan-17-one	2345	6b-Acetoxy-3b-cyano-5-bromo-5a-cholestane; 3b-Cyano-5-bromo-5a-cholestan-6b-ol acetate
1738	6a-Hydroxy-4-pregnene-3,20-dione	883	6b-Acetoxy-4-androstene-3,17-dione
871	6a-Hydroxy-5a-androst-1-ene-3,17-dione	1751	6b-Acetoxy-4-pregnene-3,20-dione
867	6a-Hydroxy-5a-androst-2-ene-1,16-dione	1561	6b-Acetoxy-5,17a-dihydroxy-5a-pregnane-3,20-dione
868	6a-Hydroxy-5a-androst-2-ene-1,17-dione	125	6b-Acetoxy-5a-androstane
		2437	6b-Acetoxy-5a-cholestan-7-one; 6b-Hydroxy-5a-cholestan-7-one acetate
		2886	6b-Benzoyloxy-5b-lumista-7,22-diene-3b,5-diol; 5b-Lumista-7,22-diene-3b,5,6b-triol 6-benzoate
		2887	6b-Benzoyloxy-9a-lumista-7,22-diene-3b,5b-diol; 9a-Lumista-7,22-diene-3b,5b,6b-triol 6-benzoate
		2563	6b-Bromo-4-cholesten-3-one
		2648	6b-Bromo-5-oxo-A-nor-3,5-seco-cholestan-3-oic acid

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2392	6b-Bromo-5a-cholestane-3b,7b-diol (crude)	2539	7b-Thiocyanocholesteryl acetate
280	6b-Fluoro-5a-androstan-17-one	1170	8,19-Epoxy-19-methoxy-3-oxoeti-4-enic acid methyl ester
1161	6b-Hydroxy-3-oxoeti-4-enic acid	2108	8,19-Epoxy-3b,5,19-trihydroxy-5b,14b,17a-card-20(22)-enolide
1031	6b-Hydroxy-3a,5a-cycloandrostan-17-one	1134	8,19-Epoxy-5-hydroxy-19-methoxy-3-oxo-5b-etianic acid methyl ester
874	6b-Hydroxy-4-androstene-3,17-dione	1121	8,19-Epoxy-5-hydroxy-3-oxo-5b-etianic acid ethyl ester
324	6b-Hydroxy-5a-androstan-17-one; 5a-Androstan-6b-ol-17-one	1173	8-Hydroxy-3-oxo-17a-eti-4-ene-19,20-dioic acid 20-amide 19,8-lactone
305	6b-Hydroxy-5a-androstan-3-one; 5a-Androstan-6b-ol-3-one	1851	8b,19-Epoxy-19,21-dihydroxy-14b,17a-pregn-4-ene-3,20-dione
2238	7,12-Dioxo-3a-hydroxy-5b-cholan-24-oic acid	1850	8b,19-Epoxy-19,21-dihydroxy-4-pregnene-3,20-dione
2241	7,12-Dioxo-3a-methoxy-5b-cholan-24-oic acid methyl ester	1816	8b,19-Epoxy-19-hydroxy-17a-pregn-4-ene-3,20-dione
2237	7,12-Dioxo-5b-cholan-24-oic acid ethyl ester	1817	8b,19-Epoxy-19-methoxy-17a-pregn-4-ene-3,20-dione
284	7,7-Difluoro-5a-androstan-17-one	1484	8b,19-Epoxy-19-methoxy-3b,5,21-trihydroxy-5b-pregnan-20-one
281	7,7-Difluoro-5a-androstan-3-one	1853	8b,19-Epoxy-21-hydroxy-19-methoxy-14b,17a-pregn-4-ene-3,20-dione
238	7,7-Difluoro-5a-androstane-3b,12b,15a-triol	1852	8b,19-Epoxy-21-hydroxy-19-methoxy-4-pregnene-3,20-dione
579	7,7-Ethylenedioxy-5a-androstan-17-one; 5a-Androstane-7,17-dione 7-ethylene ketal	1774	8b,19-Epoxy-4-pregnene-3,20-dione
1493	7,7-Ethylenedioxy-5a-pregnan-20-one; 5a-Pregnane-7,20-dione 7-ethylene ketal	1585	8b,19-Epoxy-5,19,21-trihydroxy-5b-pregnane-3,20-dione
584	7,7:17,17-Bis(ethylenedioxy)-5a-androstane; 5a-Androstane-7,17-dione bis(ethylene ketal)	1579	8b,19-Epoxy-5-hydroxy-19-methoxy-5b-pregnane-3,20-dione
2588	7-Dehydrocholesterol; Cholesta-5,7-dien-3b-ol; Provitamin D3	1547	8b,19-Epoxy-5-hydroxy-5b-pregnane-3,20-dione
2420	7-Hydroximino-5a-cholestane; 5a-Cholestan-7-one oxime	2071	8b,21-Dihydroxy-19-nor-14b,17a-pregn-4-ene-3,20-dione
2225	7a,12a-Diacetoxy-3-oxo-5b-cholan-24-oic acid methyl ester	2070	8b,21-Dihydroxy-19-nor-4-pregnene-3,20-dione
2223	7a,12a-Dihydroxy-3-oxo-5b-cholan-24-oic acid methyl ester	2066	8b-Hydroxy-19-nor-17a-pregn-4-ene-3,20-dione
2163	7a,12a-Dihydroxy-5b-cholan-24-oic acid	2065	8b-Hydroxy-19-nor-4-pregnene-3,20-dione
1352	7a,12b-Dihydroxy-5a-pregnan-20-one (crude)	1022	9(11)-Methyltestosterone; 17b-Hydroxy-17a-methyl-4,9(11)-androstadien-3-one
412	7a,16b-Dihydroxy-5a-androstan-3-one	2945	9,10-Secoergosta-5,7,22-trien-3b-ol; Dihydrotachysterol
413	7a,17b-Dihydroxy-5a-androstan-3-one	399	9a,11a-Epoxy-5a-androstan-3-one
2749	7a,8a-Epoxy-5a-ergost-22-ene	1943	9a,11b-Dichloro-17,21-dihydroxy-1,4-pregnadiene-3,20-dione
508	7a,9a,11a-Trihydroxy-5a-androstan-3-one	1770	9a,11b-Dichloro-17a-acetoxypregesterone
644	7a,9a-Dihydroxy-5a-androstane-3,17-dione	1954	9a,11b-Dichloro-6a-fluoro-21-hydroxy-16a,17a-(isopropylidenedioxy)-1,4-pregnadiene-3,20-dione; Fluchlorolone
186	7a-Benzoyloxy-5a-androstan-3b-ol	642	9a,16b-Dihydroxy-5a-androstane-3,11-dione
2326	7a-Chloro-7b-nitroso-5a-cholestane	414	9a,17b-Dihydroxy-5a-androstan-3-one
966	7a-Hydroxy-4,4-dimethyl-5a-androstane-3,6,16-trione	596	9a-Bromo-3,3-ethylenedioxy-5a-androstan-11-one; 9-Bromo-5a-androstane-3,11-dione 3-ethylene ketal
325	7a-Hydroxy-5a-androstan-17-one; 5a-Androstan-7a-ol-17-one	1596	9a-Bromo-5a-pregnane-3,11,20-trione
301	7a-Hydroxy-5a-androstan-2-one; 5a-Androstan-7a-ol-2-one	2049	9a-Chloro-16b-methyl-11b,17,21-trihydroxy-1,4-pregnadiene-3,20-dione; Beclomethasone
306	7a-Hydroxy-5a-androstan-3-one; 5a-Androstan-7a-ol-3-one	2050	9a-Chloro-17a,21-dipropionyloxy-11b-hydroxy-16b-methyl-1,4-pregnadiene-3,20-dione; Beclomethasone dipropionate
616	7a-Hydroxy-5a-androstane-3,17-dione	954	9a-Fluoro-11b,17b-dihydroxy-17-methyl-5a-androstan-3-one
1305	7a-Hydroxy-5a-pregnan-20-one	1769	9a-Fluoro-11b-hydroxy-4-pregnene-3,20-dione
2568	7a-Methoxy-4-cholesten-3-one	2039	9a-Fluoro-16b-methyl-11b,17a,21-trihydroxy-1,4-pregnadiene-3,20-dione; Betamethasone
492	7b,11a-Diacetoxy-3a-fluoro-5a-androstan-16-one	962	9a-Fluoro-17b-hydroxy-17a-methyl-5b-androstane-3,11-dione
851	7b,11a-Dihydroxy-3b-fluoro-5-androsten-17-one	588	9a-Fluoro-5a-androstane-2,17-dione
411	7b,11a-Dihydroxy-5a-androstan-3-one	618	9a-Hydroxy-5a-androstane-3,17-dione
518	7b,12b,14b-Trihydroxy-5a-androstan-15-one	603	9a-Hydroxy-5a-androstane-3,7-dione
184	7b,12b-Diacetoxy-5a,14b-androstane; 5a,14b-Androstane-7b,12b-diol diacetate	2879	9a-Lumista-7,22-diene-3b,5b,6b-triol
461	7b,12b-Dihydroxy-5a-androstan-17-one	2867	9a-Lumista-7,22-diene-3b,5b-diol
641	7b,12b-Dihydroxy-5a-androstane-3,17-dione	959	9a-Methyl-5a-androstan-3-one-2-spiro-2'-(1',3'-dithian); 9a-Methyl-5a-androstane-2,3-dione 2-trimethylenedithio ketal
1056	7b,15a-Isopropylidenedioxy-D-homo-5a-androstan-12b-ol	928	9a-Methyl-5a-androstan-3b-ol
2541	7b-Benzoyloxy-5-cholesten-3b-ol; 5-Cholestene-3b,7b-diol 7-benzoate	2922	9a-Methyl-5a-ergosta-7,22-diene-3,11-dione
360	7b-Benzoyloxy-5a-androstan-3-one; 5a-Androstan-7b-ol-3-one benzoate	2916	9a-Methyl-5a-ergosta-7,22-diene-3b,11b-diol
382	7b-Hydroxy-5a-androstan-17-one ethylene ketal; 17,17-Ethylenedioxy-5a-androstan-7b-ol	1727	9b,10a-Pregn-4-ene-3,20-dione; Retroprogesterone
326	7b-Hydroxy-5a-androstan-17-one; 5a-Androstan-7b-ol-17-one	2155	9b,11b-Epoxy-5b-cholan-24-oic acid methyl ester
307	7b-Hydroxy-5a-androstan-3-one; 5a-Androstan-7b-ol-3-one	2800	9b-Ergosta-4,22-dien-3-one
617	7b-Hydroxy-5a-androstane-3,17-dione		
2569	7b-Methoxy-4-cholesten-3-one		

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2830	9b-Ergosta-5,7,22-trien-3-one (some decomposition)	1285	b-Cortol; 5b-Pregnane-3a,11b,17a,20b,21-pentol
2880	9b-Lumista-7,22-diene-3b,5b,6b-triol	2993	b-Escin (approx.95%)
2997	a-Amyrin 3,5-dinitrobenzoate; Urs-12-en-3b-ol 3,5-dinitrobenzoate	2962	b-Sitosteryl palmitate
2995	a-Amyrin acetate; Urs-12-en-3b-ol acetate; 3b-Acetoxy-12-ursene	2045	Betamethasone 17a-valerate
2673	A-Homo-4a-cholesten-3-one	2042	Betamethasone 21-disodium phosphate
2668	A-Homo-5a-cholestan-3-one	2988	Betulin; Lup-20(30)-ene-3b,28-diol
2666	A-Homo-5a-cholestan-4x-ol	2244	Bile acids, oxidized (from ox or sheep), mostly dehydrocholic acid
2665	A-Homo-5a-cholestane	1010	Bis(17-methyltestosterone N-(3-carboxyamino-2-hydroxypropyl)carbonylmethoxime) PEG ester
1050	A-Nor(5->6b)abeo-androst-2-en-5-one	1768	Bis(17a-hydroxyprogesterone 3/20-(N-(3-carboxyamino-2-hydroxypropyl)carbonylmethoxime) PEG ester
1041	A-Nor(5->6b)abeo-androstan-5-one	950	Bis(17a-methylandrostanolone N-(3-carboxyamino-2-hydroxypropyl)carbonylmethoxime) PEG ester
1049	A-Nor(5->6b)abeo-androstane-3,5,15-trione	2606	Bis(3-oxocholesta-1,4-dien-7-yl) diselenide
1047	A-Nor(5->6b)abeo-androstane-3,5-dione	2621	Botogenin acetate; 3b-Acetoxy-25(R)-spirost-5-en-12-one
1039	A-Nor(5->6b)abeo-androstane-3b,5b-diol	2277	Bufalin; 3b,14-Dihydroxy-5b-bufa-20,22-dienolide
2932	A-Nor-19-norergosta-5,7,9-trien-3-one	2663	Calciferol 2
1040	A-Nor-5a-androstan-2-one	2664	Calciferol 3; 9,10-Seco-5,7,10(19)-cholestatrien-3b-ol
1048	A-Nor-5a-androstane-2,12,15-trione	2911	Campesta-5,7-dien-3b-ol
2632	A-Nor-5a-cholestan-2-one	2175	Chenodeoxycholic acid diacetate methyl ester; 3a,7a-Diacetoxy-5b-cholan-24-oic acid methyl ester
2631	A-Nor-5a-cholestan-2b-ol	1939	Chlormadinone acetate; 17a-Acetoxy-6-chloro-4,6-pregnadiene-3,20-dione
2634	A-Nor-5a-cholestane-2a-carboxylic acid methyl ester	2609	Cholesta-1,3,5-trien-7-one
2925	A-Nor-5a-ergost-7-en-2-one	2608	Cholesta-1,4,6-trien-3-one
2633	A-Nor-5b-cholestan-3-one	2600	Cholesta-1,4-dien-3-one
2926	A-Nor-5b-lumist-7-en-2-one	2660	Cholesta-1,4-dien-3-one 4,5-monoozonide
2928	A-Nor-5b-lumista-7,22-dien-2-one	2602	Cholesta-3,5-dien-7-one
2927	A-Nor-5b-lumista-7,22-diene	2587	Cholesta-3,5-diene
2116	Acetylstrophanthidin; 3b-Acetoxy-5a,14-dihydroxy-19-oxo-20(22)-cardenolide	2601	Cholesta-4,6-dien-3-one
2107	Acovenoside A; 3b-(6-Deoxy-3-O-methyl-a-L-talopyranosyl)-1b,14-dihydroxy-5b-card-20(22)-enolide	2603	Cholesta-4,6-dien-3-one semicarbazone
1898	Aldosterone 18,21-diacetate 3-(O-carboxymethyl)oxime:BSA conjugate	2589	Cholesta-5,22(E)-dien-3b-ol; 22(E)-Dehydrocholesterol
875	Androst-4-en-11b-ol-3,17-dione; 11b-Hydroxyandrost-4-ene-3,17-dione	2590	Cholesta-5,22(Z)-dien-3b-ol; 22(Z)-Dehydrocholesterol
878	Androst-4-en-16a-ol-3,17-dione; 16a-Hydroxyandrost-4-ene-3,17-dione	2591	Cholesta-5,24-dien-3b-ol; Desmosterol (approx.95%)
894	Androst-4-en-19-al-3,17-dione	2599	Cholesta-5,7-dien-3b-ol 3',5'-dinitrobenzoate
1140	Androst-5-en-3b-ol-17b-carboxylic acid; 3b-Hydroxyeti-5-enic acid	2399	Cholestane-3b,5a,6b-triol
734	Androst-5-ene-3b,11b,17b-triol	2396	Cholesterol-5a,6a-epoxide; 5a,6a-Epoxycholestan-3b-ol
735	Androst-5-ene-3b,16a,17a-triol	2470	Cholesterol; 5-Cholesten-3b-ol
839	Androst-5-ene-3b,16a-diol-17-one 3-acetate; 3b-Acetoxy-16a-hydroxyandrost-5-en-17-one (approx.95-98%)	2504	Cholesteryl 10-undecenoate
737	Androst-5-ene-3b,16b,17b-triol	2520	Cholesteryl 11-eicosenoate
731	Androst-5-ene-3b,17b-diol 17-benzoate; 17b-Benzoyloxyandrost-5-en-3b-ol	2477	Cholesteryl acetate; 3b-Acetoxy-5-cholestene
732	Androst-5-ene-3b,17b-diol 3-acetate 17-benzoate; 3b-Acetoxy-17b-benzoyloxyandrost-5-ene	2519	Cholesteryl arachidate; Cholesteryl eicosanoate
924	Androsta-1,4,9(11)-triene-3,17-dione	2521	Cholesteryl arachidonate; Cholesteryl eicosatetraenoate
921	Androsta-1,4-dien-11b-ol-3,17-dione; 11b-Hydroxyandrosta-1,4-diene-3,17-dione	2522	Cholesteryl behenate; Cholesteryl docosanoate
374	Androsterone glucuronide; 3a-Hydroxy-5a-androstan-17-one gluconuride	2527	Cholesteryl benzoate
333	Androsterone sulfate, sodium salt; 3a-Hydroxy-5a-androstan-17-one 3-sulfate, sodium salt	2465	Cholesteryl bromide; 3b-Bromo-5-cholestene
321	Androsterone; 3a-Hydroxy-5a-androstan-17-one	2488	Cholesteryl butyrate
2261	Apochoolic acid; 3a,12a-Dihydroxy-5b-chole-8(14)-en-24-oic acid	2464	Cholesteryl chloride; 3b-Chloro-5-cholestene
3002	b-Amyrin 3,5-dinitrobenzoate; Olean-12-en-3b-ol 3,5-dinitrobenzoate	2480	Cholesteryl chloroacetate
3000	b-Amyrin acetate; 3b-Acetoxy-12-oleanene; Olean-12-en-3b-ol acetate	2502	Cholesteryl cinnamate
2998	b-Amyrin; b-Amyrenol; Olean-12-en-3b-ol	2490	Cholesteryl crotonate; Cholesteryl 2-butenate
		2503	Cholesteryl decanoate; Cholesteryl caprate (95-97%)
		2523	Cholesteryl erucate; Cholesteryl (Z)-13-docosenoate
		2475	Cholesteryl formate
		2492	Cholesteryl hemisuccinate, monocyclohexylamine salt
		2493	Cholesteryl hemisuccinate, morpholine salt
		2494	Cholesteryl hemisuccinate, tris salt
		2511	Cholesteryl heptadecanoate
		2497	Cholesteryl hexanoate; Cholesteryl caproate
		2501	Cholesteryl hydrocinnamate; Cholesteryl 3-phenylpropionate

<u>No.</u>	<u>Compound Description</u>	<u>No.</u>	<u>Compound Description</u>
2466	Cholesteryl iodide; 3b-Iodo-5-cholestene	780	Dehydroisoandrosterone 3-sulfate, sodium salt; 3b-Hydroxyandrost-5-en-17-one 3-sulfate, sodium salt
2495	Cholesteryl isoamyl ether; 3b-(3-Methylbutoxy)-5-cholestene (approx.80%)	777	Dehydroisoandrosterone; 3b-Hydroxyandrost-5-en-17-one
2489	Cholesteryl isobutyrate	2630	Demissidine; Solanine D; 16,22,26-Nitrilocholestan-3-ol
2505	Cholesteryl laurate	2178	Deoxycholic acid 3-(b-D-glucuronide), dipotassium salt
2515	Cholesteryl linoleate; Cholesteryl (Z,Z)-9,12-octadecadienoate	2164	Deoxycholic acid methyl ester; 3a,12a-Dihydroxy-5b-cholan-24-oic acid methyl ester
2516	Cholesteryl linolelaidate; Cholesteryl 9,12-octadecadienoate	2161	Deoxycholic acid, sodium salt; 3a,12a-Dihydroxy-5b-cholan-24-oic acid, sodium salt
2517	Cholesteryl linolenate; Cholesteryl 9,12,15-octadecatrienoate	1772	Deoxycorticosterone 21-hemisuccinate:BSA conjugate
2479	Cholesteryl methyl carbonate	1759	Deoxycorticosterone 21-hemisuccinate; 21-Hydroxy-4-pregnene-3,20-dione 21-hemisuccinate
2474	Cholesteryl methyl ether; 3b-Methoxy-5-cholestene	1754	Deoxycorticosterone acetate; 21-Acetoxy-4-pregnene-3,20-dione
2524	Cholesteryl nervonate; Cholesteryl (Z)-15-tetracosenoate	1068	Des-A-10b(H)-androstane-5,17-dione
2500	Cholesteryl nonanoate; Cholesteryl pelargonate	2098	Des-A-17a-pregnane-5,20-dione
2499	Cholesteryl octanoate; Cholesteryl caprylate	2679	Des-A-cholest-6-en-5-one 2,4-dinitrophenylhydrazone
2513	Cholesteryl oleate; Cholesteryl (Z)-9-octadecenoate	2680	Des-A-cholest-9(10)-en-5-one, 2,4-dinitrophenylhydrazone
2514	Cholesteryl oleyl carbonate (approx.95%)	2947	Des-A-ergost-22-en-5-one
2508	Cholesteryl palmitate; Cholesteryl hexadecanoate	2948	Des-A-lumist-22-en-5-one
2510	Cholesteryl palmitelaidate; Cholesteryl (E)-9-hexadecenoate	2949	Des-A-lumist-22-en-5-one oxime
2509	Cholesteryl palmitoleate; Cholesteryl (Z)-9-hexadecenoate	2096	Des-A-pregnan-20b-hydroxy-5-one
2507	Cholesteryl pentadecanoate	2097	Des-A-pregnane-5,20-dione
2486	Cholesteryl propyl carbonate	2099	Des-A-pregnane-5,20-dione bis(2,4-dinitrophenylhydrazone)
2472	Cholesteryl sulfate, potassium salt	2038	Dexamethasone; 9a-Fluoro-16a-methyl-11b,17a,21-trihydroxy-1,4-pregnadiene-3,20-dione
2473	Cholesteryl sulfate, sodium salt; 5-Cholesten-3b-ol sulfate, sodium salt	3007	Diethylstilbestrol
2506	Cholesteryl tridecanoate	2106	Digitoxigenin; 3,14-Dihydroxy-5b-card-20(22)-enolide
2496	Cholesteryl valerate; Cholesteryl pentanoate	2103	Digitoxin; 3b-(Digitoxosyl(1->4)digitoxosyl(1->4)digitoxosyloxy)-14-hydroxy-5b-card-20(22)-enolide
2469	Cholesterylacetanilide; N-Acetyl-N-phenyl-5-cholesten-3b-amine	2329	Dihydrocholesterol; 5a-Cholestan-3b-ol (95-97%)
2468	Cholesterylaniline; N-Phenyl-5-cholesten-3b-amine	2617	Diosgenin; 25R-Spirost-5-en-3b-ol
2118	Convallatoxin; 3b-(6-Deoxy-a-L-mannopyranosyloxy)-5a,14-dihydroxy-19-oxo-5b-card-20(22)-enolide; Strophanthidin a-L-rhamnoside	3008	dl-2,3-Bis(3'-hydroxybenzyl)butane-1,4-diol; Enterodiol
2411	Coprostan-3-one; 5b-Cholestan-3-one	27	ent-Estradiol
2323	Coprostane; 5b-Cholestane	2996	Epi-a-amyrin acetate; Urs-12-en-3a-ol acetate; 3a-Acetoxy-12-ursene
1787	Cortisolone; 17a,21-Dihydroxy-4-pregnene-3,20-dione; Reichstein's S	2994	Epi-a-amyrin; Urs-12-en-3a-ol
1896	Corticosteron-18-oic acid 11,18-lactone; 11b,21-Dihydroxy-3,20-dioxopregn-4-en-18-oic acid 11,18-lactone	3001	Epi-b-amyrin acetate; 3a-Acetoxy-12-oleanene; Olean-12-en-3a-ol acetate
1893	Cortisone 21-glucuronide methyl ester	2999	Epi-b-amyrin; Olean-12-en-3a-ol
1894	Cortisone 21-glucuronide methyl ester 3,20-disemicarbazone	335	Epiandrosterone sulfate, sodium salt; 3b-Hydroxy-5a-androstan-17-one 3-sulfate, sodium salt
1895	Cortisone 21-glucuronide triacetyl methyl ester 3,20-disemicarbazone	323	Epiandrosterone; 3b-Hydroxy-5a-androstan-17-one
1892	Cortisone 21-glucuronide, barium salt	100	Equilin methyl ether; 3-Methoxy-1,3,5(10),7-estratetraen-17-one
1882	Cortisone 21-sulfate, triethylammonium salt	99	Equilin; 3-Hydroxy-1,3,5(10),7-estratetraen-17-one
1884	Cortisone acetate; 21-Acetoxy-17a-hydroxypregn-4-ene-3,11,20-triene; Kendall's E; Reichstein's Fa	3010	Equol diacetate
1897	d-Aldosterone 21-hemisuccinate:BSA conjugate	3011	Equol dibenzoate
1888	d-Aldosterone 21-hemisuccinate; 11b,18-Epoxy-18,20-dihydroxypregn-4-ene-3,20-dione 21-hemisuccinate	2946	Ergocalciferol; Vitamin D2; Calciferol; 9,10-Secoergosta-5,7,10(19),22-tetraen-3b-ol
1169	d-Aldosterone g-lactone; 11b,18-Epoxy-18-hydroxy-3-oxoandrost-4-ene-17b-carboxylic acid g-lactone	2799	Ergosta-4,22-dien-3-one
1880	d-Aldosterone; 11b,18-Epoxy-18,21-dihydroxypregn-4-ene-3,20-dione; Reichstein's X	2829	Ergosta-5,7,22-trien-3-one; Ergosterone (crude)
101	d-Equilenin; 3-Hydroxy-1,3,5(10),6,8-estrápentaen-17-one	2816	Ergosta-7,9(11),22-trien-5a-ol
1057	D-Homo-5a-androstan-11-one	2831	Ergosterone 22,23-dibromide; 22,23-Dibromoergosta-5,7-dien-3-one
2057	Danazol; 17b-Hydroxy-17a-pregna-2,4-dien-20-yno[2,3-d]isoxazole	42	Estradiol phenylpropionate (assumed 17b)
2245	Dehydrocholic acid, sodium salt; 3,7,12-Trioxo-5b-cholan-24-oic acid, sodium salt	63	Estrilol 16a-glucuronide
2243	Dehydrocholic acid; 3,7,12-Trioxo-5b-cholan-24-oic acid	66	Estrilol 17b-(b-D-glucuronide), sodium salt; 1,3,5(10)-Estratriene-3,16a,17b-triol 17-glucuronide, sodium salt
785	Dehydroisoandrosterone 3-acetate; 3b-Acetoxyandrost-5-en-17-one	65	Estrilol 17b-glucuronide
		68	Estrilol 3-benzyl ether; 3-Benzyloxy-1,3,5(10)-estratriene-16a,17b-diol
		67	Estrilol 3-glucuronide pentaacetate, sodium salt

<u>No.</u>	<u>Compound Description</u>
54	Estriol 3-methyl ether; 3-Methoxy-1,3,5(10)-estratriene-16a,17b-diol
50	Estriol 3-sulfate, sodium salt + 33% methylglucamine

<u>No.</u>	<u>Compound Description</u>
60	Estriol triacetate; 3,16a,17b-Triacetoxy-1,3,5(10)-estratriene
51	Estriol trisulfate, trisodium salt (approx.95%)