

Inner Diameter x Outer Diameter Sizes Available	Pump Size*	Default Silicone	Clear C-Flex	Opaque C-Flex	ADCF C-Flex	PVC	PharMed	Pharma-Pure
1/16" x 3/16" (1.6 mm x 4.8 mm)	14	X	X	X	X	X	X	
1/8" x 1/4" (3.2 mm x 6.4 mm)	16	X	X	X	X	X	X	X
1/8" x 3/16" (3.2 mm x 4.8 mm)						X		
3/16" x 5/16" (4.8 mm x 7.9 mm)	25	X	X	X	X	X	X	X
3/16" x 3/8" (4.8 mm x 9.5 mm)	15	X	X	X	X	X	X	
1/4" x 3/8" (6.4 mm x 9.5 mm)	17	X	X		X	X	X	X
1/4" x 7/16" (6.4 mm x 11.1 mm)	24	X	X	X	X	X	X	
1/4" x 1/2" (6.4 mm x 12.7 mm)	26	X	X	X	X	X	X	X
5/16" x 7/16" (7.9 mm x 11.1 mm)	18	X	X	X	X	X	X	X
5/16" x 1/2" (7.9 mm x 12.7 mm)		X						X
3/8" x 1/2" (9.5 mm x 12.7 mm)		X	X	X	X	X		X
3/8" x 5/8" (9.5 mm x 15.9 mm)	73	X	X	X	X	X	X	X
3/8" x 3/16" (9.5 mm x 14.3 mm)		X						
1/2" x 5/8" (12.7 mm x 15.9 mm)				X		X		
1/2" x 3/4" (12.7 mm x 19 mm)	82	X	X		X	X	X	X
1/2" x 11/16" (12.7 mm x 17.5 mm)	19	X	X	X	X	X		X
1/2" x 7/8" (12.7 mm x 22.2 mm)		X						
5/8" x 7/8" (15.9 mm x 22.2 mm)			X		X	X		X
3/4" x 1" (19 mm x 25.4 mm)	90	X	X		X	X	X	X
3/4" x 1 1/8" (19 mm x 28.6 mm)		X	X		X			
3/4" x 1 1/4" (19 mm x 31.8 mm)			X	X	X			
7/8" x 1 1/8" (22.2 mm x 28.6 mm)		X	X		X			
7/8" x 1 1/4" (22.2 mm x 31.8 mm)		X						
1" x 1 1/4" (25.4 mm x 31.8 mm)		X	X		X	X		
1" x 1 3/8" (25.4 mm x 34.9 mm)	92	X					X	

\*この表のポンプサイズはMasterflex™チュービングポンプシステムのもので、内径/外径寸法はMasterflexブランドのものとは異なる場合があります。試験を行い、実際のアプリケーションに適していることを確認してください。MasterflexはCole-Parmer Instrument Companyの登録商標です。